Mini-preparation of plasmid DNA by alkaline lysis

Unit 1.6, Ausubel et al. 1999 (NUF 00.5.2; modified j.a.m. 07/07)

- 1. Inoculate 1.5-5 ml LB amended with appropriate antibiotic(s) with a colony taken from a selective plate. Grow until late exponential phase (until dense), 7-12 hours depending on antibiotic and growth temperature.
- 2. Microcentrifuge 1.5 ml of cell culture for 30 s 1 min. in an eppendorf tube and remove supernatant. The pellet may be stored at -20°C .
- 3. Resuspend the pellet in 100 µl of Solution I by vigorous vortexing (make sure no clumps of cells are visible) and let stand for 5 min. at room temperature.
- 4. Add 200 μl lysis solution (Solution II), mix gently (by inversion), and incubate for 10 min. at room temperature.
- 5. Add 150 μl of potassium acetate solution (Solution III), vortex briefly, and incubate 10 min. on ice.
- 6. Microcentrifuge for 4 min. at room temperature and transfer the supernatant to a fresh tube.
- 7. Add 2-4 µl of 10 mg/ml RNAse A and incubate at 37°C for 30-45 min.
- 8. Add 400 µl chloroform-isoamyl alcohol (24:1 by vol) and vortex. Microcentrifuge for 5 min. at room temperature and transfer the upper (aqueous) phase to a fresh tube.
- 9. Add 800 µl of 95% ethanol or 400 µl isopropanol, vortex, and microcentrifuge for 10 min at room temperature. Discard the supernatant.
- 10. Wash pellet with 700 μl of 70% ethanol (do not resuspend), microcentrifuge 3 min. at room temperature, and discard the supernatant.
- 11. Dry the DNA.
- 12. Resuspend in 20-50 µl TE or sddH2O.

The typical yield of high-copy-number plasmids such as pUC is about 3-5 µg DNA per ml of original cell culture.

Birnboim & Doly I

50 mM glucose 25 mM Tris-HCl 10 mM EDTA pH 8.0

Autoclave and store at room temperature.

Birnboim & Doly II

Prepare 2 solutions:

0.4 M NaOH

2% w/v SDS

Store at room temperature. Mix 1:1 immediately before use.

Birnboim & Doly III

294 g potassium acetate (3M) 50 ml 90% formic acid OR 115 ml glacial acetic acid. ddH₂O to 1 L pH should be \sim 5.5

Store at room temperature or at 4°C.