Protocol: Hygromycin B Selection (FAK null cells)

Materials:

- Hygromycin B (Invitrogen cat# 10687-010) 50 mg/ml
- Complete medium: aMEM w/out nucleosides (Invitrogen cat# 12561) + 15 % Fetal Bovine Serum + 1% Pen/Strep
- Syringe (sterile) / 0.22 um syringe filter (sterile)
- Sterile centrifuge tubes
- Tube rack, gloves, 70 % ethanol, etc. other general items used while working in the sterile hood

Procedure:

- 1. Sterile filter and aliquot Hygromycin B stock solution
- 2. Prepare medium containing Hygromycin B prior to each use (personal preference)
 - add 350 ul of Hygromycin B (50 mg/ml) to 49.65 mL of complete media (final concentration of 350 ug/ml FAK -/- cell specific)
- 3. Mix by pipetting up and down before adding to cells.
- 4. Change media on cells under hygromycin B selection every 3 days (hygromycin B degrades @ 37°C incubator)

Note: Hygromycin B and alph-MEM without nucleosides may be purchased from BioStores on the Stanford campus