Protocol: Actin Staining

Materials:

General:

Cells: NIH/3T3 fibroblasts seeded @1,000 to 3,000 cells/cm² 24 to 48 hrs prior to staining

Coverslips Timer
Kimwipes Waste beaker

Aluminum Foil 50 mL centrifuge tube ("formaldehyde waste")

Pipets/Pipet aid Markers
Pipet tips/pipetters Gloves

Aspirator

Reagents/ Solutions:

Phosphate-buffered saline (PBS), pH 7.4

3.7 % Formaldehyde in PBS (use methanol-free formaldehyde)

0.1 % Triton X-100 solution (in PBS) PAP pen (Sigma Aldrich Cat # Z377821)

Vectashield mounting media with DAPI – (Vector Laboratories Cat # H-1200)

Clear fingernail polish

Primary Blocking Solution:

PBS + 1% (w/v) BSA + 0.1% (v/v) Nonidet P40

Alexa Fluor 488 Phalloidin Solution: **protect from light – (**Invitrogen Cat # A12379)

5 uL Alexa Fluor 488 Phalloidin in 200 uL Primary Blocking Solution (for each coverslip area)

Procedure:

- 1. Wash cells 2X with PBS (10 mL per wash) remove all PBS from around slide after final wash
- 2. Fix cells: apply 2 mL of 3.7% formaldehyde solution on top of slide for 10 min. @ room temp.

 Note: formaldehyde is toxic; avoid contact with skin, eyes, etc.; dispose of as hazardous waste, do not pour down sink
- 3. Wash cells 2X with PBS (10mL per wash) place waste in the "formaldehyde waste" tube
- 4. Permeablize cells: apply 2 mL of 0.1% Triton X-100 solution on top of slide for 5 min. @ room temp.
- 5. Wash cells 2X with PBS (10 mL per wash)
- 6. Make a PAP pen region the on the slide see instructions/pattern below
- 7. Wash coverslip area 1X with PBS (~200 uL should cover area)
- 8. Actin staining: apply 200 uL of Alexa Fluor 488 Phalloidin Solution for 20 min. @ room temp
 - * **Protect from light** for all subsequent steps place aluminum foil over culture dish to prevent photobleaching of the *Alexa Fluor 488*
- 9. Remove the *Phalloidin Solution* and wash the coverslip area 3X with *Primary Blocking Solution*
- 10. Remove the *Blocking Solution* and add Vectashield mounting medium (25uL per coverslip area)
- 11. Place a coverslip over the stained area and seal the coverslip by pipetting a bead of clear fingernail polish along the edge of the coverslip.
- 12. Allow the nail polish to dry before viewing on a fluorescent microscope.

- i. Remove the slide from the culture dish and dry the bottom of the slide with a Kimwipe
- ii. Place the slide on the pattern
- iii. Use a Kimwipe to wipe dry the area in gray
- do not disturb the white coverslip area you will be staining
- iv. Apply the PAP pen along the black square around the coverslip area
- v. Dry culture dish, and place slide back in the dish

