Protocol: Transfecting NIH 3T3s with GFP-Actin

Materials:

- 3T3 Fibroblasts (seeded in 35 mm dishes or a 6-well plate)
- FuGENE 6 transfection reagent (Roche Applied Science Cat# 11 815 091 001)
- pAcGFP1-Actin vector (Clontech Cat.# 632453)
- DMEM (Invitrogen #11995-065)
- Water bath
- Serological Pipets
- Pipet aid
- Pipette tips

- 20 μL, 200 μL, 1000 uL pipetter
- Microscope
- Sterile polystyrene tube
- Timer
- Waste beaker
- Biohazard Bag
- 70% ethanol
- Kimwipes
- Markers
- Gloves

Procedure:

Adding serum-free DMEM to cells

- 1. Heat DMEM to 37°C in a water bath
- 2. Prepare the hood:
 - spray and wipe the hood surface with 70% ethanol
 - place the following materials into the hood: waste beaker, tube rack, polystyrene tube, pipetting materials, pre-warmed DMEM
 - tape a biohazard bag to the front of the hood
- 3. Observe the cells under microscope check for contamination and note cell confluence
- 4. Place the cells in the hood
- 5. Remove the cell culture media, and add 2 mL of DMEM (without serum or Pen/Strep)
- 6. Label the dish (initials) and place the dish back in the incubator

Diluting FuGENE 6 with DMEM

- 7. Aseptically place FuGENE 6 and GFP-Actin DNA in the hood
- 8. Pipette 91 uL of DMEM into a sterile polystyrene tube
- 9. Pipette 9 uL of FuGENE 6 directly into the DMEM
 - Avoid direct contact between FuGENE and the tube wall
- 10. Tap to mix
- 11. Incubate for 5 min @ room temp

Forming FuGENE/DNA complex

- 12. Add 5.3 uL of pAcGFP1-Actin vector (1 ug) to the dilute FuGENE
- 13. Tap to mix
- 14. Incubate for 15 min @ room temp
- 15. Remove the dish of cells from the incubator
- 16. Add the FuGENE/DNA complex to the cells in a drop-wise fashion
- 17. Gently swirl dish to mix
- 18. Label the dishes (GFP-Actin, date)
- 19. Place the cells back in the incubator (add serum containing medium in 6-24 hrs)
- 20. Clean the hood and place all material that contacted cells/medium in a biohazard waste bag; clean the waste beaker by adding bleach for ~15 min. before washing