## Useful Numbers for Cell Culture

There are various sizes of dishes and flasks used for cell culture. Some useful numbers such as surface area and volumes of dissociation solutions are given below for various size culture vessels.

|                | Surface Area<br>(mm²) | Seeding<br>Density | Cells at<br>Confluency¹ | Versene<br>(ml of 0.53 mM<br>EDTA) | Trypsin<br>(ml of 0.05%<br>tryspin, 0.53 mM<br>EDTA) | Growth<br>Medium<br>(ml) |
|----------------|-----------------------|--------------------|-------------------------|------------------------------------|--|--------------------------|
| Dishes         |                       |                    |                         |                                    |  |                          |
| 35 mm          | 962                   | $0.3 	imes 10^6$   | $1.2 \times 10^6$       | 1                                  | 1  | 2                        |
| 60 mm          | 2,827                 | $0.8 	imes 10^6$   | $3.2 	imes 10^6$        | 3                                  | 2  | 3                        |
| 100 mm         | 7,854                 | $2.2 \times 10^6$  | $8.8 	imes 10^6$        | 5                                  | 3  | 10                       |
| 150 mm         | 17,671                | $5.0 	imes 10^6$   | $20.0 	imes 10^6$       | 10                                 | 8  | 20                       |
| Cluster Plates |                       |                    |                         |                                    |  |                          |
| 6-well         | 962                   | $0.3 	imes 10^6$   | $1.2 	imes 10^6$        | 2                                  | 2  | 3–5                      |
| 12-well        | 401                   | $0.1 	imes 10^6$   | $0.4 	imes 10^6$        | 1                                  | 1  | 1–2                      |
| 24-well        | 200                   | $0.05 	imes 10^6$  | $0.2 	imes 10^6$        | 0.5                                | 0.5  | 0.5-1.0                  |
| Flasks         |                       |                    |                         |                                    |  |                          |
| T-25           | 2,500                 | $0.7 	imes 10^6$   | $2.8 	imes 10^6$        | 3                                  | 3  | 3–5                      |
| T-75           | 7,500                 | $2.1 \times 10^6$  | $8.4 	imes 10^6$        | 5                                  | 5  | 8–15                     |
| T-160          | 16,000                | $4.6 \times 10^6$  | $18.4 \times 10^6$      | 10                                 | 10   | 15–30                    |

 $<sup>^1</sup>$ The number of cells on a confluent plate, dish, or flask will vary with cell type. For this table, HeLa cells were used.