Review for the "ga5-1 complementation" section.

A. Choose all the options that you think are correct. You may want to choose more than one option for some questions.

| 1. Gibberellins are: |
|--|
| Encoded in the genome of plants Synthesized via the terpenoid pathway Produced by certain bacteria Produced by certain fungi Produced by plants. |
| 2. What do you need to generate cDNA from Arabidopsis? |
| Total DNA Total RNA Messenger RNA Genomic DNA |
| 3. What template would you use to amplify the CDS of GA20OX1 gene of Arabidopsis by PCR? |
| RNAGenomic DNAcDNAThe CDS of GA20OX1 |
| 4. You amplified the GA20OX1 gene from Arabidopsis using primers that introduced AscI and PacI restriction sites to the PCR product. After digesting the PCR product using AscI and PacI you run an agarose gel. |
| You run this restriction gel to see whether the PCR product was properly cut. You run the gel to see whether you had a PCR product. You run the gel to see whether the PCR product was the correct size. You run the PCR product on the gel in order to remove the restriction enzymes. |
| 5. What is the best way of making more copies of a plasmid? |
| Use PCR to make more copies of the DNA Transform it into E. coli Transform it into plants. |
| 6. What do you use agarose gel electrophoresis for? |
| Visualize the DNA your are working on. Analyze the length of the DNA you are working with by running it beside a standard of known size. Cut your DNA using restriction enzymes. Separate fragments of DNA of different sizes. Amplify a fragment of DNA by PCR. |

| 7. To sequence a fragment of DNA from a plasmid you need to: |
|--|
| Know some of the sequence of the plasmid in order to design primers flanking the region that needs to be sequenced. Cut the DNA using restriction enzymes. Purify your plasmid. Know what is the concentration of your plasmid preparation. Previously run your plasmid on an agarose gel. |
| 1 Teviously run your plasmid on an agarose get. |
| 8. In class you transformed your ligation reaction into E. coli cells. The transformed cells were able to grow on medium with the antibiotic kanamycin. Which of the following sentences explains why this was the case? |
| The transformation corrected a mutation in the E. coli genome that made the bacteria susceptible to the antibiotic The pMDC32 vector encodes for an enzyme that degrades kanamycin A gene that encodes for an enzyme that degrades kanamycin was inserted into the bacterial genome. |
| 9. We transformed the <i>ga5-1 Arabidopsis</i> mutant using a binary vector containing a wild type GA20OX1 gene to |
| Replace the mutated GA20OX gene found in the <i>ga5-1</i> mutant with the wild type GA20OX1 using Agrobacterium. Introduce the wild type GA20OX1 in the genome of <i>ga5-1</i> . Test whether this wild type GA20OX1 could complement the <i>ga5-1</i> semidwarf phenotype. |
| 10. What substrates does the DNA polymerase use? |
| Nucleotide triphosphates Deoxy nucleotide triphosphates Nucleosides Deoxynucleotide monophosphates |
| 11. During the "Green Revolution" high yield varieties of wheat and rice where developed. The researchers |
| Realized that varieties with shorter stems were less prone to lodging. Were trying to understand the role of gibberellins in plant growth. Were looking to optimize the flowering time of these plants. Were looking for varieties that produced more seeds. |
| 12. You used the vector pMDC32 in class. pMDC32 is a binary vector. This vector |
| Can replicate in both E. coli and Agrobacterium. Encodes for the proteins needed for plant transformation. Has two origins of replication. Is inserted twice in the Agrobacterium genome. |

____ Is part of the binary vector system in Agrobacterium that is formed of two plasmids, one that contains the T-DNA, and one, the helper plasmid that encodes for the VIR genes needed to transfer the T-DNA to the plant.

B. You want to clone the genomic fragment of your favorite gene (YFG) from Arabidopsis into the binary vector pMDC32. Below is a list of the steps you need to do. Organize them in the right order. Some of the procedures described here might have to occur at the same time, you can label those 2a and 2b for example. Remember that you have done most of the steps yourself in class.

a-Sequence the insert in your newly generated plasmid to test whether its sequence matches that of YFG. b-Check whether the E. coli cells transformed with the ligation reaction containing both the cut pMDC32 vector and your cut PCR product generated more colonies than the cells transformed with the ligation reaction that contained only your cut vector.

c-Set up a restriction digestion for both the pMDC32 vector and your purified PCR product.

d-Purify the PCR product to remove the DNA polymerase and the nucleotides.

e-Set up a ligation reaction using only your cut pMDC32 vector.

f-Plan your cloning strategy by looking for appropriate restriction sites in the pMDC32 vector

g-Design primers to amplify the genomic section of YFG and insert restriction sites to allow cloning into the pMDC32 vector.

h-Transform your ligation reaction into competent E. coli cells.

i-Run both your cut pMDC32 vector and your cut PCR product on an agarose gel.

j-Use the E. coli colonies on the LB medium containing kanamycin to inoculate a 3 ml overnight culture.

k-Set up a restriction digestion reaction of your new plasmids to test whether any of them contains your cut PCR product.

l-If you identify a new plasmid containing an insert of the correct size you will prepare this plasmid for sequencing.

m-See whether any of the new plasmids contain an insert of the size of your PCR product.

n-PCR YFG using your total DNA preparation and the primers you designed.

o-See whether your PCR product has the correct size.

p-Cut both the cut vector backbone and your cut PCR product from the gel and purify it using the gel purification kit.

q- Set up a ligation reaction using your cut pMDC32 vector and your cut PCR product.

r-Plate the transformed E. coli cells on LB medium containing kanamycin.

s-Run your restriction digestion on an agarose gel.

t-Prepare total DNA from a wild type Arabidopsis plant.

u-See whether your pMDC32 vector has been completely cut.

v-Use the overnight 3 ml E. coli culture to extract the plasmid product of your ligation reaction using a plasmid miniprep.