

PCR PROTOCOL

Purpose: 'Full length' trnL chloroplast DNA sequencing.

Original Source: Taberlet *et al.* 1991

Document By: TRK (3/6/2020)

Notes: This is the standard PCR protocol to develop a reference library to which dietary P6 sequences may be matched (P6 is a short loop within *trnL*). It is often used for intraspecific population genetic or phylogenetic studies requiring variation in chloroplast DNA haplotypes. Target size varies ~600bp.

Primers:

trnL(UAA)c	CGAAATCGGTAGACGCTACG
trnL(UAA)d	GGGGATAGAGGGACTTGAAC

Reaction (12.5 µL) using standard NEB Taq buffer without MgCl₂.

REAGENTS (Final concentration)	1 X (µL)
10X Buffer without Mg	1.25
25 mM MgCl ₂ (2.5 mM final)	1.25
dNTPs mix -mM EACH (200µM EACH final)	1.00
PRIMER1 – 10uM (0.2µM final)	0.25
PRIMER2 – 10uM (0.2µM final)	0.25
BSA - mg/ml (0.1 mg/ml final)	0.10
DMSO (4% final)	0.05
MOLECULAR WATER	6.28
NEB TAQ	0.06
DNA TEMPLATE	2.0

***Make extra master mix and include a negative control.**

Program: *cp50* on Eppendorf Mastercycler Nexus GSX1

Process	Temperature (°C)	Time
Initial denaturing	95	4 min
35 cycles: Denaturing	94	30 s
35 cycles: Annealing	50	30 s
35 cycles: Extension	72	1:15 min
Final Extension	72	5 min
Hold	10	Hold

Approximate total time: 2.0 hours.