METABOLISME GLUKOSA, UREA, DAN TRIGLISERIDA (TEKNIK SPEKTROFOTOMETRI)

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Pemeriksaan Kadar Gula Darah

	BLANKO	STANDARD	SAMPLE
reagensia kit	1000 µl	1000 µl	1000 µl
standard	-	10 µl	-
Sample	-	-	10 µl
Mix, incubate at 37°C, 10 min, Measure absorbance with λ 500 nm			

Pemeriksaan Kadar Total Protein Plasma

	BLANKO	STANDARD	SAMPLE
reagensia kit	1000 µl	1000 µl	1000 µl
standard	-	10 µl	-
Sample	-	-	10 µl
Mix, incubate at 37°C, 10 min, Measure absorbance with λ 530 nm			

Pemeriksaan Kadar Urea Plasma

	BLANKO	STANDARD	SAMPLE
reagensia A	1000 µl	1000 µl	1000 µl
reagensia B	1000 µl	1000 µl	1000 µl
Mix, incubate at 25°C, 5 min, Measure absorbance with λ 530 nm			
standard	-	10 µl	-
Sample	-	-	10 µl
Mix, incubate at 25°C, 5 min, Measure absorbance with λ 600 nm			

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Glucose GOD FS*

Diagnostic reagent for quantitative in vitro determination of glucose in serum or plasma on photometric systems

Order Information

Cat. No.	Kit size
1 2500 99 10 021	R 5 x 25 mL + 1 x 3 mL Standard
1 2500 99 10 026	R 6 x 100 mL
1 2500 99 10 023	R 1 x 1000 mL
1 2500 99 10 704	R 8 x 50 mL
1 2500 99 10 717	R 6 x 100 mL
1 2500 99 10 917	R 10 x 60 mL
1 2500 99 10 192	R 4 x 60 mL
1 2500 99 10 952	6150 Tests on ADVIA 1650/1800
1 2500 99 10 030	6 x 3 mL Standard

Summary [1,2]

Measurement of glucose concentration in serum or plasma is mainly used in diagnosis and monitoring of treatment in diabetes mellitus. Other applications are the detection of neonatal hypoglycemia, the exclusion of pancreatic islet cell carcinoma as well as the evaluation of carbohydrate metabolism in various diseases.

Method

"GOD-PAP": enzymatic photometric test

Principle

Determination of glucose after enzymatic oxidation by glucose oxidase. The colorimetric indicator is quinoneimine, which is generated from 4-aminoantipyrine and phenol by hydrogen peroxide under the catalytic action of peroxidase (Trinder's reaction) [3].

Reagents

Components and Concentrations

Phosphate buffer	pH 7.5	250 mmol/L
Phenol		5 mmol/L
4-Aminoantipyrine		0.5 mmol/L
Glucose oxidase (GOD)		≥ 10 kU/L
Peroxidase (POD)		≥ 1 kU/L

Standard: 100 mg/dL (5.55 mmol/L)

Storage Instructions and Reagent Stability

The reagent is stable up to the end of the indicated month of expiry, if stored at 2 - 8 °C, protected from light and contamination is avoided. Do not freeze the reagents! The standard is stable up to the end of the indicated month of expiry, if stored at 2 - 25 °C.

Note: It has to be mentioned, that the measurement is not influenced by occasionally occurring color changes, as long as the absorbance of the reagent is < 0.3 at 546 nm.

Warnings and Precautions

- The reagent contains sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
- Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents.

Waste Management

Please refer to local legal requirements.

Reagent Preparation

Reagent and standard are ready to use.

Materials required but not provided

NaCl solution 9 g/L
General laboratory equipment

Specimen

Serum, heparin plasma or EDTA plasma
Separate at the latest 1h after blood collection from cellular contents.

Stability in plasma after addition of a glycolytic inhibitor (Fluoride, moniodacetate, mannose) [4]:

2 days at 20 - 25 °C

7 days at 4 - 8 °C

1 day at -20 °C

Stability in serum (separated from cellular contents, hemolysis free) without adding a glycolytic inhibitor [2,5]:

8 h at 25 °C

72 h at 4 °C

Discard contaminated specimens!

Assay Procedure

Application sheets for automated systems are available on request.

Wavelength	500 nm, Hg 546 nm
Optical path	1 cm
Temperature	20 - 25 °C/37 °C
Measurement	Against reagent blank

	Blank	Sample or standard
Sample or standard	-	10 µL
Dist. water	10 µL	-
Reagent	1000 µL	1000 µL
Mix, incubate 20 min. at 20 - 25 °C or 10 min. at 37 °C. Read absorbance against the blank within 60 min.		

Calculation

With standard or calibrator

$$\text{Glucose [mg/dL]} = \frac{\Delta A_{\text{Sample}}}{\Delta A_{\text{Std/Cal}}} \times \text{Conc. Std/Cal [mg/dL]}$$

Conversion factor

$$\text{Glucose [mg/dL]} \times 0.05551 = \text{Glucose [mmol/L]}$$

TRIGLISERIDA

N.S. BIO-TEC

TRIGLYCERIDES (GPO/PAP)

Enzymatic Colorimetric Determination of Serum Triglycerides

Ref. 5 X 30 ml

INTENDED USE

NS Biotec triglycerides reagent is intended for the in vitro quantitative determination of triglycerides in serum and plasma on both automated and manual systems.


CLINICAL SIGNIFICANCE

Triglycerides are esters of the trihydric alcohol glycerol with 3 long chain fatty acids. They are the main lipids present in human plasma; the others are cholesterol, phospholipids, and non-esterified fatty acids. Triglycerides are synthesized in the intestinal mucosa by the esterification of glycerol and free fatty acids. They are then released into the mesenteric lymphatics and distributed to most tissues for storage. Triglycerides are the main storage lipids in humans, where they constitute about 95% of adipose tissue lipids. Elevated levels of triglycerides have been associated with high risk in severe atherosclerosis. High triglycerides levels and hyperlipidemia in general can be an inherited trait or can be secondary to disorders including diabetes mellitus, nephrosis, biliary obstruction, and metabolic disorders associated with endocrine disturbances¹⁻³.

ASSAY PRINCIPLE

Triglycerides are generally determined by a combination of hydrolysis to glycerol and free fatty acids and measurement of the amount of glycerol released. The most commonly used methods involve alkaline hydrolysis and either chemical or enzymatic measurement of glycerol. Chemical means of analysis generally rely on measurement of the product of periodate oxidation of glycerol. Eggstein and Kreutz developed an enzymatic method for measuring glycerol released from triglycerides by alkaline hydrolysis⁴. This method was based on the coupled reaction sequence catalyzed by glycerol kinase, pyruvate kinase, and lactate dehydrogenase. A method for complete enzymatic hydrolysis to triglycerides avoiding the need for serum pretreatment was described by Bucolo and David, using a combination of lipase and at least one proteolytic enzyme⁵. Wahlefeld reported that certain esterases could be combined with a lipase to achieve complete triglycerides hydrolysis⁶. Both methods employed a coupled enzymatic reaction sequence⁷ to measure glycerol. NS Biotec triglycerides reagent combines the use of lipoproteinlipase, glycerol kinase, and glycerol phosphate oxidase with the peroxidase/4-chlorophenol/4-aminoantipyrine system of Trinder⁸ for the measurement of triglycerides in human serum. The series of reactions involved in the assay system are as follows:

1. Triglycerides are hydrolyzed by lipoprotein lipase (LPL) to glycerol and fatty acids.
2. Glycerol is then phosphorylated to glycerol-3-phosphate by ATP in a reaction catalyzed by glycerol kinase (GK).
3. The oxidation of glycerol-3-phosphate is catalyzed by glycerol phosphate oxidase (GPO) to form dihydroxyacetone phosphate and hydrogen peroxide (H_2O_2).
4. In presence of peroxidase (POD), the hydrogen peroxide (H_2O_2) formed effects the oxidative coupling of 4 - chlorophenol and 4- aminoantipyrine (4-AAP) to form a red-colored quinoneimine dye.

The intensity of the color produced is directly proportional to triglycerides concentration. It is determined by measuring the increase in absorbance at 500–550 nm.

EXPECTED VALUES

Males:	40 – 180 mg/dl (0.45 – 1.82 mmol/l)
Females:	35 – 135 mg/dl (0.4 – 1.54 mmol/l)

For the recognition of the risk factor hyper-triglyceridemia. The

following limits are recommended:

Suspicious	>150 mg/dl (1.71 mmol/l)
Elevated	>200 mg/dl (2.28 mmol/l)

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference range. For diagnostic purposes, the triglycerides results should always be assessed in conjunction with the patient's medical history, clinical examination, and other findings.

REAGENTS

R ₁	Triglycerides standard	200 mg/dl
	Pipes buffer, pH 7.8	50 mmol/l
	p-Chlorophenol	2.0 mmol/l
R ₂	Lipoprotein lipase	1500 U/l
	Glycerol kinase	800 U/l
	Glycerol phosphate oxidase	4000 U/l
	Peroxidase	440 U/l
	4-Aminoantipyrine	0.4 mmol/l
	ATP	0.3 mmol/l
	Mg ²⁺	40 mmol/l
	Sodium cholate	0.2 mmol/l

Reagent Preparation & Stability

All reagents are ready for use and stable up to the expiry date given on label when stored at 2–8°C.

SPECIMEN

Serum or plasma* from fasting patients.

- The only accepted anticoagulants are heparin and EDTA.

Specimen Preparation & Stability

Patients should refrain from eating for 10 to 14 hours before blood is drawn. Samples must be drawn in a soap and glycerol free collection device.

Blood should be collected by venipuncture, after the patient has been in a seated position for at least 5 minutes. Tourniquet usage should be kept to a minimum and the specimen should be allowed to clot for 30 minutes at room temperature⁸.

The best specimen is unhemolysed serum, and should be analyzed on the day of collection. Specimens are stable for 7 days when stored at 4°C; several months at -20°C and for years at -70°C¹.

PROCEDURE

- **Manual Procedure**

Wavelength	500 - 550 nm
Cuvette	1 cm light path
Temperature	20-25 or 37 °C
Zero adjustment	against reagent blank
Specimen	Serum or plasma

	Blank	Standard	Specimen
R ₂	1.0 ml	1.0 ml	1.0 ml
Standard	—	10 μ l	—
Specimen	—	—	10 μ l

Mix, incubate for 5 minutes at 37°C or 10 minutes at 20-25°C. Measure the absorbance of specimen (A_{specimen}) and standard (A_{standard}) against reagent blank.

The color is stable for 60 minutes.

- Automated Procedure
User defined parameters for different auto analyzers are available upon request.
- CALCULATION
Calculate the triglycerides concentration by using the following formulae:
$$\text{Triglycerides Concentration} = \frac{\text{Absorbance of Specimen}}{\text{Absorbance of Standard}} \times \text{Standard}$$
- Unit conversion
$$\text{mg/dl} \times 0.0114 = \text{mmol/l}$$

LINEARITY
When run as recommended, the assay is linear up to 900 mg/dl (20.7 mmol/l).
If result exceeds 900 mg/dl (10.26 mmol/l), specimen should be diluted with 0.9% NaCl solution and reassayed. Multiply the result by the dilution factor.

SENSITIVITY
The sensitivity is defined as the change of analytical response per unit change in analyte concentration at a path length of 1 cm.

When run as recommended the sensitivity of this assay is 3.0 mg/dl (0.034 mmol/l).

QUALITY CONTROL
It is recommended that controls (normal and abnormal) be included in:

- Each set of assays, or
- At least once a shift, or
- When a new bottle of reagent is used, or
- After preventive maintenance is performed or a clinical component is replaced.

Commercially available control material with established triglycerides values may be routinely used for quality control.
Failure to obtain the proper range of values in the assay of control material may indicate:

- Reagent deterioration,
- Instrument malfunction, or
- Procedure errors.

The following corrective actions are recommended in such situations:

- Repeat the same controls.
- If repeated control results are outside the limits, prepare fresh control serum and repeat the test.
- If results on fresh control material still remain outside the limits, then repeat the test with fresh reagent.
- If results are still out of control, contact NS Biotec Technical Services.

- INTERFERING SUBSTANCES**
- Anticoagulants:
The only acceptable anticoagulants are heparin and EDTA.
 - Bilirubin:
No interference from free bilirubin up to level of 10 mg/dl and from conjugated bilirubin up to a level of 12 mg/dl.
 - Drugs:
Methyldopa and noramidopyrine causes artificially low triglycerides values at the tested drug level. For a more comprehensive review of drugs affecting triglycerides assays refer to the publication by Young⁹.
 - Haemoglobin:
No interference from haemoglobin up to a level of 600 mg/dl.
 - Lipemia:
Extremely lipemic samples can produce a normal triglycerides result (triglycerides grater than 3000 mg/dl).

- Others:
Ascorbic acid levels higher than 2.0 mg/dl decrease the apparent triglycerides concentration significantly.
- WARNING & PRECAUTIONS**
- NS Biotec triglycerides reagent is for in vitro diagnostic use only. Normal precautions exercised in handling laboratory reagents should be followed.
 - Warm up working solution to the corresponding temperature before use.
 - The reagent and sample volumes may be altered proportionally to accommodate different spectrophotometer requirements.
 - Valid results depend on an accurately calibrated instrument, timing, and temperature control.
 - The reagent blank will not exceed an absorbance of 0.06 but don't use the reagent if it is turbid or if the absorbance is greater than 0.2 at 500 nm.
 - Extremely lipemic specimens can produce a normal result. Dilute specimens' 1+4 with saline and reassayed. Multiply the result by 5.

BIBLIOGRAPHY

- Stein, EA (1987): Lipids, lipoproteins, and apolipoproteins. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders 448-481.
- Naito, HK (1984): Disorders of lipid metabolism. In: Kaplan LA, Pesce AJ, eds. Clinical Chemistry, theory, analysis, and correlation. St. Louis: Mosby Company. 550-593.
- Fredrickson, DS, Levy, RI and Lees, RS (1967): Fat transport in lipoproteins an integrated approach to mechanisms and disorders. N Engl J Med. 276: 34-43.
- Eggstein, M & Kreutz,
- F (1966): Klin. Wschr. 44: 262-267.

	Consult Instruction for Use
	Caution Consult Accompanying Documents
	In Vitro Diagnostic Medical Device
	Temperature Limitation
	Manufacturer
	Authorized Representative In The European Community
	Catalogue Number
	Batch Code
	Use By


NS BIOTEC
MEDICAL EQUIPMENTS
7 Sad Zaghloul sq. Mehatef Elraml,
Alexandria – Egypt
Tele: 002 03 486 5140
Fax: 002 03 485 7597
Website : www.nsbiotec.com
E- mail : info@nsbiotec.com

Wellkang Ltd
sulte B , 29 Harley
street LONDON ,
W 1 G 9 QR , U.K.
www.CE-marking .EU

UREA

Principle

Reagents

Kit 5 x 100 ml (Ref.99 36 48). Contents:
A. 5 x 100 ml. Urease / Salicylate. (Ref. 99 21 04)
B. 1 x 15 ml Alkaline hypochlorite. (Ref. 99 14 75)
C. 1 x 5 ml Standard. (Ref. 99 02 41)
Aqueous solution of Urea equivalent to 40 mg / dl. (6.6 mmol/L). Ready-to-use.

Working reagent

A. Dissolve the contents of the Urease/Salicylate vial with the volume of deionized water stated on the label.

The concentrations in the working reagent are:

Phosphate buffer pH 6.8	20 mM
Sodium salicylate	61 mM
Sodium nitroprusiate	3.4 mM
EDTA-Na ₂	1.34 mM
Urease	≥ 23 U/ml
Stabilizers	

B. Dilute the contents of the Alkaline hypochlorite vial up to 500 ml. of deionized water.

Concentrations of reagent solution are:

Alkaline hypochlorite	7.5 mM
NaOH	160 mM

Storage and stability

The components of the kit, stored at 2-8°C, will remain stable until the expiration date stated on the label. Once the Urease/Salicylate vial has been dissolved, will remain stable for 3 weeks, if stored in amber bottle at 2-8° C. The Alkaline hypochlorite solution will remain stable for 8 months, if stored in the same way.

Sample

Serum, plasma and urine. Urea will remain stable in serum for at least 1 day at room temperature (≤ 25°C), 5 days at 2-8°C and 6 months when frozen (-20° C). In urine, urea will remain stable, when kept at 2-8°C, for 5 days, provided that the pH value be lower than 4.
If a urine sample is to be assayed, it should be previously diluted 1/100 with deionized water.
Multiply the final result by 100.

Caution

Reagent B: In case of contact with the skin, mucose or eyes, wash thoroughly with water and ask the physician.
The reagent A contains Sodium azide at 0.09%. Handle with care. The disposal of the residues has to be made according to legal local regulations.

Procedure	BL ml	SA ml	ST ml
Standard	—	—	0.01
Sample	—	0.01	—
Reagent A	1.00	1.00	1.00

Mix and incubate 3 min. at 37°C. or 5 min. at room temperature (≤ 25°C).

Reagent B	1.00	1.00	1.00
-----------	------	------	------

Mix and incubate again 3 min. at 37°C. or 5 min. at room temperature (≤ 25°C).

Reading
Wavelength: Hg 578 nm; 600 nm.
Blank: The contents of the BL tube.
Colour stability: 4 hours.

Calculations

SA O.D. _____
ST O.D. _____

_____ x 40 = mg urea / dl

S.I. Units
(mg/dl) x 0.1665 = mmol/L

Normal values
Serum, plasma: 15 - 45 mg/dl.
Urine: 20 - 35 g / L.

Performance Characteristics

Linearity: Up to 400 mg/dl of Urea. For higher values, dilute the sample 1/2 in deionized water and assay once again. Multiply the final result by 2.
The analytical performance characteristics of the product depend both of the reagent and the reading system used, manual or automatic. The following data have been obtained manually.

Intraseres Variation Coefficient: 1.66%
Interseries Variation Coefficient: 2.05%
Recovery: 97.9 %.

Any glassware contamination by ammonium salts or ammonia should be avoided. Serum samples should be free from hemolysis and turbidity. Fluoride as well as ammonium heparinate inhibit the reaction.

Quality control

Seriscann Normal (Normal Control Serum)
(Ref. 99 41 48) and Seriscann Anormal (Abnormal Control Serum) (Ref. 99 46 85).

References

Foster, L.B., Hochholzer, J.M. (1971), Clin. Chem., 17, 921-925.
Wilcox, A., Wallace, E.C., Sterling, R.E., David, H.A., Ware, A.G. (1966), Clin. Chem. 12, 151-157