

Pichia Expression Kit

For Expression of Recombinant Proteins in *Pichia pastoris*

Catalog no. K1710-01

Revision date: 21 December 2009

Manual part no. 25-0043

MAN0000012

User Manual

Table of Contents

	Kit Contents and Storage	iv
	Required Equipment and Supplies	V
n	ntroduction	1
	Pichia pastoris Expression System	
	Experimental Outline	
	•	
VI	lethods	
	Pichia Strains	
	E. coli Strains	
	Selecting a <i>Pichia</i> Expression Vector	
	pHIL-D2	
	pPIC3.5	
	pHIL-S1	
	pPIC9	
	Signal Sequence Processing	
	Cloning into the <i>Pichia</i> Expression Vectors	
	Preparing Transforming DNA	
	Growing Pichia for Spheroplasting	
	Preparing Spheroplasts	
	Transforming <i>Pichia</i>	
	Screening for Mut ⁺ and Mut ^S Transformants	
	PCR Analysis of <i>Pichia</i> Integrants	
	Expression of Recombinant <i>Pichia</i> Strains	
	Analyzing Samples by SDS-Polyacrylamide Gel Electrophoresis	
	Optimizing <i>Pichia</i> Protein Expression	
	Scaling Up Expression	
	Protein Purification and Glycosylation	
^	ppendix	50
_	• •	
	E. coli Media Recipes	
	Pichia Media Recipes	
	Recombination and Integration in <i>Pichia</i>	
	Electroporation of <i>Pichia</i>	
	PEG 1000 Transformation Method for <i>Pichia</i>	
	Lithium Chloride Transformation Method	
	Direct PCR Screening of <i>Pichia</i> Clones	
	Isolating Total DNA from <i>Pichia</i>	
	Detecting Multiple Integration Events	
	Procedure for Total RNA Isolation from <i>Pichia</i>	
	β-Galactosidase Assay	
	Accessory Products	
	Technical Support	
	Purchaser Notification	
	References	94

Kit Contents and Storage

Kit Contents

The *Pichia* **Expression Kit** is shipped at room temperature and contains the following components.

Spheroplast Module (Box 1). Store at room temperature.

Reagent	Amount	Components	
SOS medium	20 mL	1 M Sorbitol	
		0.3X YPD	
		10 mM CaCl ₂	
Sterile Water	$2 \times 125 \text{ mL}$	Autoclaved, deionized water	
SE	2 × 125 mL	1 M Sorbitol	
		25 mM EDTA, pH 8.0	
SCE	2 × 125 mL	1 M Sorbitol	
		10 mM Sodium citrate buffer, pH 5.8	
		1 mM EDTA	
1 M Sorbitol	2 × 125 mL		
CaS	$2 \times 60 \text{ mL}$	1 M Sorbitol	
		10 mM Tris-HCl, pH 7.5;	
		10 mM CaCl ₂	
40% PEG	25 mL	40% (w/v) PEG 3350 (Reagent grade) in water	
СаТ	25 mL	20 mM Tris-HCl, pH 7.5	
		20 mM CaCl ₂	

Spheroplast Module (Box 2). Store at -20°C.

Reagent	Amount	Components
Zymolyase	$10 \times 20 \ \mu L$	3 mg/mL Zymolyase in water
		(100,000 units/g lytic activity)
1 M DTT	$10 \times 1 \text{ mL}$	1 M dithiothreitol in water

Stab Vials: Pichia and E. coli stabs. Store at 4°C.

Strain	Amount	Genotype	Phenotype (<i>Pichia</i> only)
GS115	1 stab	his4	Mut ⁺
KM71	1 stab	arg4 his4 aox1::ARG4	Mut ^S , Arg ⁺
GS115 Albumin	1 stab	HIS4	Mut ^s
GS115 β-Gal	1 stab	HIS4	Mut ⁺
TOP10F′	1 stab	F' {proAB, lacI ^q , lacZΔM15, Tn10 (Tet ^R)} mcrA, Δ(mrr-hsdRMS-mcrBC), φ80lacZΔM15, ΔlacX74, deoR, recA1, λ ⁻ araD139, Δ(ara-leu)7697, galU, galK, rpsL(Str ^R), endA1, nupG	

Kit Contents and Storage, continued

Kit Contents, continued

Vectors. Store at −20°C.

Reagent	Description
pHIL-D2	Vector for intracellular expression in <i>Pichia</i> .
10 μg, 20 μL at 0.5 μg/μL in TE buffer, pH 8.0*	
pPIC3.5	Vector for intracellular expression in <i>Pichia</i> .
10 μg, 20 μL at 0.5 μg/μL in TE buffer, pH 8.0	
pHIL-S1	Vector for secreted expression in Pichia.
10 μg, 20 μL at 0.5 μg/μL in TE buffer, pH 8.0	Uses the <i>PHO1</i> signal sequence.
pPIC9	Vector for secreted expression in Pichia.
10 μg, 20 μL at 0.5 μg/μL in TE buffer, pH 8.0	Uses the α -factor signal sequence.

^{*}TE buffer, pH 8.0: 10 mM Tris-HCl, 1 mM EDTA, pH 8.0

Primers. Store at -20°C.

5' AOX1 sequencing primer	5'-GACTGGTTCCAATTGACAAGC-3'
2 μg (312 pmoles), lyophilized	
3´ AOX1 sequencing primer	5′-GCAAATGGCATTCTGACATCC-3′
2 μg (314 pmoles), lyophilized	
α-Factor sequencing primer	5′-TACTATTGCCAGCATTGCTGC-3′
2 μg (315 pmoles), lyophilized	

Media

The following prepackaged media is included for your convenience. Instructions for use are provided on the package. Store at room temperature.

Media	Amount	Yield
YP Base Medium	2 pouches	2 liters of YP medium
YP Base Agar Medium	2 pouches	2 liters of YP medium
Yeast Nitrogen Base	1 pouch	500 mL of 10X YNB



The *Pichia* Spheroplast Module for transforming *Pichia* by spheroplasting is available separately from Invitrogen (see page 88 for ordering information).

Intended Use

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

Required Equipment and Supplies

Required Equipment and Supplies (not provided)

- 30°C rotary shaking incubator
- Water baths capable of 37°C, 45°C, and 100°C
- Centrifuge suitable for 50 mL conical tubes (floor or table-top)
- Baffled culture flasks with metal covers (50 mL, 250 mL, 500 mL, 1000 mL, and 3 L)
- 50 mL sterile, conical tubes
- 6 mL and 15 mL sterile snap-top tubes (Falcon 2059 or similar)
- UV Spectrophotometer
- Mini agarose gel apparatus and buffers
- Agarose and low-melt agarose
- Polyacrylamide gel electrophoresis apparatus and buffers
- Media for transformation, growth, screening, and expression (see Recipes, pages 59–66)
- 5% SDS solution (10 mL per transformation)
- Sterile cheesecloth or gauze
- Breaking Buffer (see **Recipes**, page 66)
- Acid-washed glass beads (available from Sigma)
- Replica-plating equipment (optional)
- Bead Beater[™] (optional, available from Biospec)

Introduction

Pichia pastoris Expression System

Review Articles

The information presented here is designed to give you a concise overview of the *Pichia pastoris* expression system. It is by no means exhaustive. For further information, read the articles cited in the text along with the following review articles (Buckholz & Gleeson, 1991; Cregg & Higgins, 1995; Cregg *et al.*, 1993; Nico-Farber *et al.*, 1995; Romanos, 1995; Sreekrishna *et al.*, 1988; Wegner, 1990). A general review of foreign gene expression in yeast is also available (Romanos *et al.*, 1992).

General Characteristics of *Pichia* pastoris

As a eukaryote, *Pichia pastoris* has many of the advantages of higher eukaryotic expression systems such as protein processing, protein folding, and posttranslational modification, while being as easy to manipulate as *E. coli* or *Saccharomyces cerevisiae*. It is faster, easier, and less expensive to use than other eukaryotic expression systems such as baculovirus or mammalian tissue culture, and generally gives higher expression levels. As a yeast, it shares the advantages of molecular and genetic manipulations with *Saccharomyces*, and it has the added advantage of 10- to 100-fold higher heterologous protein expression levels. These features make *Pichia* very useful as a protein expression system.

Similarity to Saccharomyces

Many of the techniques developed for *Saccharomyces* may be applied to *Pichia*. These include:

- Transformation by complementation
- Gene disruption
- Gene replacement

In addition, the genetic nomenclature used for *Saccharomyces* has been applied to *Pichia*. For example, the *HIS4* gene in both Saccharomyces and *Pichia* encodes histidinol dehydrogenase. There is also cross-complementation between gene products in both *Saccharomyces* and *Pichia*. Several wild-type genes from *Saccharomyces* complement comparable mutant genes in *Pichia*. Genes such as *HIS4*, *LEU2*, *ARG4*, *TRP1*, and *URA3* all complement their respective mutant genes in *Pichia*.

Pichia pastoris as a Methylotrophic Yeast

Pichia pastoris is a methylotrophic yeast, capable of metabolizing methanol as its sole carbon source. The first step in the metabolism of methanol is the oxidation of methanol to formaldehyde using molecular oxygen by the enzyme alcohol oxidase. In addition to formaldehyde, this reaction generates hydrogen peroxide. To avoid hydrogen peroxide toxicity, methanol metabolism takes place within a specialized cell organelle, called the peroxisome, which sequesters toxic by-products away from the rest of the cell. Alcohol oxidase has a poor affinity for O₂, and *Pichia pastoris* compensates by generating large amounts of the enzyme. The promoter regulating the production of alcohol oxidase is the one used to drive heterologous protein expression in *Pichia*.

Pichia pastoris Expression System, continued

Two Alcohol Oxidase Proteins

Two genes in *Pichia pastoris* code for alcohol oxidase–AOX1 and AOX2. The AOX1 gene product accounts for the majority of alcohol oxidase activity in the cell. Expression of the AOX1 gene is tightly regulated and induced by methanol to very high levels, typically $\geq 30\%$ of the total soluble protein in cells grown on methanol. The AOX1 gene has been isolated and a plasmid-borne version of the AOX1 promoter is used to drive expression of the gene of interest encoding the desired heterologous protein (Ellis *et al.*, 1985; Koutz *et al.*, 1989; Tschopp *et al.*, 1987a). While AOX2 is about 97% homologous to AOX1, growth on methanol is much slower than with AOX1. This slow growth on methanol allows isolation of Mut^S strains (aox1) (Cregg *et al.*, 1989; Koutz *et al.*, 1989).

Expression

Expression of the *AOX1* gene is controlled at the level of transcription. In methanol-grown cells approximately 5% of the polyA+ RNA is from the *AOX1* gene. The regulation of the *AOX1* gene is a two step process: a repression/derepression mechanism plus an induction mechanism (*e.g.*, *GAL1* gene in *Saccharomyces* (Johnston, 1987)). Briefly, growth on glucose represses transcription, even in the presence of the inducer methanol. For this reason, growth on glycerol is recommended for optimal induction with methanol. Note that growth on glycerol only (derepression) is not sufficient to generate even minute levels of expression from the *AOX1* gene. The inducer, methanol, is necessary for even detectable levels of *AOX1* expression (Ellis *et al.*, 1985; Koutz *et al.*, 1989; Tschopp *et al.*, 1987a).

Phenotype of aox1 mutants

Loss of the *AOX1* gene, and thus a loss of most of the cell's alcohol oxidase activity, results in a strain that is phenotypically Mut^S (Methanol utilization slow). This has in the past been referred to as Mut⁻. The Mut^S designation has been chosen to accurately describe the phenotype of these mutants. This results in a reduction in the cells' ability to metabolize methanol. The cells, therefore, exhibit poor growth on methanol medium. Mut⁺ (Methanol utilization plus) refers to the wild type ability of strains to metabolize methanol as the sole carbon source. These two phenotypes are used when evaluating *Pichia* transformants for integration of your gene (Experimental Outline, page 4).

Intracellular and Secretory Protein Expression

Heterologous expression in *Pichia pastoris* can be intracellular or secreted. Secretion requires the presence of a signal sequence on the expressed protein to target it to the secretory pathway. While several different secretion signal sequences have been used successfully, including the native secretion signal present on some heterologous proteins, success has been variable. The secretion signal sequence from the *Saccharomyces cerevisiae* factor prepro peptide has been used with the most success (Cregg *et al.*, 1993; Scorer *et al.*, 1993).

The major advantage of expressing heterologous proteins as secreted proteins is that *Pichia pastoris* secretes very low levels of native proteins. Since there is very low amount of protein in the minimal *Pichia* growth medium, this means that the secreted heterologous protein comprises the vast majority of the total protein in the medium and serves as the first step in purification of the protein (Barr *et al.*, 1992) . However, that if there are recognized glycosylation sites (Asn-X-Ser/Thr) in your protein's primary sequence, glycosylation may occur at these sites.

Pichia pastoris Expression System, continued

Posttranslational Modifications

In comparison to *Saccharomyces cerevisiae*, *Pichia* may have an advantage in the glycosylation of secreted proteins because it may not hyperglycosylate. Both *Saccharomyces cerevisiae* and *Pichia pastoris* have a majority of N-linked glycosylation of the high-mannose type; however, the length of the oligosaccharide chains added posttranslationally to proteins in *Pichia* (average 8–14 mannose residues per side chain) is much shorter than those in *Saccharomyces cerevisiae* (50–150 mannose residues) (Grinna and Tschopp, 1989; Tschopp *et al.*, 1987b). Very little O-linked glycosylation has been observed in *Pichia*.

In addition, *Saccharomyces cerevisiae* core oligosaccharides have terminal α 1,3 glycan linkages whereas *Pichia pastoris* does not. It is believed that the α 1,3 glycan linkages in glycosylated proteins produced from *Saccharomyces cerevisiae* are primarily responsible for the hyper-antigenic nature of these proteins making them particularly unsuitable for therapeutic use. Although not yet proven, this is predicted to be less of a problem for glycoproteins generated in *Pichia pastoris*, because it may resemble the glycoprotein structure of higher eukaryotes (Cregg *et al.*, 1993).

Experimental Outline

Selection of Vector and Cloning

To utilize the strong, highly inducible P_{AOXI} promoter for expressing your protein, four expression vectors are included in this kit. pHIL-D2 and pPIC3.5 are used for intracellular expression, and pHIL-S1 and pPIC9 are used for secreted expression (see pages 14–17 for more information). Before cloning your insert, you must:

- decide whether you want intracellular or secreted expression.
- analyze your insert for the following restriction sites: Sac I, Stu I, Sal I, Not I, and Bgl II. We recommend these sites for linearizing your construct prior to Pichia transformation. If your insert has all of these sites, refer to pages 29–30 for alternate sites.

Transformation and Integration

Two different phenotypic classes of His⁺ recombinant strains can be generated: Mut⁺ and Mut^S. Mut^S refers to the "Methanol utilization slow" phenotype caused by the loss of alcohol oxidase activity encoded by the *AOX1* gene. A strain with a Mut^S phenotype has a mutant *aox1* locus, but is wild type for *AOX2*. This results in a slow growth phenotype on methanol medium. Transformation of strain GS115 can yield both classes of transformants, His⁺ Mut⁺ and His⁺ Mut^S, while KM71 yields only His⁺ Mut^S, because the strain itself is Mut^S. Both Mut⁺ and Mut^S recombinants are useful to have, because one phenotype may favor better expression of your protein than the other. Because of clonal variation, you should test 6–10 recombinants per phenotype. There is no way to predict beforehand which construct or isolate will better express your protein. We strongly recommend that you analyze *Pichia* recombinants by PCR to confirm the integration of your construct (see page 43).

After you have successfully cloned your gene, you will linearize your plasmid to stimulate recombination when the plasmid is transformed into *Pichia*. The table below describes the types of recombinants you will get by selective digestion of your plasmid.

Restriction Enzyme	Integration Event	GS115 Phenotype	KM71 Phenotype
Sal I or Stu I	Insertion at his4	His+ Mut+	His+ Mut ^S
Sac I	Insertion at 5′ AOX1 region	His+ Mut+	His+ Mut ^S
Not I or Bgl II	Replacement at AOX1 locus	His+ Mut ^S His+ Mut+	His+ Mut ^S (not recommended, see page 7)

Experimental Outline, continued

Expression and Scale-up

After confirming your *Pichia* recombinants by PCR, you will test expression of both His⁺ Mut⁺ and His⁺ Mut^S recombinants. This procedure involves growing a small culture of each recombinant, inducing them with methanol, and taking time points. If looking for intracellular expression, analyze the cell pellet from each time point by SDS polyacrylamide gel electrophoresis (SDS-PAGE). If looking for secreted expression, analyze both the cell pellet and supernatant from each time point. We recommend that you analyze your SDS-PAGE gels by Coomassie staining and, if you have an antibody to your protein, by western blot. We also suggest checking for protein activity by an activity assay, if one is available. Not all proteins express to the level of grams per liter, so it is advisable to check by western blot or activity assay, and not just by Coomassie staining of SDS-PAGE gels for production of your protein.

Choose the *Pichia* recombinant strain that best expresses your protein and optimize induction based on the suggestions on pages 52–53. After you optimize expression, scale-up your expression protocol to produce more protein.

Experimental Process

The overall experimental process is divided into two major sections: **Generating Recombinant Strain** and **Induction (Mut**⁺ **and/or Mut**⁵). Each section contains a table outlining the major steps of the experimental process. Each step is discussed in detail further in the manual. Refer to the indicated pages to read about particular steps of interest. The discussion about recombination and integration in *Pichia* will help you choose the right vector. For more information, refer to the review by Higgins (Higgins, 1995).

Generating Recombinant Strain

The goal of this section is to create a *Pichia pastoris* strain containing your integrated gene of interest. Before starting your experiments, determine which vector to use.

Step	Procedure	Page
1	Select the appropriate expression vector (For more information, refer to Recombination and Integration in <i>Pichia</i> , pages 69–72)	11–17
2	Clone your gene of interest into selected vector	19–24
3	Transform <i>E. coli</i> , select ampicillin-resistant transformants, and confirm the presence and orientation of your gene of interest	25
4	Linearize the constructs with appropriate restriction enzymes to generate His ⁺ Mut ^s and His ⁺ Mut ⁺ recombinant strains	27–30
5	Transform and select His ⁺ transformants (GS115 recombinants, His ⁺ Mut ⁺ ; KM71 recombinants, His ⁺ Mut ^S)	31–37
6	Screen His ⁺ transformants for Mut ⁺ and Mut ^S strains (6–10 recombinants of each phenotype)	38–42
7	Confirm the integration of your gene of interest in Mut ⁺ and Mut ^S recombinants by PCR	43–44

Experimental Outline, continued

Mut⁺ Induction

The method of induction depends on whether the recombinant is Mut⁺ or Mut^S. The differences primarily occur in the culture volumes and the time of induction (see below). Refer to the following pages for more detailed instructions.

Step	Procedure	Page
1	Guidelines for expression of recombinant proteins in Pichia	45–46
2	Grow His ⁺ Mut ⁺ recombinants in 25 mL of buffered glycerol medium to a final OD_{600} = 2–6	47
3	Harvest the cells and resuspend them to an OD_{600} of 1.0 (~100–200 mL) with methanol medium. Place the cell suspension in a 1 liter baffled flask	47
4	Incubate the culture at 30°C with shaking and take samples for analysis at 0, 6, 12, 24, 36, 48, 60, 72, 84, and 96 hours	47
5	Analyze the medium (if protein of interest is targeted for secretion) and the cell lysates (for intracellular and secreted expression) for protein via PAGE/Coomassie Blue staining, western blot, activity, ELISA, or immunoprecipitation	49–51
6	Optimize expression of your His+ Mut+ recombinant	52-53
7	Scale-up your expression for protein purification	54–56

Mut^S Induction

This is very similar to Mut⁺ induction except that Mut^S grow very slowly on methanol. To compensate, cells are concentrated to increase cell mass before induction.

Step	Procedure	Page
1	Guidelines for expression of recombinant proteins in Pichia	45–46
2	Grow His ⁺ Mut ^S recombinants in 100–200 mL of buffered glycerol medium to a final $OD_{600} = 2-6$	48
3	Harvest the cells and resuspend them to an OD_{600} of 10.0 (~10–20 mL) with methanol medium. Place the cell suspension in a 100 mL or 250 mL baffled flask.	48
4	Incubate the culture at 30°C with shaking and take samples for analysis at 0, 24, 48, 72, 96, 120, and 144 hours	48
5	Analyze the medium (if protein of interest is targeted for secretion) and the cell lysates (for intracellular and secreted expression) for protein via PAGE/Coomassie Blue staining, western blot, activity, ELISA, or immunoprecipitation	49–51
6	Optimize expression of your His+ Mut+ recombinant	52-53
7	Scale-up your expression for protein purification	54–56

Methods

Pichia Strains

Introduction

Pichia pastoris is quite similar to Saccharomyces cerevisiae as far as general growth conditions and handling. You should be familiar with basic microbiological and sterile techniques before attempting to grow and manipulate any microorganism. You should also be familiar with basic molecular biology and protein chemistry. Some general references to consult are Guide to Yeast Genetics and Molecular Biology (Guthrie & Fink, 1991), Current Protocols in Molecular Biology (Ausubel et al., 1994), Molecular Cloning: A Laboratory Manual (Sambrook et al., 1989), Protein Methods (Bollag et al., 1996), and Guide to Protein Purification (Deutscher, 1990).

Genotype of Pichia Strain

The *Pichia* host strains GS115 and KM71 have a mutation in the histidinol dehydrogenase gene (*his4*) that prevents them from synthesizing histidine. All expression plasmids carry the *HIS4* gene that complements *his4* in the host, so transformants are selected for their ability to grow on histidine-deficient medium. Spontaneous reversion of GS115 and KM71 to His⁺ prototrophy is less than 1 out of 10⁸.

The parent strain of KM71 has a mutation in the argininosuccinate lyase gene (*arg4*) that prevents the strain from growing in the absence of arginine. The wild-type *ARG4* gene was used to disrupt *AOX1*, creating KM71, a Mut^S, Arg⁺, His⁻ strain.

Both GS115 and KM71 will grow on complex medium such as YPD (also known as YEPD) and on minimal media supplemented with histidine. Until transformed, neither GS115 nor KM71 will grow on minimal medium alone as they are His⁻.

Note: Mut^S (\underline{M} ethanol \underline{ut} ilization slow) phenotype has in the past been referred to as Mut⁻. The Mut^S designation has been chosen to accurately describe the phenotype of these mutants.

Construction of KM71

The ARG4 gene (~2 kb) was inserted into the cloned, wild-type AOX1 gene between the BamH I site (codons 15/16 of AOX1) and the Sal I site (codons 227/228 of AOX1). ARG4 replaces codons 16 through 227 of AOX1. This construct was transformed into the parent strain of KM71 (arg4 his4) and Arg^+ transformants were isolated and analyzed for the Mut^S phenotype. Genetic analysis of Arg^+ transformants showed that the wild-type AOX1 gene was replaced by the aox1::ARG4 construct.



The advantage of using KM71 is that there is no need to screen for the Mut phenotype on methanol minimal medium. All transformants will be Mut^S. Secondly, since the *AOX1* locus was not completely deleted, it is theoretically possible to replace *aox1::ARG4* with your construct by gene replacement. The phenotype of this strain would be His⁺ Mut^S Arg⁻. This means the recombinant strain would require arginine in the medium to grow. Unfortunately, simple inclusion of arginine does not totally alleviate the effects of the *arg4* mutation, and *arg4* strains do not grow well on minimal medium supplemented with arginine. Therefore, we do **not** recommend that you generate His⁺ transformants in KM71 by replacing the *aox1::ARG4* construct.

Pichia Strains, continued

Control Expression Strains

GS115/His⁺ **Mut**^S **Albumin**: This strain is a control for secreted expression and the Mut^S phenotype when screening *Pichia* transformants (page 38). The gene for serum albumin was cloned with its native secretion signal, then integrated into *Pichia* at the *AOX1* locus. This strain secretes albumin (67 kDa) into the medium at levels > 1 gram/liter.

GS115/His⁺ **Mut**⁺ **β-galactosidase**: This strain is a control for intracellular expression and the Mut⁺ phenotype when screening *Pichia* transformants (page 38). The gene for β-galactosidase (*lacZ*) was integrated into *Pichia* at the *his4* locus. This strain expresses β-galactosidase (117 kDa) at levels that can be detected on Coomassie-stained SDS-PAGE (see pages 49–51) or assayed using ONPG (see page 86–87).

Growth of *Pichia* Strains

The growth temperature of *Pichia pastoris* is 28–30°C for liquid cultures, plates, and slants. Growth above 32°C during induction can be detrimental to protein expression and can even lead to cell death. Other important facts:

- Doubling time of log phase Mut⁺ or Mut^S Pichia in YPD is ~2 hours
- Mut⁺ and Mut^S strains do not differ in growth rates unless grown on methanol
- Doubling time of log phase Mut⁺ Pichia in methanol medium (MM) is 4–6 hours
- Doubling time of log phase Mut^S Pichia in MM is ~18 hours
- One $OD_{600} = \sim 5 \times 10^7 \text{ cells/mL}$

Note that growth characteristics may vary depending on the recombinant strain.

Growth on Methanol

When plates or medium containing methanol are used as growth medium, it is advisable to add methanol every day to compensate for loss due to evaporation or consumption.

- For plates add 100 μL of 100% methanol to the lid of the inverted plate.
- For liquid medium add 100% methanol to a final concentration of 0.5%.

Some researchers have had success adding methanol to 1% every day for Mut^s strains and up to 3% for Mut⁺ without any negative effect to their liquid culture.

Pichia Strains, continued



Storing *Pichia* Strains

Make frozen stocks for long-term storage of all three *Pichia* strains included in this kit (see below).

To store cells for weeks to months, use YPD medium or YPD agar slants (see page 61).

- 1. Streak for single colonies of the desired strain on YPD.
- 2. Transfer one colony to a YPD stab and grow for 2 days at 30°C.
- 3. You can store the cells on YPD for several weeks at 4°C.

To store cells for months to years, store frozen at -80°C.

- 1. Culture a single colony of the desired strain overnight in YPD.
- 2. Harvest the cells and suspend in YPD containing 15% glycerol at a final OD₆₀₀ of 50–100 (approximately 2.5×10^9 – 5.0×10^9 cells/mL).
- 3. Freeze the cells in liquid nitrogen or a dry ice/ethanol bath, and store at -80° C.



After extended storage at 4°C or –80°C, we recommend checking the His⁺ transformants for correct genotype and viability by streaking on MM, MD or MGY plates before using again.

E. coli Strains

Genotype of E. coli Strain

The *E. coli* strain, TOP10F′ is provided in case no suitable *E. coli* strain is available. Other strains which may be suitable are TOP10, DH5 α F′, JM109, or any other strain which is recombination deficient (recA) and deficient in endonuclease A (endA).

F´ {proAB, lacI¹, lacZΔM15, Tn10 (TetR)} mcrA, $\Delta(mrr-hsdRMS-mcrBC)$, $\phi80lacZ\DeltaM15$, $\Delta lacX74$, recA1, λ^- araD139, $\Delta (ara-leu)7697$, galU, galK, $rpsL(Str^R)$, endA1, nupG

Note: If you do not plan to perform single-stranded DNA rescue, *E. coli* strains that do not carry the F´ episome are also suitable for use.




We recommend that you make a frozen stock of TOP10F' to keep on hand.

- 1. Culture TOP10F′ in 5 mL LB with 10 μg/mL tetracycline. Grow overnight.
- 2. Mix thoroughly 0.85 mL of culture with 0.15 mL sterile glycerol.
- 3. Transfer to a freezer vial and freeze in liquid nitrogen or a dry ice/ethanol bath.
- 4. Store at -80°C.

Selecting a Pichia Expression Vector

Generic Structure

All the vectors included in this kit share several general features shown in black, while some of the vectors also have signal sequences (Sig) and/or an f1 bacteriophage origin. For details of each individual plasmid refer to pages 14–17.





There is no yeast origin of replication in any of the *Pichia* expression vectors included in this kit. His⁺ transformants can only be isolated if recombination occurs between the plasmid and the *Pichia* genome.

Selecting a Pichia Expression Vector, continued

Features

The table below describes the general and optional features of the ${\it Pichia}$ expression vectors.

Feature	Description	Benefit
5' AOX1	An ~1000 bp fragment containing the <i>AOX1</i> promoter	Allows methanol-inducible high level expression in <i>Pichia</i>
		Targets plasmid integration to the <i>AOX1</i> locus.
Sig	DNA sequence coding for an N-terminal protein secretion signal	Targets desired protein for secretion
MCS	Multiple Cloning Site	Allows insertion of your gene into the expression vector
TT	Native transcription termination and polyadenylation signal from <i>AOX1</i> gene (~260 bp)	Permits efficient transcription termination and polyadenylation of the mRNA
HIS4	Pichia wild-type gene coding for histidinol dehydrogenase (~2.4 kb) and used to complement Pichia his4 strains	Provides a selectable marker to isolate <i>Pichia</i> recombinant strains
3´ AOX1	Sequences from the <i>AOX1</i> gene that are further 3′ to the TT sequences (~650 bp)	Targets plasmid integration at the <i>AOX1</i> gene
Amp pBR322 origin	Ampicillin resistance gene <i>E. coli</i> origin of replication	Allows selection, replication, and maintenance in <i>E. coli</i>
f1 origin	Bacteriophage f1 origin of replication (458 bp)	Permits generation of single- stranded DNA for mutagenesis
Not I	Unique restriction sites	Permits linearization of vector for
Bgl II		efficient integration into the <i>Pichia</i>
Sac I		genome
Sal I		
Stu I		

Selecting a Pichia Expression Vector, continued

Selecting a Vector

If your protein is cytosolic and non-glycosylated, you may elect to express the protein intracellularly. However, there is evidence of a non-glycosylated protein being secreted without extensive modification (Despreaux and Manning, 1993). Note that the protein in question was a secreted, bacterial protein with one N-glycosylation site. Check your protein sequence for possible N-glycosylation sites (Asn-X-Ser/Thr) before cloning a cytosolic protein into a secretion vector.

If your protein is normally secreted, glycosylated, or directed to an intracellular organelle, you may wish to try secreting your protein. We recommend that you try both the native secretion signal and the α -factor signal sequence (in pPIC9) to secrete your protein. There has been better success reported with the α -factor signal sequence than with the *PHO1* signal sequence in pHIL-S1. This may be due to the lack of *KEX2*-like processing signals in the *PHO1* signal sequence (Laroche *et al.*, 1994).

pHIL-D2

Description

The details of pHIL-D2 are listed below:

- 8,209 bp nonfusion vector
- One unique *Eco*R I site
- For intracellular expression of your gene
- Requires an initiating ATG codon in a Kozak consensus sequence for proper translation initiation of your gene (Cavener and Stuart, 1991; Kozak, 1987; Kozak, 1990)
- HIS4 selection in Pichia
- For insertion at *AOX1* in GS115 or KM71, linearize with *Sac* I (generates His⁺ Mut⁺ in GS115 and His+ Mut^S in KM71)
- For insertion at *HIS4*, linearize with *Sal* I or *Stu* I (generates His⁺ Mut⁺ in GS115 and His⁺ Mut^S in KM71)
- For a gene replacement at *AOX1* in GS115, linearize with *Not* I (generates His⁺ Mut^S)

Refer to page 29 for alternate restriction sites if your insert DNA has a *Not* I, *Sac* I, *Sal* I, or *Stu* I site.

Map of pHIL-D2

The map below shows the location and size of each feature of pHIL-D2. For the details of the multiple cloning site refer to page 21. The complete sequence of pHIL-D2 is available at www.invitrogen.com or from Technical Support (page 89).

Comments for pHIL-D2: 8209 nucleotides

5' *AOX1* promoter fragment: bases 14-941 5' *AOX1* primer site: bases 868-888

EcoR I Site: bases 956-961


3' AOX1 primer site: bases 1036-1056

3' AOX1 transcription

termination (TT) fragment: bases 963-1295

HIS4 ORF: bases 4223-1689
3' AOX1 fragment: bases 4578-5334
Ampicillin resistance gene: bases 5686-6546
f1 origin of replication: bases 7043-6588

pBR322 origin: bases 7138-7757



pPIC3.5

Description

The details of pHIL-D2 are listed below:

- 7,751 bp nonfusion vector
- BamH I, SnaB I, EcoR I, Avr II, Not I unique sites
- Intracellular expression of your gene
- Requires an initiating ATG codon in a Kozak consensus sequence for proper translation initiation of your gene (Cavener and Stuart, 1991; Kozak, 1987; Kozak, 1990)
- HIS4 selection in Pichia
- For insertion at *AOX1* in GS115 or KM71, linearize with *Sac* I (generates His⁺ Mut⁺ in GS115 and His⁺ Mut^S in KM71)
- For insertion at *HIS4*, linearize with *Sal* I or *Stu* I (generates His⁺ Mut⁺ in GS115 and His⁺ Mut^S in KM71)
- For a gene replacement at *AOX1* in GS115, linearize with *Bgl* II (generates His⁺ Mut^s)

Refer to page 29 for alternate restriction sites if your insert DNA has a *Not* I, *Sac* I, *Sal* I, or *Stu* I site.

Map of pPIC3.5

The map below shows the location and size of each feature of pPIC3.5. For the details of the multiple cloning site refer to page 22. The complete sequence of pPIC3.5 is available at www.invitrogen.com or from Technical Support (page 89).

Comments for pPIC3.5: 7751 nucleotides

5' AOX1 promoter fragment: bases 1-937 5' AOX1 primer site: bases 855-875 Multiple Cloning Site: bases 938-968


3' AOX1 primer site: bases 1055-1075

3' AOX1 transcription termination (TT) fragment: bases 981-1314

HIS4 ORF: bases 4242-1708

3' AOX1 fragment: bases 4598-5354 pBR322 origin: bases 6436-5764

Ampicillin resistance gene: bases 7442-6582



pHIL-S1

Description

The details of pHIL-S1 are listed below:

- 8,260 bp fusion vector
- Xho I, EcoR I, Sma I, BamH I unique sites
- Secreted expression using the PHO1 secretion signal
- For expression, your gene must be cloned in frame with the initiation codon of the signal sequence.
- HIS4 selection in Pichia
- For insertion at AOX1 in GS115 or KM71, linearize with Sac I (generates His⁺ Mut⁺ in GS115 or His⁺ Mut^S in KM71)
- For insertion at *HIS4* in GS115 or KM71, linearize with *Sal* I or *Stu* I (generates His⁺ Mut⁺ in GS115 or His⁺ Mut^S in KM71)
- For gene replacement at AOX1 in GS115, linearize with Bgl II (generates His⁺ Mut^S)

Refer to page 30 for alternate restriction sites if your insert DNA has a *Bgl* II, *Sac* I, *Sal* I, or *Stu* I site.

Map of pHIL-S1

The map below shows the location and size of each feature of pHIL-S1. For the details of the multiple cloning site, refer to page 23. The complete sequence of pHIL-S1 is available at www.invitrogen.com or from Technical Support (page 89).


Comments for pHIL-S1: 8260 nucleotides

5' *AOX1* promoter fragment: bases 1-941 5' *AOX1* primer site: bases 856-876 PHO1 secretion signal (S): bases 942-1007 Multiple Cloning Site Region: bases 1006-1026 3' *AOX1* primer site: bases 1099-1119

3' AOX1 transcription

termination (TT) fragment: bases 1025-1190

HIS4 ORF: bases 4286-1753 3' AOX1 fragment: bases 4641-5397 pBR322 origin: bases 6556-5937 f1 origin of replication: bases 6651-7106 Ampicillin resistance gene: bases 7922-7062



pPIC9

Description

The details of pHIL-S1 are listed below:

- 8,023 bp fusion vector
- Xho I, SnaB I, EcoR I, Avr II, Not I unique sites
- Secreted expression of your gene using the -factor secretion signal
- For expression, your gene must be cloned in frame with the initiation codon of the signal sequence.
- HIS4 selection in Pichia
- For insertion at AOX1 in GS115 or KM71, linearize with Sac I (generates His⁺ Mut⁺ in GS115 and His⁺ Mut^S in KM71)
- For insertion at *HIS4*, linearize with *Sal* I or *Stu* I (generates His⁺ Mut⁺ in GS115 and His⁺ Mut^S in KM71)
- For gene replacement at AOX1 in GS115, linearize with Bgl II (generates His⁺ Mut^S)

Refer to page 30 for alternate restriction sites if your insert DNA has a *Bgl* II, *Sac* I, *Sal* I, or *Stu* I site.

Map of pPIC9K

The figure below shows the map of pPIC9. Details of the multiple cloning site are provided on page 24. The sequence of pPIC9 is available at www.invitrogen.com or from Technical Support (page 89).


Comments for pPIC9: 8023 nucleotides

5' AOX1 promoter fragment: bases 1-948 5' AOX1 primer site: bases 855-875 α -Factor secretion signal(s): bases 949-1215 α -Factor primer site: bases 1152-1172 Multiple Cloning Site: bases 1192-1241 3' AOX1 primer site: bases 1327-1347 3' AOX1 transcription

termination (TT): bases 1253-1586 *HIS4* ORF: bases 4514-1980

3' AOX1 fragment: bases 4870-5626 pBR322 origin: bases 6708-6034

Ampicillin resistance gene: bases 7713-6853



Signal Sequence Processing

Signal Sequence Processing

When cloning into the *Xho* I site of pPIC9, the secretion signal sequence between the *Xho* I site and *Sna*B I may need to be regenerated.

The processing of the α -factor mating signal sequence in pPIC9 occurs in two steps:

- 1. The preliminary cleavage of the signal sequence by the *KEX2* gene product, with the final *KEX2* cleavage occurring between arginine and glutamine in the sequence Glu-Lys-Arg * Glu-Ala-Glu-Ala, where * is the site of cleavage.
- 2. The STE13 gene product further cleaves the Glu-Ala repeats.

Optimizing Signal Cleavage

In *Saccharomyces cerevisiae*, the Glu-Ala repeats are not necessary for cleavage by *KEX2*, but the *KEX2* cleavage after Glu-Lys-Arg may be more efficient when followed by Glu-Ala repeats. A number of amino acids are tolerated at site X instead of Glu in the sequence Glu-Lys-Arg-X. These amino acids include the aromatic amino acids, small amino acids, and histidine. Proline, however, inhibits *KEX2* cleavage. For more information on *KEX2* cleavage, see (Brake *et al.*, 1984).

There are some cases where *STE13* cleavage of Glu-Ala repeats is not efficient, and Glu-Ala repeats are left on the N-terminus of the expressed protein. This is generally dependent on the protein of interest.



The *PHO1* signal sequence is atypical of signal sequences even though it is a native *Pichia* secretion signal. If cloning into the *Xho* I site, we recommend regenerating the full signal sequence between the *Xho* I and *Eco*R I sites (see page 23). However, recent evidence suggests that the *PHO1* signal sequence might have to be modified to include *KEX2*-like processing sites for efficient cleavage to occur (Laroche *et al.*, 1994).

Cloning into the Pichia Expression Vectors

Introduction

After selecting a vector into which to clone your gene of interest (see pages 11–17), develop a cloning strategy. The *AOX1* promoter and the multiple cloning site are presented on the following pages for each vector along with a summary of considerations for each vector to help you decide on a strategy.



We recommend that you transform the three supercoiled *Pichia* expression vectors into *E. coli* to prepare permanent stock.

- Resuspend each vector in 10 μ L sterile water to prepare a 1 μ g/ μ L solution. Store the stock solution at -20° C.
- Use the stock solution to transform competent *E. coli* and select transformants on LB agar plates containing 50–100 μg/mL ampicillin (LB-Amp).

General Considerations

The following are some general considerations applicable to all vectors.

- The codon usage in *Pichia* is believed to be the same as *Saccharomyces cerevisiae* because many genes have proven to be cross-functional.
- Maintain plasmid constructions in a *rec*A mutant *E. coli* strain such as the TOP10F′ strain provided in the kit.
- The native 5′ end of the *AOX1* mRNA is noted in each multiple cloning site. This information is necessary to calculate the size of the expressed mRNA of the gene of interest.
- Translation termination is determined by either stop codons in the gene of interest or in the 3′ *AOX1* sequence. The stop codons in the 3′ *AOX1* sequence are noted in each figure on the following pages.
- The premature termination of transcripts due to "AT rich regions" has been observed in *Pichia* and other eukaryotic systems (Henikoff and Cohen, 1984; Irniger *et al.*, 1991; Scorer *et al.*, 1993; Zaret and Sherman, 1984). If you are expressing a gene with high AT content, refer to page 53.
- The predicted protease cleavage sites for the *PHO1* and α -factor signal sequences are indicated in each figure.
- If you are attempting to secrete a protein using its native secretion signal, we recommend that you also try pPIC9 in parallel. When using pPIC9, the open reading frame (ORF) of the mature gene of interest is cloned in frame and downstream of the α -factor.

General Cloning Strategies

Strategies generally fall into three different categories:

- 1. Ligation of a compatible restriction fragment:
 - a. Forced (directional) insertion involving the use of two different sites in the multiple cloning site (for pPIC3.5, pHIL-S1, or pPIC9 vectors).
 - b. Ligation of the fragment with the same restriction end on both ends into a single, compatible site (e.g. *Eco*R I cloning in pHIL-D2).
- 2. PCR amplification of the fragment containing the gene of interest in such a way that compatible restriction ends are generated for ligation into the appropriate vector.
- 3. Direct cloning of an amplified fragment containing the gene of interest via the TA Cloning[®] Kit (see page 88 for ordering), followed by subcloning of a compatible fragment into the appropriate *Pichia* expression vector.

Cloning Procedures

Refer to (Ausubel *et al.*, 1994), pages 3.16.1 to 3.17.3. or (Sambrook *et al.*, 1989), pages 5.10 to 5.13. for help with cloning.


Bacterial Transformation

Once you have decided on a cloning strategy, you will need to prepare competent *E. coli* cells for transformation before setting up your ligation reactions. See *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994) or *Molecular Biology: A Laboratory Manual* (Sambrook *et al.*, 1989) for preparation of electrocompetent or chemically competent *E. coli* or use your laboratory's procedure.

For a high-efficiency and convenient transformation, we recommend One Shot[®] TOP10 Chemically Competent *E. coli* or One Shot[®] TOP10 Electrocompetent Cells, which are available separately from Invitrogen (see page 88 for ordering information).

P_{AOX1} and Multiple Cloning Site of pHIL-D2

The figure below shows the detail of the multiple cloning site and surrounding sequences.



Special Considerations

For pHIL-D2, the fragment containing the gene of interest should have a yeast consensus sequence (Romanos *et al.*, 1992). An example of a yeast consensus sequence is provided below. The ATG initiation codon is shown underlined.

(A/Y)A(A/T)AATGTCT

Note that other sequences are also possible. Although not as strong as the mammalian Kozak translation initiation sequence, the yeast consensus sequence is thought to have a 2- to 3-fold effect on the efficiency of translation initiation.

- Shorter, 5' untranslated leaders reportedly work better in AOX1 expression.
 In pHIL-D2, make the untranslated region as short as possible when cloning your gene.
- If your insert has a *Not* I site, refer to page 29 for alternate restriction sites to linearize your plasmid for *Pichia* transformation.

P_{AOX1} and Multiple Cloning Site of pPIC3.5

The sequence below shows the detail of the multiple cloning site and surrounding sequences.

AOX1 mRNA 5'end (824)		5' AOX	1 Primer Site (855-8	75)		
TTATCATCAT TATTAGCTTA	CTTTCATAAT	TGCGACTGGT	TCCAATTGAC	AAGCTTTTGA	TTTTAACGAC	
			BamH I	SnaB I EcoR	I Avr II	
TTTTAACGAC AACTTGAGAA	GATCAAAAAA	CAACTAATTA	TTCGAAGGAT	CCTACG <u>TAG</u> A	ATTCCCTAGG	
Not I						
GCGGCCGCGA ATTAATTCGC	$\mathtt{CT}\underline{\mathtt{TAG}}\mathtt{ACATG}$	ACTGTTCCTC	AGTTCAAGTT	GGGCACTTAC	GAGAAGACCG	
3' AOX 1 Primer Site (1055-1075)						
GTCTTGC <u>TAG</u> ATTCTAATCA	AGAGGATGTC	AGAATGCCAT	TTGCCTGAGA	GATGCAGGCT	TCATTTTGA	
AOX1 mRNA 3' end (1146) ↓						
TACTTTTTTA TTTGTAACCT	ATATAGTATA	GGATTTTTT	TGTCATTTTG	TTTCTTC		

Special Considerations

For pPIC3.5, the fragment containing the gene of interest should have a yeast consensus sequence (Romanos *et al.*, 1992). An example of a yeast consensus sequence is provided below. The ATG initiation codon is shown underlined.

(A/Y)A(A/T)AATGTCT

Note that other sequences are also possible. Although not as strong as the mammalian Kozak translation initiation sequence, the yeast consensus sequence is thought to have a 2- to 3-fold effect on the efficiency of translation initiation.

- Shorter, 5' untranslated leaders reportedly work better in AOX1 expression.
 In pPIC3.5, make the untranslated region as short as possible when cloning your gene.
- If you are digesting with BamH I and SnaB I or SnaB I and EcoR I, digest with SnaB I first. If you digest with BamH I or EcoR I first, the SnaB I site will be too close to the end of the DNA and will not digest properly.
- If your insert has a *Bgl* II, *Sac* I, *Sal* I, *or Stu* I site, refer to page 29 for alternate restriction sites to linearize your plasmid for *Pichia* transformation.

P_{AOX1} and Multiple Cloning Site of pHIL-S1

The sequence below shows the detail of the multiple cloning site and surrounding sequences.

Special Considerations

- The fragment containing the gene of interest must be cloned in frame with the secretion signal open reading frame.
- If the *Xho* I site is used for cloning, it must be recreated for efficient cleavage of the fusion protein to occur. It is part of the *PHO1* signal peptide sequence.
- An initiating ATG is provided by the signal sequence. Translation will initiate at the ATG closest to the 5′ end of the mRNA.
- If your insert has a *Not* I site, refer to page 29 for alternate restriction sites to linearize your plasmid for *Pichia* transformation.
- The *PHO1* cleavage site has been confirmed for several different fusion proteins by N-terminal peptide sequencing.
- In -factor (pPIC9) or native general, more success has been reported with the secretion signals than with pHIL-S1. This may be due to the lack of *KEX2*-like processing signals (Laroche *et al.*, 1994).

^{*} If the Xho I site (which is part of the PHO1 signal cleavage sequence) is used for cloning, it must be recreated in order for efficient cleavage of the fusion protein to occur.

P_{AOX1} and Multiple Cloning Site of pPIC9

The sequence below shows the detail of the multiple cloning site and surrounding sequences.

AOX1 mRNA 5' end (824) ACAGCAATAT ATAAACAGAA GGAAGCTGCC CTGTCTTAAA CCTTTTTTTT TATCATCATT ATTAGCTTAC 5' AOX1 Primer Site (855-875) TTTCATAATT GCGACTGGTT CCAATTGACA AGCTTTTGAT TTTAACGACT TTTAACGACA ACTTGAGAAG α-Factor (949-1215) ATCAAAAAAC AACTAATTAT TCGAAGGATC CAAACG ATG AGA TTT CCT TCA ATT TTT ACT GCA Met Arg Phe Pro Ser Ile Phe Thr Ala GTT TTA TTC GCA GCA TCC TCC GCA TTA GCT GCT CCA GTC AAC ACT ACA ACA GAA GAT Val Leu Phe Ala Ala Ser Ser Ala Leu Ala Ala Pro Val Asn Thr Thr Thr Glu Asp GAA ACG GCA CAA ATT CCG GCT GAA GCT GTC ATC GGT TAC TCA GAT TTA GAA GGG GAT Glu Thr Ala Gln Ile Pro Ala Glu Ala Val Ile Gly Tyr Ser Asp Leu Glu Gly Asp TTC GAT GTT GCT GTT TTG CCA TTT TCC AAC AGC ACA AAT AAC GGG TTA TTG TTT ATA Phe Asp Val Ala Val Leu Pro Phe Ser Asn Ser Thr Asn Asn Gly Leu Leu Phe Ile Xho I α-Factor Primer Site (1152-1172) AAT ACT ACT ACT GCC AGC ATT GCT GCT AAA GAA GAA GGG GTA TCT CTC GAG AAA AGA Asn Thr Thr Ile Ala Ser Ile Ala Ala Lys Glu Glu Gly Val Ser Leu Glu Lys Arg SnaB I EcoR I GAG GCT GAA GCT TAC GTA GAA TTC CCT AGG GCC GCC AAT TAA TTCGCCT<u>TAG</u> Glu Ala Glu Ala Tyr Val Glu Phe Pro Arg Ala Ala Ala Asn *** ACATGACTGT TCCTCAGTTC AAGTTGGGCA CTTACGAGAA GACCGGTCTT GCTAGATTCT AATCAAGAGG 3' AOX1 Primer Site (1327-1347) ATGTCAGAAT GCCATTTGCC TGAGAGATGC AGGCTTCATT TTTGATACTT TTTTATTTGT AACCTATATA AOX1 mRNA 3' end (1418) GTATAGGATT TTTTTTGTCA

Special Considerations

- The fragment containing the gene of interest must be cloned in frame with the secretion signal open reading frame.
- If the *Xho* I site is used for cloning, the sequence between the *Xho* I site and the *Sna*B I site encoding the KEX2 site (Glu-Lys-Arg-X) must be recreated for efficient cleavage of the fusion protein to occur. It is part of the α -factor signal peptide sequence. Refer to the discussion on page 18.
- An initiating ATG is provided by the signal sequence. Translation will initiate at the ATG closest to the 5′ end of the mRNA.
- If your insert has a *Bgl* II site, refer to page 30 for alternate restriction sites to linearize your plasmid for *Pichia* transformation.

^{*} If cloning into the Xho I site, the sequence between the Xho I site and SnaB I site (underlined) must be recreated in order for efficient cleavage of the fusion protein to occur.

Transformation into E. coli

Introduction

At this point you have ligation reactions that you will transform by chemical means or electroporation into competent *E. coli* cells (TOP10F´ or equivalent). For procedures to prepare competent cells, refer to *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994) or *Molecular Biology: A Laboratory Manual* (Sambrook *et al.*, 1989).

Analyzing Transformants

- 1. After transformation, plate the transformation mix onto LB plates with $50-100 \, \mu g/mL$ ampicillin (see **Appendix**, page 59), and select ampicillin resistant colonies.
- 2. Pick 10 ampicillin resistant transformants and inoculate into LB medium with $50-100 \mu g/mL$ ampicillin. Grow overnight at $37^{\circ}C$ with shaking.
- 3. Isolate plasmid DNA by miniprep for restriction analysis and sequencing (see below). To sequence the *Pichia* expression vectors, use the primers provided.
- 4. Make a glycerol stock of your desired clone for safekeeping by combining 0.85 mL of an overnight bacterial culture with 0.15 mL of sterile glycerol. Mix by vortexing and transfer to a labeled storage tube. Freeze the tube in liquid nitrogen or a dry ice/ethanol bath and store at -70°C.
- 5. After you confirm your construct by sequencing, proceed to **Preparing Transforming DNA**, page 27.

Sequencing Recombinant Clones

We strongly recommend that you sequence your construct before transforming into *Pichia* to confirm the following:

- The correct reading frame (for secretion)
- An ATG in the proper context for eukaryotic translation initiation

Use the primers listed on the next page to sequence your constructs. Resuspend each primer in 20 μ L sterile water to prepare a stock solution of 0.1 μ g/ μ L. For the location of the priming sites, see pages 21–24.

For sequencing protocols, refer to Unit 7 in *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994) or Chapter 13 in *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989).

Transformation into E. coli, continued

5´ *AOX1* Sequencing Primer

- is 5' GACTGGTTCCAATTGACAAGC 3'
- hybridizes 5' of the MCS in the AOX1 promoter region
- allows the determination of the 5′ AOX1-gene of interest junction
- confirms that the ORFs are preserved where necessary

3´ *AOX1* Sequencing Primer

- is 5′ GCAAATGGCATTCTGACATCC 3′
- hybridizes just 3′ of the MCS in the 3′ AOX1 (TT) region
- allows the determination of the 3′ AOX1-gene of interest junction

α-Factor Sequencing Primer

- is 5'-TACTATTGCCAGCATTGCTGC-3'
- hybridizes within the -factor leader region in pPIC9
- allows the determination of the 5' end of the gene of interest
- confirms that the ORFs are preserved where necessary

Preparing Transforming DNA

Introduction

You should have a *Pichia* multi-copy expression vector with your gene of interest cloned in the correct orientation for expression. The table below describes what you will be doing in the next few sections.

Step	Action	Pages
1	Prepare your DNA for transformation	28
2	Grow GS115 or KM71 to prepare spheroplasts	32
3	Prepare spheroplasts for transformation	33
4	Transform GS115 or KM71 with your DNA	35
5	Select His ⁺ transformants and characterize for Mut ⁺ /Mut ^S phenotype	38
6	Test 10 His ⁺ Mut ⁺ and 10 His ⁺ Mut ⁵ by PCR for integration of your gene	43



We recommend isolating both His⁺ Mut⁺ and His⁺ Mut^S *Pichia* transformants as it is difficult to predict beforehand what construct will best express your protein (see pages 69 and 70). By linearizing your construct DNA in the 5′ AOX1 region or in the HIS4 gene and using GS115 (Mut⁺) and KM71 (Mut^S), you can easily isolate Mut⁺ and Mut^S recombinants. Plan on using ~10 µg digested DNA for each transformation.

Preparing Plasmid DNA

Plasmid DNA for *Pichia* transformation should be at least pure enough for restriction digestion; however, the cleaner the DNA, the more efficient the transformation. We recommend the PureLink™ HiPure Plasmid Miniprep Kit (see page 88) to prepare plasmid DNA for routine *Pichia* transformations. Refer to www.invitrogen.com or contact Technical Support for more information on a large selection of plasmid purification columns. You may prepare plasmid DNA using alkaline lysis, phenol:chloroform extraction, and ethanol precipitation.

Preparing Transforming DNA, continued

Linearizing Plasmid DNA

We recommend that you linearize your vector in such a manner to generate both Mut⁺ and Mut⁵ recombinants. It is possible that one phenotype will express your protein of interest better than the other will.

- To isolate His⁺ Mut⁺ transformants of GS115, linearize all constructs with *Sal* I, *Stu* I, or *Sac* I.
- To isolate His⁺ Mut^S transformants of KM71, linearize plasmid constructs with *Sal* I, *Stu* I, or *Sac* I.

Note: If your insert DNA has all three sites, please see the next page.

• To isolate His⁺ Mut^S transformants of GS115, linearize plasmid constructs with *Not* I (pHIL-D2) or *Bgl* II (pPIC3.5, pHIL-S1, and pPIC9)

Note: If you wish to generate recombinants that are Mut^S, use KM71 because it is much easier and more efficient to generate Mut^S recombinant strains using single crossover events than double crossover events (e.g. insertions at *AOX1* or *his4* as opposed to gene replacement at *AOX1*).

Procedure

- 1. Digest both your construct and the parent vector. You will transform GS115 and/or KM71 with the parent vector as a background control for expression.
- 2. Analyze a small portion of your digest by agarose gel electrophoresis to confirm complete digestion of your fragment. The number of transformants and frequency of targeting will be reduced if digestion is not complete.
- 3. Extract the digest with phenol:chloroform:isoamyl alcohol (25:24:1) and ethanol precipitate the digested DNA. Resuspend DNA pellet in 10–20 μ L of TE buffer. It is not necessary to purify the fragment containing your gene away from the rest of the plasmid.
- 4. Store at –20°C until ready to transform.

Preparing Transforming DNA, continued

Alternate Restriction Sites

If your insert DNA contains *Sac* I, *Sal* I, and *Stu* I sites, you need to linearize your construct with another enzyme. Use the following table to select another enzyme. A single digestion which linearizes the vector in either of the *AOX1* recombination sequences allows integration, but at lower efficiencies. Remember to digest the parent vector with the same enzyme when preparing your DNA samples for transformation.

pPHIL-D2

Restriction Enzyme	5´ <i>AOX1</i> 14–940 bp	3´ AOX1 4,577–5,333 bp	Vector backbone 5,333+ bp	HIS4 gene 1,688–4,222 bp
Sac I	221	_	_	_
Pme I	424	_	_	_
Ври 1102 І	599	_	_	_
Nsi I	689	_	_	_
Xcm I	711	_	_	_
Not I	8	5,337	_	_
Dra I*	424	5,169; 5,311	5,896; 6,588	_
Sal I	_	_	_	2,887
Stu I	_	_	_	2,972
BspE I	_	_	_	3,554

pPIC3.5

Restriction Enzyme	5´ AOX1 1–937 bp	3´ AOX1 4,616–5,393 bp	Vector backbone 5,393+ bp	HIS4 gene 1,715–4,249 bp
Sac I	209	_	_	_
Pme I	414	_	_	-
Ври 1102 І	589	_	_	_
Nsi I	678	_	_	-
Xcm I	699	_	_	_
Bgl II	2	5363	_	_
Dra I*	414	5,201; 5,343	6,534; 6,553; 7,245	_
Sal I	_	_	_	2,919
Stu I	_	_	_	3,004
BspE I	_	_	_	3,586

Preparing Transforming DNA, continued

Alternate Restriction Sites, continued

pHIL-S1

Restrictio n Enzyme	5′ <i>AOX1</i> 1–940 bp	3′ AOX1 4,639–5,395 bp	Vector backbone 5,395+ bp	HIS4 gene 1,750–4,284 bp
Sac I	209	_	_	_
Pme I	412	_	_	_
Ври 1102 І	587	_	_	_
Nsi I	677	_	_	_
Xcm I	699	_	_	_
Bgl II	2	5,394	_	_
Dra I*	412	5,232; 5,374	7,021; 7,713	_
Sal I	_	_	_	2,950

pPIC9

Restriction Enzyme	5´ AOX1 1–948 bp	3´ AOX1 4,881–5,638 bp	Vector backbone 5,638+ bp	HIS4 gene 1,980–4,514 bp
Sac I	209	_	_	_
Pme I	414	_	_	_
Ври 1102 І	589	_	_	_
Nsi I	678	_	_	_
Xcm I	699	_	_	_
Bgl II	2	5,622	_	_
Dra I*	414	5,460; 5,602	6,793; 6,812; 7,504	-
Sal I	_	_	_	3,178
Stu I	_	_	_	3,263
BspE I	_	_	_	3,845

^{*}Restriction sites are used to generate gene replacements at *AOX1* in GS115 only.

Growing Pichia for Spheroplasting

Introduction

In general, spheroplasting and electroporation (page 73) provide the highest efficiency of transformation for most researchers (10³ to 10⁴ transformants per μg DNA). *Pichia* can also be transformed using PEG 1000 (page 74) or lithium chloride (page 76). These two protocols, particularly lithium chloride, do **not** perform as well as spheroplasting or electroporation. If you do not have an electroporation device, we recommend spheroplasting or using the PEG 1000 method. Transformation in *Pichia* is less efficient than for *Saccharomyces*. For references on general yeast transformation, see (Cregg *et al.*, 1985; Hinnen *et al.*, 1978).

Spheroplasting

The cell wall of yeast prevents uptake of DNA. To enable yeast to take up DNA, it is necessary to partially remove the cell wall. Zymolyase is a β -glucanase that hydrolyzes the glucose polymers with $\alpha 1,\! 3$ linkages in the cell wall. Addition of Zymolyase partially digests the cell wall. It is critical not to overdigest the cell wall as doing so will cause the cells to die. Zymolyase digestion is monitored by the sensitivity of the cells to SDS. Aliquots of cells are added to SDS, lysing the spheroplasts. This causes a clearing of the solution that is monitored by the absorbance (light-scattering) at 800 nm. It has been empirically determined that when 70% spheroplasting has been achieved, digestion is optimal. Cells are then washed with an isotonic solution to remove the enzyme and incubated with DNA. The cells are resuspended in sorbitol to facilitate cell wall regeneration and plated.

Preparing Media

Prepare the following media several days in advance and store at 4°C (see **Appendix**, pages 61–64 for details):

YPD (Yeast extract Peptone Dextrose) medium, 1 liter

YPD plates, 1 liter

RDB (Regeneration Dextrose Base) plates, 1 liter

RDHB (Regeneration Dextrose Histidine Base) plates, 1 liter

Prepare the following solution on the day of transformation and maintain at 45°C: 5% SDS solution in water

RD (Regeneration Dextrose), molten agarose, 100 mL

Growing Pichia for Spheroplasting, continued

Solutions

Spheroplasting and Transformation Reagents

Provided:

1 M Sorbitol

SE: 1 M sorbitol, 25 mM EDTA, pH 8.0

DTT: 1 M DTT in water

SCE: 1 M sorbitol, 1 mM EDTA and 10 mM sodium citrate buffer, pH

5.8

CaS: 1 M sorbitol, 10 mM Tris-HCl, pH 7.5, 10 mM CaCl₂

Zymolyase: 3 mg/mL in water

40% PEG: 40% (w/v) PEG 3350 (Reagent grade) in water

CaT: 20 mM Tris, pH 7.5 and 20 mM CaCl₂ SOS: 1 M sorbitol, 0.3X YPD, 10 mM CaCl₂

Prepared fresh for each transformation:

SED: 19 mL of SE and 1 mL of 1 M DTT (see page 33) PEG/CaT: 1:1 mixture of 40% PEG and CaT (see page 35)

Procedure

- 1. Streak GS115 or KM71 onto a YPD plate such that isolated, single colonies will grow. Incubate the plate at 28–30°C for 2 days.
- 2. Inoculate 10 mL of YPD in a 50 mL conical tube or 100 mL shake flask with a single colony of GS115 or KM71 from the YPD plate and grow overnight at 28–30°C with vigorous shaking (250–300 rpm). You may store this culture at 4°C for several days.
- 3. Place 200 mL of YPD in each of three 500 mL culture flasks. Inoculate the flasks with 5, 10, and 20 μ L of cells from the culture made in Step 2 and incubate them overnight with vigorous shaking (250–300 rpm) at 28–30°C.
- 4. The next morning, bring the transformation solutions (SE, SCE, Sterile Water, SOS, PEG, CaS, CaT, 1 M sorbitol) provided in the kit, the RDB plates (for plating transformants), and the RDHB plates (for viability control) to room temperature.
- 5. Check the OD_{600} of each of the three culture flasks. Harvest the cells from the culture that has an OD_{600} between 0.2 and 0.3. Centrifuge the cells at room temperature for 5–10 minutes at 1,500 × g. Decant the supernatant and discard the other cultures. Proceed to **Preparing Spheroplasts**, page 33.

Note: If the cultures are all over 0.3, choose one of the cultures and dilute (1:4) with fresh medium and incubate at $28-30^{\circ}$ C until the OD₆₀₀ is between 0.2 and 0.3 (2–4 hours). Harvest the cells and proceed as in Step 5, above.

Preparing Spheroplasts

Before Starting

You should have a cell pellet from Step 5, page 32.

- Prepare 100 mL of molten RD agarose and keep at 45°C (see Appendix, page 64)
- Thaw one tube of 1 M DTT (provided in the kit)
- Prepare fresh SED for one batch of spheroplasts as follows:

Using sterile technique, transfer 19 mL of SE (provided) to an appropriate sterile container (e.g., 50 mL conical tube). Add 1 mL of 1 M DTT and mix well. For best results this solution of SED should be made and used immediately.



The quality and freshness of DTT is critical for a successful spheroplast preparation. The 1 M DTT provided is analytical reagent grade and must be stored at -20° C.

Washing the Cells

- 1. Wash the cells from Step 5, page 32 by resuspending the pellet in 20 mL of sterile water (provided). Resuspend the pellet by swirling the tube. Transfer to a sterile, 50 mL conical tube.
- 2. Pellet the cells by centrifugation at $1,500 \times g$ for 5 minutes at room temperature. Decant and discard the supernatant. The cell pellet will be used to prepare spheroplasts.
- 3. Wash the cell pellet once by resuspending in 20 mL of fresh SED, prepared above and centrifuge at $1,500 \times g$ for 5 minutes at room temperature.
- 4. Wash the cells once with 20 mL of 1 M sorbitol and centrifuge as described in Step 2.
- 5. Resuspend the cells by swirling in 20 mL of SCE buffer and divide the suspension into two 50 mL conical tubes (~10 mL each).
- 6. Remove one tube of Zymolyase from -20°C and place it on ice. Mix well by flicking the tube several times. Zymolyase is provided as a slurry and does not go into solution. It is important to mix the slurry thoroughly before each use to ensure addition of a consistent amount of Zymolyase.

Adding Zymolyase

Use one tube of cells prepared above to determine the optimal time of digestion with Zymolyase to make spheroplasts. After you determine the optimal, use the other tube of to make spheroplasts.

Zymolyase digests the cell wall and makes the cells extremely fragile. **Handle the sample gently**. The moment after adding Zymolyase, the cell wall begins to get digested.

- Prepare at least 20 mL of a 5% SDS solution (not provided) for use below.
- Set your UV-Vis spectrophotometer to 800 nm and blank with 800 μ L 5% SDS and 200 μ L SCE.
- Set up 17 sterile microcentrifuge tubes and label them 0, 2, 4, 5, 6, 7, 8, 9, 10, 15, 20, 25, 30, 35, 40, 45, and 50. Add 800 μL of 5% SDS to each tube.

Preparing Spheroplasts, continued

Adding Zymolyase, continued

- 1. From one tube of cells (Step 5, page 33), withdraw 200 µL cells and add to the tube marked "0". This is your zero time point. Set the tube aside on ice.
- 2. Add 7.5 μL of Zymolyase to the same tube of cells, mix it gently by inversion, and incubate the cells at 30°C. **Do not shake the sample.** This sample will be used to establish the incubation time for optimal spheroplasting as described below. Keep the second tube of cells at room temperature for use in Step 6 below. Keep the remainder of the Zymolyase on ice.
- 3. Monitor the formation of the spheroplasts as follows: At time 2 minutes, withdraw 200 μ L of cells (from the suspension in Step 2) and add to the tube marked "2". Repeat at time t= 4, 5, 6, 7, 8, 9, 10, 15, 20, 25, 30, 35, 40, 45, and 50 minutes after adding Zymolyase. Read the OD₈₀₀ for all samples.
- 4. Determine the percent of spheroplasting for each time point using the equation:

```
% Spheroplasting = 100 - [(OD_{800} \text{ at time } t/OD_{800} \text{ at time } 0) \times 100]
```

For example: At time t = 0, the $OD_{800} = 0.256$

At time t = 15, the $OD_{800} = 0.032$

Calculation: % spheroplasting = $100 - [(0.032/0.256) \times 100]$

 $= 100 - [(0.125) \times 100]$

= 100 - 12.5

= 87.5%

5. Determine the time of incubation that results in approximately 70% spheroplasting. This time of incubation is variable due to differences in lots of Zymolyase. In Invitrogen labs, it takes approximately 15–40 minutes of Zymolyase treatment to achieve optimal spheroplasting.

Note: It is important to establish the minimum time required for the desired amount of spheroplasting. Prolonged incubation with Zymolyase is deleterious to spheroplasts and will result in lower transformation efficiency.

- 6. Add $7.5 \,\mu$ L Zymolyase to the remaining tube of cells as described in Step 1 above. Incubate the tube at 30°C for the time that was established in Step 5 to obtain the optimal level (70%) of spheroplasting.
- 7. Harvest the spheroplasts by centrifugation at $750 \times g$ for 10 minutes at room temperature. Decant and discard the supernatant.
- 8. Wash the spheroplasts once with 10 mL of 1 M sorbitol (**gently** disperse the pellet by tapping the tube, **do not vortex**). Collect the spheroplasts by centrifugation at $750 \times g$ for 10 minutes at room temperature.
- 9. Wash the spheroplasts once with 10 mL of CaS and centrifuge as in Step 7. **Gently** resuspend the spheroplasts in 0.6 mL of CaS. The spheroplasts must be used immediately (up to 30 minutes) for transformation (page 35). They cannot be stored for much longer. This preparation yields enough spheroplasts for six transformations.

Transforming Pichia

Before Starting

Make sure your RDB plates are at room temperature and that you have molten RD top agarose available. Thaw your linearized DNA and keep on ice. You should have the following:

- Your construct linearized with Sal I, Stu I, or Sac I to favor isolation of His⁺ Mut⁺ recombinants in GS115
- Your construct linearized with Sal I, Stu I, or Sac I to favor isolation of His⁺ Mut^S recombinants in KM71
- Your construct linearized with Not I, Bgl II, or equivalent to favor isolation of His⁺ Mut^S recombinants in GS115
- Parent plasmid linearized with same restriction enzyme

Controls should include no DNA or linearized pBR322 DNA and plasmid only (no cells) to check for contamination.

Procedure

- 1. For each transformation, dispense $100~\mu L$ of the spheroplast preparation from Step 9 (previous page) into a sterile 15 mL snap-top Falcon 2059 tube (or equivalent).
- 2. Add 10 µg of DNA and incubate the tube at room temperature for 10 minutes.
- 3. During the 10 minute incubation, make a fresh PEG/CaT solution. Since each transformation requires 1.0 mL of the PEG/CaT solution, calculate the amount you need and prepare this volume by adding together equal volumes of 40% PEG and CaT (a 1:1 solution).
- 4. Add 1.0 mL of fresh PEG/CaT solution to the cells and DNA, mix gently, and incubate at room temperature for 10 minutes.
- 5. Centrifuge the tube at $750 \times g$ for 10 minutes at room temperature and carefully aspirate the PEG/CaT solution. Invert the tube and tap it gently to drain the excess PEG/CaT solution.
- 6. Resuspend the pellet of transformed cells in 150 μL of SOS medium and incubate it at room temperature for 20 minutes.
- 7. Add 850 µL of 1 M sorbitol. Proceed to **Plating**, next page.

Transforming Pichia, continued

Plating

Pichia spheroplasts need to be plated in top agarose or agar to protect them from lysis prior to selection.

- 1. Mix together $100-300 \,\mu\text{L}$ of each spheroplast-DNA solution from Step 7, previous page, with $10 \,\text{mL}$ of molten RD agarose and pour on RDB plates. Allow the top agarose to harden.
 - **Note:** There is enough of the spheroplast-DNA solution to plate duplicate and triplicate plates.
- 2. Invert plates and incubate at 28–30°C. Transformants should appear in 4–6 days.
- 3. For cell viability: Mix 100 μ L of spheroplasts with 900 μ L of 1 M sorbitol.
- Mix 100 μL of this diluted sample with 10 mL of molten RDH and pour on a RDHB plate. Allow top agarose to harden.
- 5. Invert the plates and incubate at 28–30°C. Appearance of colonies after 4–6 days demonstrates that the spheroplasts can regenerate into dividing cells.

Evaluating Your Transformation Experiment

After 4–6 days, His $^+$ transformants on your sample plates will become apparent. Transformation efficiency using the spheroplast method is generally 10^3 to 10^4 His $^+$ transformants/µg of DNA. There should be no colonies on the "No DNA", pBR322 plate, or the plasmid only (no cells) plate

Optional Method

Plating in top agarose can cause the transformants to be on top or be imbedded in the top agarose, making it difficult to pick and patch colonies as described in the next section. The following protocol allows you to collect the transformants and re-plate them directly onto plates without using top agarose.

- 1. Scrape the agarose containing the His⁺ transformants with a sterile spreader into a sterile, 50 mL, conical centrifuge tube and mix with 20 mL sterile deionized water. Vortex the suspension vigorously to separate the cells from the agarose.
- 2. Filter the suspension through 4 folds of sterile cheesecloth. Centrifuge the filtrate at $1,500 \times g$ for 5 minutes at room temperature. Centrifugation pellets the cells on the bottom of the tube and any remaining agarose on top of the cells.
- 3. Remove the agarose pellet carefully from the top of the cells by gently shaking the tube to disperse only the agarose pellet into the water. Decant the supernatant with the agarose pellet.
- 4. Resuspend the cell pellet in 5 mL of sterile deionized water and sonicate for 10 seconds using a microtip and 20–30% power. Sonicate to get the cells into solution and not to lyse the cells.
- 5. Dilute cells by 10⁴ and plate 50 μL and 100 μL onto MD plates. Incubate the plates overnight at 30°C. Proceed to **Screening for Mut**⁺ **and Mut**^S **Transformants,** page 38.

Transforming Pichia, continued

Chemically Competent *Pichia* Cells

The *Pichia* EasyComp[™] Kit (see page 88) provides enough reagents to produce 6 preparations of competent cells; each preparation yields enough competent cells for 20 transformations. You may use these cells immediately or store them frozen for future use. Each 50 μ L aliquot of competent *Pichia* cells with 3 μ g linearized plasmid DNA yields approximately 50 colonies on selective medium. Contact Technical Support for more information (page 89).

Screening by Functional Assay

Some researchers have used a functional assay to directly test for high expressing *Pichia* recombinant clones without first screening for Mut^S or Mut⁺ phenotypes. After testing for high expression, be sure to also check the Mut phenotype. This will help you optimize expression of your recombinant clone.

Screening for Mut⁺ and Mut^S Transformants

Introduction

After colonies appear on plates of His⁺ GS115 transformants, score for the Mut⁺ and Mut^S phenotypes. The kit contains two strains that will provide examples of Mut⁺ and Mut^S phenotypes. GS115 Albumin is Mut^S and GS115 β -Gal is Mut⁺. His⁺ KM71 recombinants do not need to be screened for their Mut phenotype as they all will be Mut^S.

Remember also to isolate two control strains for background protein expression in *Pichia*. One control is the parent plasmid linearized in such a way to generate His⁺ Mut^S transformants. The other control is the parent plasmid linearized to generate His⁺ Mut⁺ transformants.

Screening for His⁺ Mut⁺ in GS115

Transformation of GS115 with *Sal* I- or *Stu* I-linearized constructs favor recombination at the *HIS4* locus. Most of the transformants should be Mut+; however, with the presence of the *AOX1* sequences in the plasmid, there is a chance that recombination will occur at the *AOX1* locus, disrupting the wild-type *AOX1* gene and creating His⁺ Mut^S transformants. Testing on Minimal Dextrose (MD) and Minimal Methanol (MM) plates allows you to isolate His⁺ Mut⁺ transformants (see next page).

His+ Mut^S in KM71

There is no need to test recombinants for the Mut phenotype in KM71; all His⁺ transformants in KM71 will be Mut^S because of the disruption of the *AOX1* gene (*aox1::ARG4*). Transforming KM71 with *Sal* I- or *Stu* I- linearized plasmid constructs favor recombination at the *HIS4* locus while *Sac* I-linearized plasmid constructs favor recombination at the 5′ region of the *AOX1* gene. Purify His⁺ transformants on minimal plates without histidine to ensure pure clonal isolates before testing for expression (see page 45) or confirming integration by PCR (page 43).

Screening for His⁺ Mut^S in GS115

Transforming GS115 with *Not* I-linearized pHIL-D2 or *Bgl* II-linearized pPIC3.5, pHIL-S1, and pPIC9 constructs favors recombination at the *AOX1* locus. Displacement of the alcohol oxidase (*AOX1*) structural gene occurs at a frequency of 5–35% of the His⁺ transformants. Patching or replica-plating on Minimal Dextrose (MD) versus Minimal Methanol (MM) plates can readily distinguish Mut⁺ and Mut^S transformants.

Because Mut^S transformants do not produce alcohol oxidase (the product of the AOX1 gene), they cannot efficiently metabolize methanol as a carbon source; therefore, they grow poorly on minimal methanol (MM) medium. This slow growth on methanol can be used to distinguish His⁺ transformants, in which the AOX1 gene has been disrupted (His⁺ Mut^S), from His⁺ transformants with an intact AOX1 gene (His⁺ Mut⁺).

Screening for Mut⁺ and Mut^S Transformants, continued

Materials Needed

You can prepare the following media (see page 64) and materials several days in advance, and store at 4°C:

- Minimal Dextrose (MD) agar plates, 1 liter
- Minimal Methanol (MM) agar plates, 1 liter
- Sterile toothpicks and Scoring Templates (see page 42)
- Streak out the strains GS115 Albumin (His $^+$ Mut s) and GS115 β -Gal (His $^+$ Mut $^+$) on an MD or MGY plate as controls for Mut $^+$ and Mut s growth on MD and MM plates.



In contrast to His⁺ Mut^S transformants generated using a construct linearized with *Not* I or *Bgl* II, most of the His⁺ transformants generated by the *Sac* I, *Sal* I ,or *Stu* I-digested construct should be Mut⁺ arising from gene insertion events at the *his4* or *AOX1* loci, leaving an intact *AOX1* locus.

His⁺ Mut^S or His⁺ Mut⁺ in GS115

Use the plates containing His⁺ transformants, and screen for the Mut⁺ and Mut^S phenotypes as described below.

- Using a sterile toothpick, pick one His⁺ colony, and streak or patch it in a regular pattern on an MM plate and an MD plate.
 Note: Make sure to patch the MM plate first.
- 2. Use a new toothpick for each transformant, and continue picking transformants until you have patched 100 transformants (2–3 plates).
- 3. To differentiate Mut⁺ from Mut^S, make one patch for each of the controls (GS115/His⁺ Mut^S Albumin and GS115/His⁺ Mut⁺ β -gal) onto the MD and MM plates.
- 4. Incubate the plates at 30°C for 2 days.
- 5. After 2 days or longer at 30°C, score the plates. Mut⁺ transformants grow well on both MD and MM plates. Mut⁵ transformants grow well only on MD plates, and show little or no growth on the MM plates.



We recommend that you purify your His⁺ transformants to ensure pure clonal isolates. You may do this before or after testing for the Mut phenotype.

Screening for Mut* and MutS Transformants, continued

Replica-Plating Procedure

The following procedure gives a lower rate of misclassifications, but it increases the overall Mut⁺/Mut^S screening procedure by 2 days. For the procedure, you need replica-plating equipment.

- 1. Using sterile toothpicks, patch 100 His⁺ transformant on MD plates (2–3 plates). For controls, make one patch from each of the strains GS115/His⁺ Mut^s Albumin and GS115/His⁺ Mut⁺ β -gal onto the MD plates.
- 2. Incubate the plates at 28–30°C for 2 days.
- 3. After 2 days, replica-plate the patches from the MD plates onto fresh MM and MD plates to screen for Mut^S transformants.
- 4. Incubate the replica plates at 28–30°C for 2 days.
- 5. After 2 days at 28–30°C, score the replica plates. Look for patches that grow normally on the MD replica plates, but show little or no growth on the MM replica plates. Including His⁺ Mut⁺ and His⁺ Mut^S control patches on each plate will provide examples of Mut⁺ and Mut^S phenotypes.

Screening by Functional Assay

Some researchers have used a functional assay to directly screen for high expressing *Pichia* recombinant clones without first screening for Mut⁵ or Mut⁺ phenotypes. If you elect to screen directly for high-expressing recombinants, be sure to also check the Mut phenotype to help you optimize expression of your recombinant clone.

Multiple Integration Events

Pichia pastoris is capable of integrating multiple copies of transforming DNA via recombination into the genome at sites of sequence homology (see page 72 for figure). Although the exact mechanism of multiple integration events is not fully understood, such events are reasonably common among selected transformants (in this case, His⁺ transformants).

Successful expression of the gene of interest to useful levels may depend upon the generation of a recombinant strain that contains multiple copies integrated at the *AOX1* or *HIS4* loci. In addition to simply screening expression levels among several His⁺ Mut^S or His⁺ Mut⁺ recombinants via SDS-PAGE analysis, it may be desirable to determine the existence of strains that have multiple integrants in the His⁺ Mut^S or His⁺ Mut⁺ recombinant strain.

Refer to the **Appendix**, page 81, for methods to detect multiple integration events.

Screening for Mut⁺ and Mut^S Transformants, continued

Vectors for Multiple Integration

Three vectors for isolating or generating *Pichia* recombinants containing multiple integrations of your desired gene are available separately from Invitrogen (see page 88 for ordering information). Two of the vectors, pPIC3.5K and pPIC9K, are used *in vivo* to identify possible transformants with multiple copies of your gene. The other vector, pAO815, is used for creating tandem copies of your gene *in vitro* before transforming into *Pichia*.

pPIC3.5K and pPIC9K


The vectors pPIC3.5K (for intracellular expression) and pPIC9K (for secreted expression) contain the bacterial kanamycin resistance gene cloned between the *HIS4* gene and the 3′ *AOX1* region. The kanamycin resistance gene confers resistance to G418 in *Pichia*. Multiple insertions of the kanamycin gene into the *Pichia* chromosome increase the resistance to G418. Because the kanamycin resistance gene is linked to your gene, isolation of hyper-resistant G418 transformants may also indicate that your gene is present in multiple copies.

pAO815

The vector pAO815 (for intracellular expression) allows you to generate multiple copies of your gene *in vitro* by creating an expression cassette (e.g., P_{AOXI}–your gene of interest–HIS4), and cloning multiple copies in tandem in the vector. The vector is then transformed into *Pichia* and transformants are selected and tested for increased expression of the desired protein.

Screening for Mut⁺ and Mut^S Transformants, continued

Scoring Templates



PCR Analysis of Pichia Integrants

Introduction

Use the following protocol to analyze *Pichia* integrants to determine if the gene of interest has integrated into the *Pichia* genome. Isolate genomic DNA from 6 to 10 Mut^S or Mut⁺ *Pichia* clones and the strain transformed with the parent plasmid using the protocol on page 79. After isolating your DNA, use the procedure below to identify integrants. Amplify the gene of interest using the α -factor primer (for pPIC9 only) or the 5′ *AOX1* primer paired with the 3′ *AOX1* primer included in the kit. This protocol can confirm the integration of the gene of interest but does not provide information on the site of integration.

Note: A more direct method for PCR screening is available in the Appendix, page 78.

Isolating Genomic DNA from *Pichia*

Invitrogen's Easy-DNATM Kit provides a fast and easy method to isolate genomic DNA from *Pichia pastoris*. See page 88 for ordering information.

Analysis by PCR

1. Set up PCR reactions as follows:

 $\begin{array}{lll} 10 \text{X PCR Buffer} & 5 \ \mu\text{L} \\ \text{Genomic DNA (\sim1 μg)} & 5 \ \mu\text{L} \\ 100 \ \text{mM dNTPs ($25 \ \text{mM each})} & 1 \ \mu\text{L} \\ 5' \ AOX1 \ \text{Primer ($0.1 \ \mu\text{g}/\mu\text{L})} & 5 \ \mu\text{L*} \\ 3' \ AOX1 \ \text{Primer ($0.1 \ \mu\text{g}/\mu\text{L})} & 5 \ \mu\text{L*} \\ \text{Sterile water} & \text{to } 50 \ \mu\text{L} \\ Taq \ \text{Polymerase ($5 \ \text{U}/\mu\text{L})} & 0.25 \ \mu\text{L} \end{array}$

For amplification controls, use 100 ng of recombinant plasmid (positive control) and 100 ng of the appropriate plasmid without insert (negative control).

2. Load thermocycler and run the following program:

Step	Temperature	Time	Cycle
Hot Start	94°C	2 minutes	1X
Denaturation	94°C	1 minute	
Annealing	55°C	1 minute	25–35X
Extension	72°C	1 minute	
Final Extension	72°C	7 minutes	1X

3. Analyze 10 µL on a 1X TAE, 0.8 % agarose gel.

^{*}Resuspend the primers in 20 μ L sterile water to prepare a 0.1 μ g/ μ L solution. You may decrease the amount of primer. For ~20 pmoles primer, use 2 μ L of each resuspended primer.

PCR Analysis of Pichia Integrants, continued

Interpreting PCR

If screening Mut⁺ integrants, you should see two bands, one corresponding to the size of your gene of interest, the other to the *AOX1* gene (approximately 2.2 kb). If screening Mut^S integrants in GS115, you should see only the band that corresponds to the gene of interest. In KM71, because of the *ARG4* insert in *AOX1*, the PCR product is 3.6 kb. Parent plasmids produce the following sized PCR products. Add these fragments to the size of your insert to interpret your PCR results.

Vector	PCR Product
pHIL-D2	188 bp
pPIC3.5	214 bp
pHIL-S1	262 bp
pPIC9 (using the 5' AOX1 primer)	492 bp
pPIC9 (using the α-Factor primer)	195 bp




If you use the α -factor primer as a PCR primer, you will not see a band with GS115 or KM71. This is because there is no α -factor signal associated with the chromosomal *AOX1* gene.

Sometimes there will be ghost bands appearing in your PCR. These do not seem to be significant as they have not been shown to be a problem.

Example of PCR Analysis

The figure below shows the results of a typical PCR analysis using the procedure on the previous page. Genomic DNA was isolated from *Pichia* recombinants and from appropriate controls. Ten microliter samples from each PCR were run on a 0.8% agarose gel. Lanes 1 and 8 contain markers for a 1 kb ladder; Lanes 2–4 are *Pichia* recombinants; Lane 5 is pHIL-D2 with the gene of interest; Lane 6 is GS115/pHIL-D2 (no insert); and Lane 7 is pHIL-D2 alone.



Discussion

Lane 7 shows the 188 bp PCR product made from pHIL-D2 by priming with the 5′ and 3′ AOX1 primers (see page 26). Lane 6 shows the 188 bp product and the wild-type AOX1 gene (2.2 kb) from GS115/pHIL-D2. Lane 5 shows the expected size of our gene of interest cloned into pHIL-D2 (650 bp + 188 bp = 838 bp). Analysis of the *Pichia* recombinants in lanes 2–4 reveal that lanes 2 and 3 contain the insert, and that the recombinant in lane 2 may be a Mut^S because there is no wild-type AOX1. Lane 4, although from a His⁺ transformant, does not contain the gene of interest.

Expression of Recombinant Pichia Strains

Introduction

The purpose of this section is to determine the optimal method and conditions for expressing your gene. Below are some factors and guidelines that need to be considered before starting expression in *Pichia pastoris*. As with any expression system, optimal expression conditions are dependent on the characteristics of the protein being expressed.

Media

You need BMGY/BMMY (buffered complex glycerol or methanol medium), BMG/BMM (buffered minimal glycerol or methanol medium), or MGY/MM (minimal glycerol or minimal methanol medium) for expression (see **Appendix**, pages 63–65). BMG, BMM, BMGY, and BMMY are usually used for the expression of secreted proteins, particularly if pH is important for the activity of your protein. Unlike MGY and MM, they are all buffered media. Because these media are buffered with phosphate buffer, you may use a wide range of pH values to optimize production of your protein. BMGY/BMMY contain yeast extract and peptone which may help stabilize secreted proteins and prevent or decrease proteolysis of secreted proteins. Inclusion of yeast extract and peptone act as a "mixed feed" allowing better growth and biomass accumulation.

Proteases

There are some proteins specifically susceptible to proteases that have optimal activity at neutral pH. If this is the case, expression using MGY and MM media may be indicated. As *Pichia* expression progresses in an unbuffered medium such as MM, the pH drops to 3 or below, inactivating many neutral pH proteases (Brierley *et al.*, 1994). *Pichia* is resistant to low pH, so the low pH will not affect growth. In contrast, it has been reported that by including 1% Casamino acids (Difco) and buffering the medium at pH 6.0, extracellular proteases were inhibited, increasing the yield of mouse epidermal growth factor (Clare *et al.*, 1991b).

If you know your protein of interest is especially susceptible to neutral pH proteases, you may want to do your expressions in an unbuffered medium (MM). If there is no evidence that your secreted protein of interest is susceptible to proteases at neutral pH, we recommend you do your initial expressions in BMMY. If the expressed protein is degraded, you may then try expression in an unbuffered medium.

If the above options fail to protect your protein from degradation, you may wish to transform your gene into SMD1168. This strain genotype is *his4 pep4* and is deficient in proteinase A activity. You can use the same procedures for transforming and expressing GS115 with SMD1168. You can also use this strain in large-scale fermentation.

Aeration

The most important parameter for efficient expression in *Pichia* is adequate aeration during methanol induction. As a general rule when inducing expression, never allow cultures to be more than 10–30% of your total flask volume. We strongly recommend that you use baffled flasks. See **Appendix**, page 66 for suppliers of baffled flasks. Cover the flasks with cheesecloth (2–3 layers) or another loose fitting cover. Never use tight fitting covers. (Aeration is not as critical when generating biomass before induction.)

Expression of Recombinant Pichia Strains, continued

Kinetics of Growth

Note that while Mut⁺ and Mut^S strains grow at essentially the same rate in YPD or glycerol media, Mut⁺ grows faster than Mut^S when both are grown on methanol because of the presence of the *AOX1* gene product.

Temperature and Shaking

Perform expression at 30°C, in a shaking incubator. It is critical that the temperature does not exceed 30°C. If your incubator temperature fluctuates, set the temperature at 28°C. If using a floor shaking incubator, shake at 225–250 rpm. If using a table-top shaker that sits inside an incubator, shake at 250–300 rpm.

Before Starting

After you have verified recombinants in GS115 or KM71 as well as a control recombinant of GS115 or KM71/Vector (no insert), proceed with expression. When performing your expression, it is important to run the proper controls so that you will be able to interpret your expression results. You should use the following expression controls:

GS115/His⁺ Mut^S albumin Mut^S- Secretion control GS115/His⁺ Mut⁺ β-Gal Mut⁺- Intracellular control GS115 or KM71/Vector (no insert) Background control

Recombination can occur in many different ways that can affect expression. We recommend that you screen 6–10 verified recombinant clones for expression levels. Start with colonies from the freshest plates available. Colony viability drops over time, so if you have any doubts, it is better to streak out your strain. (You may also start the cultures with a small sample from a frozen glycerol stock that was generated from a single colony.)

Guidelines for Expression

The following information is designed to get you started with expression. You may have to change the conditions to optimize expression for your particular protein. Use bottom or side baffled flasks whenever possible (see page 66). If you are analyzing a number of recombinants, you can try 50 mL conical tubes. Make sure that the medium is well-aerated by increasing the rate of shaking or placing the tubes at an angle in the shaker.

Expression of Recombinant Pichia Strains, continued

Mut⁺ Intracellular or Secreted

Test the effectiveness of your expression conditions by growing GS115 β -Gal, which is Mut⁺ and expresses β -galactosidase intracellularly. Include GS115 or KM71 transformed with the parent vector as a control for background intracellular expression.

- 1. Using a single colony, inoculate 25 mL of MGY, BMG, or BMGY medium in a 250 mL baffled flask. Grow at 28–30°C in a shaking incubator (250–300 rpm) until the culture reaches an OD_{600} = 2–6 (log-phase growth, approximately 16–18 hours).
- 2. Harvest the cells by centrifuging at 1,500–3,000 \times g for 5 minutes at room temperature. Decant the supernatant and resuspend the cell pellet to an OD₆₀₀ of 1.0 in MM, BMM, or BMMY medium (approximately 100–200 mL) to induce expression.
- 3. Place the culture in a 1 liter baffled flask. Cover the flask with 2 layers of sterile gauze or cheesecloth, and return it to the incubator to continue growth.
- 4. Add 100% methanol to a final concentration of 0.5% methanol every 24 hours to maintain induction. Be sure to check the volume of the culture, and add methanol accordingly. Evaporation may reduce the culture volume.
- 5. At each of the times indicated below, transfer 1 mL of the expression culture to a 1.5 mL microcentrifuge tube. These samples will be used to analyze expression levels and to determine the optimal time post-induction to harvest. Centrifuge the samples at maximum speed in a tabletop microcentrifuge at room temperature for 2–3 minutes.

 Time points (hours): 0, 6, 12, 24 (1 day), 36, 48 (2 days), 60, 72 (3 days), 84, and 96 (4 days).
- 6. For secreted expression, transfer the supernatant to a separate tube. Freeze the supernatant and the cell pellets quickly in liquid N_2 or a dry ice/alcohol bath, and store them at -80° C until ready to assay.. For intracellular expression, decant the supernatant and store just the cell pellets at -80° C until ready to assay. Freeze quickly in liquid N_2 or a dry ice/alcohol bath.
- 7. Analyze the supernatants and the cell pellets for protein expression by Coomassie-stained SDS-PAGE and western blot, or functional assay (see **Analysis by SDS-Polyacrylamide Gel Electrophoresis**, page 49).

Expression of Recombinant Pichia Strains, continued

Mut^S Intracellular or Secreted

You can test the effectiveness of your expression conditions by growing GS115, which is Mut^S and secretes albumin to the medium. Remember to include GS115 or KM71 transformed with the parent vector as a control for background intracellular expression.

- 1. Using a single colony, inoculate 100 mL of MGY, BMG, or BMGY medium in a 1 liter baffled flask. Grow the culture at $28-30^{\circ}$ C in a shaking incubator (250–300 rpm) until the culture reaches an $OD_{600} = 2-6$ (approximately 16–18 hours.).
- 2. Harvest the cells by centrifuging at $1,500-3,000 \times g$ for 5 minutes at room temperature. To induce expression, decant the supernatant and resuspend the cell pellet in MM, BMM, or BMMY medium using 1/5 to 1/10 of the original culture volume (approximately 10-20 mL).
- 3. Place the cells in a 100 mL baffled flask. Cover the flask with 2 layers of sterile gauze or cheesecloth, and return them to the incubator to continue to grow.
- 4. Add 100% methanol to a final concentration of 0.5% every 24 hours to maintain induction.
- 5. At each of the times indicated below, transfer 1 mL of the expression culture to a 1.5 mL microcentrifuge tube. These samples will be used to analyze expression levels and determine the optimal time post-induction to harvest. Centrifuge at maximum speed in a tabletop microcentrifuge for 2–3 minutes at room temperature.
 - Time points (hours): 0, 24 (1 day), 48 (2 days), 72 (3 days), 96 (4 days), 120 (5 days), and 144 (6 days).
- 6. For secreted expression, **transfer the supernatant to a separate tube.** Store the supernatant and the cell pellets at –80°C until ready to assay. Freeze quickly in liquid N₂ or a dry ice/alcohol bath.
 - For intracellular expression, decant the supernatant and store just the cell pellets at -80° C until ready to assay. Freeze quickly in liquid N_2 or a dry ice/alcohol bath.
- 7. Analyze the cell pellets for protein expression by Coomassie-stained SDS-PAGE and western blot or functional assay (see **Analyzing Samples by SDS-Polyacrylamide Gel Electrophoresis**, next page).

Analyzing Samples by SDS-Polyacrylamide Gel Electrophoresis

Polyacrylamide Gel Electrophoresis

Invitrogen offers a wide range of pre-cast NuPAGE® and Tris-Glycine polyacrylamide gels and electrophoresis apparatus. The patented NuPAGE® Gel System avoids the protein modifications associated with LaemmLi-type SDS-PAGE, ensuring optimal separation for protein analysis. In addition, Invitrogen carries a large selection of molecular weight protein standards and staining kits. For more information about the appropriate gels, standards, and stains to use, refer to www.invitrogen.com or contact Technical Support (page 89).

If you are pouring your own gels, note that any standard SDS-polyacrylamide gel apparatus and protocol will work, for example, a 12% polyacrylamide gel with a 5% stacking gel is recommended for proteins ranging in size from 40–100 kDa. For other recommendations, see standard texts such as *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994) or *Molecular Biology: A Laboratory Manual* (Sambrook *et al.*, 1989), *Guide to Protein Purification* (Deutscher, 1990), or *Protein Methods* (Bollag *et al.*, 1996).

Preparing Samples

You will need to prepare Breaking Buffer (see page 66) and have acid-washed 0.5 mm glass beads on hand.

Preparing of cell pellets (Intracellular and Secreted Expression):

- 1. Thaw cell pellets quickly and place on ice.
- 2. For each 1 mL sample, add 100 μ L Breaking Buffer to the cell pellet and resuspend.
- 1. Add an equal volume of acid-washed glass beads (size 0.5 mm). Estimate equal volume by displacement.
- 4. Vortex 30 seconds, then incubate on ice for 30 seconds. Repeat for a total of 8 cycles.
- 5. Centrifuge at maximum speed for 10 minutes at 4°C. Transfer the clear supernatant to a fresh microcentrifuge tube.
- 6. Take 50 μL of supernatant (cell lysate) and mix with an appropriate volume of denaturing PAGE Gel Loading buffer (Sample Buffer).
- 7. Heat the sample as recommended and load 10– $20\,\mu L$ per well. Thickness of the gel and number of wells will determine volume loaded. You may store the remaining sample at $-20\,^{\circ}\text{C}$ for western blots, if necessary. You may store the cell lysates at $-80\,^{\circ}\text{C}$ for further analysis.

Analyzing Samples by SDS-Polyacrylamide Gel Electrophoresis, continued

Preparing Samples, continued

Preparing of supernatant (Secreted Expression only):

- 1. Thaw supernatants and place on ice.
- 2. Mix $50~\mu L$ of the supernatant with an appropriate volume of denaturing PAGE Gel Loading buffer.
- 3. Heat the sample as recommended; then load 10– $30~\mu L$ onto the gel. Remaining sample may be stored at -20°C for western blots, if necessary. Supernatants may be stored at -80°C for further analysis.
- 4. Analyze Coomassie-stained gel and western blot (if necessary) for your protein.
- 5. If no protein is seen by Coomassie or by western blot, then concentrate the supernatant 5–10 fold and analyze samples again by western blot. Centricon and Centriprep filters (Millipore) are very useful for this purpose.

Protein Concentration

Lowry, BCA (Pierce) or Bradford protein determinations can be performed to quantify the amounts of protein in the cell lysates and medium supernatants. In general, *Pichia* cell lysates contain 5–10 µg/µL protein. *Pichia* medium supernatants will vary in protein concentration primarily due to the amount of your secreted protein. *Pichia* secretes very few native proteins. If the protein concentration of the medium is >50 µg/mL, 10 µL of medium will give a faint band on a Coomassie-stained SDS-PAGE gel.

Analyzing Samples by SDS-Polyacrylamide Gel Electrophoresis, continued

Controls

Include the following samples as controls on your SDS-PAGE:

- Molecular weight standards appropriate for your desired protein
- A sample of your protein as a standard (if available)
- A sample of GS115 or KM71 with the parent plasmid transformed into it. This control shows the background of native *Pichia* proteins that are present intracellularly. Inclusion of this sample will help you differentiate your protein from background if you express it intracellularly.
- Analyze the GS115 β -Gal and Albumin controls also as they should indicate any problems with the media or expression conditions



In addition to Coomassie-stained SDS-PAGE, we strongly recommend performing a western blot or another more sensitive assay to detect your protein. Visualization of the expressed protein depends on several factors, including its expression level, its solubility, its molecular weight, and whether an abundant cellular protein of the same size masks it. Western blot analysis, enzymatic activities, or a defined purification profile, if available, may help to identify the expressed protein among the native *Pichia* cellular proteins.

Analyzing Protein Expression

Inspection of your Coomassie-stained SDS-PAGE should reveal the induction over time of your protein co-migrating with your standard. If there is no recombinant protein visible, then perform a western blot or a functional assay if you have one.

If you detect low expression of your recombinant protein, see **Optimizing** *Pichia* **Protein Expression**, page 52.

Test your expression conditions with the one of the two control strains included in the kit (GS115 β -Gal or Albumin).

If there is no indication of expression at all, perform a northern analysis to see if and how much full-length mRNA is induced. See page 84 for an RNA isolation protocol.

Optimizing Pichia Protein Expression

Introduction

Based on available data, there is approximately a 75% chance of expressing your protein of interest in *Pichia pastoris* at reasonable levels. The biggest hurdle seems to be generating initial success—*i.e.*, expression of your protein at **any** level. While there are relatively few examples of expression of ≥ 10 grams/liter, there are many examples of expression in the ≥ 1 gram/liter range, making the *Pichia pastoris* expression system one of the most productive eukaryotic expression systems available. Likewise, there are several examples of proteins that have been successfully expressed in *Pichia pastoris* that were completely unsuccessful in baculovirus or *Saccharomyces cerevisiae*, suggesting that the *Pichia pastoris* system is an important alternative to have available. If you obtain no or low protein expression in your initial expression experiment, use the following guidelines to optimize expression.

Proteolysis or Degradation

- Do a time course study of expression. Check to see if there is a time point that yields a larger percentage of full-length protein.
- If secreting your protein, check to see if your protein is susceptible to neutral pH proteases by expressing in unbuffered medium (MM). In addition, try 1% Casamino acids with buffered medium to inhibit extracellular proteases.
- Try using SMD1168 (proteinase A-deficient) for expression (see **Proteases**, page 45).

Low Secreted Expression Levels

- Check cell pellet to see if overall expression is low or if the protein did not secrete. If it did not secrete, try a different signal sequence (e.g., a native or α-factor signal sequence).
- Concentrate your supernatant by ammonium sulfate precipitation or ultrafiltration (see page 56).
- For Mut⁺, induce expression with a higher density culture.

Low Expression Levels

- Check both Mut⁺ and Mut^S recombinants for increased expression. Some proteins express better in one type of genetic background than another.
- If secreting your protein, try intracellular expression. The protein may not be processed correctly and fail to secrete. Be sure you check your cell pellets for evidence of expression. If you are having problems with intracellular expression, try secreting your protein. It probably will glycosylate, which may be desirable or not. If glycosylation is undesirable, oligosaccharides can be removed with Peptide:N-Glycosidase F (New England Biolabs, page 58).
- Scale up to fermentation (page 56). *Pichia* is a yeast and is particularly well suited to fermentation.

Optimizing Pichia Protein Expression, continued

No Expression

Be sure to try some of the easier things solutions on the previous page because very low expression may be perceived as no expression. If none of those improve protein expression, perform a northern blot analysis to check for transcription of your gene. There is a protocol in the **Appendix** for RNA isolation from *Pichia* (see page 84).

Analyze for the presence of your insert by PCR (see page 43). You can reasonably analyze 12–20 transformants by PCR. Remember to include the vector only and original (one copy) construct controls in order to analyze your PCR experiment.

If you see premature transcriptional termination, check the AT content of your gene. In *Saccharomyces*, there are a few consensus sequences that promote premature termination. One of these, TTTTTATA, resembles a sequence in HIV-1 gp120, ATTATTTTAT AAA, which when expressed in *Pichia* gave premature termination of the mRNA. When this sequence was changed, longer transcripts were found (Scorer *et al.*, 1993).

Hyperglycosylation

If your protein is hyperglycosylated:

- Try intracellular expression to bypass the secretion pathway and to avoid protein modification.
- Try deglycosylating the protein with Peptide:N-Glycosidase F or other enzymes (see page 57).
- Engineer your gene to remove N-linked glycosylation sites.

Scaling Up Expression

Guidelines for Expression

Once expression is optimized, scale-up your expression protocol to produce more protein by increasing the culture volume using larger baffled flasks (below) or by fermentation. Use the guidelines below to scale-up your expression protocol. To purify your protein, see the references listed on page 57.

Mut⁺ Intracellular or Secreted

- 1. Using a single colony, inoculate 25 mL of MGYH, BMGH, or BMGY in a 250 mL baffled flask. Grow at 28–30°C in a shaking incubator (250–300 rpm) until culture reaches an OD_{600} = 2–6 (approximately 16–18 hours).
- 2. Use this 25 mL culture to inoculate 1 liter of MGYH, BMGH, or BMGY in a 3 or 4 liter baffled flask and grow at $28-30^{\circ}$ C with vigorous shaking (250–300 rpm) until the culture reaches log phase growth (OD₆₀₀ = 2–6).
- 3. Harvest the cells using sterile centrifuge bottles by centrifuging at $1,500-3,000 \times g$ for 5 minutes at room temperature. To induce expression, decant the supernatant and resuspend cell pellet to an $OD_{600} = 1.0$ (2–6 liters) in MMH, BMMH, or BMMY medium to start induction.
- 4. Aliquot the culture between several 3 or 4 liter baffled flask. Cover the flasks with 2 layers of sterile gauze or cheesecloth and return to incubator. Continue to grow at 28–30°C with shaking.
- 5. Add 100% methanol to 0.5% every 24 hours until the optimal time of induction is reached as determined from the time course study.
- 6. Harvest cells by centrifuging at $1,500-3,000 \times g$ for 5 minutes at room temperature.
- 7. For intracellular expression, decant the supernatant and store the cell pellets at -80°C until ready to process.

For secreted expression, **save the supernatant, chill it to 4°C, and concentrate it, if desired** (see page 56). Proceed directly to purification (page 57) or store the supernatant at –80°C until ready to process further.

Scaling Up Expression, continued

Mut^S Intracellular or Secreted

- 1. Using a single colony, inoculate 10 mL of MGYH, BMGH, or BMGY in a 100 mL baffled flask. Grow at 28–30°C in a shaking incubator (250–300 rpm) until the culture reaches an $OD_{600} = 2-6$ (approximately 16–18 hours).
- 2. Use this 10 mL culture to inoculate 1 liter of MGYH, BMGH, or BMGY in a 3 or 4 liter baffled flask and grow at 28–30°C with vigorous shaking (250–300 rpm) until the culture reaches log phase growth (OD₆₀₀ = 2–6).
- 3. Harvest the cells by centrifuging at $1,500-3,000 \times g$ for 5 minutes at room temperature. To induce expression, decant the supernatant and resuspend cell pellet in 1/5 to 1/10 of the original culture volume of MMH, BMMH, or BMMY medium (approximately 100-200 mL).
- 4. Place the culture in a 1 liter baffled flask. Cover the flask with 2 layers of sterile gauze or cheesecloth and return to incubator. Continue to grow at 28–30°C with shaking.
- 5. Add 100% methanol to 0.5% every 24 hours until the optimal time of induction is reached.
- 6. Harvest cells by centrifuging at $1,500-3,000 \times g$ for 5 minutes at room temperature.
- 7. For intracellular expression, decant the supernatant and store the cell pellets at –80°C until ready to process.

For secreted expression, **save the supernatant, chill it to 4°C, and concentrate it, if desired** (see next page). Proceed directly to purification (page 57) or store the supernatant at –80°C until ready to process further.



To increase the amount of cells for Mut^S recombinants, increase the number of flasks, put 200–300 mL in a 3 liter flask, or try fermentation.

Scaling Up Expression, continued

Concentrating Proteins

Proteins secreted into the media are usually >50% homogeneous and require some additional purification. There are several general methods to concentrate proteins secreted from *Pichia*. These general methods include:

- Ammonium sulfate precipitation
- Dialysis
- Centrifuge concentrator for small volumes (e.g., Centricon or Centriprep devices available from Millipore)
- Pressurized cell concentrators for large volumes (e.g., Amicon ultrafiltration devices available from Millipore)
- Lyophilization

A general guide to protein techniques is *Protein Methods* (Bollag et al., 1996).

Cell Lysis

A general procedure for cell lysis using glass beads is provided on the next page. There is also a cell lysis protocol in *Current Protocols in Molecular Biology*, page 13.13.4. (Ausubel *et al.*, 1994) and in *Guide to Protein Purification* (Deutscher, 1990). We also recommend lysis by French Press (follow the manufacturer's suggestions for yeast).

Fermentation

Basic guidelines are available for fermentation of *Pichia* from Invitrogen. We recommend that only those with fermentation experience or those who have access to people with experience attempt fermentation. Contact Technical Support for more information (page 89).

Protein Purification and Glycosylation

Introduction

At this point, you have an optimized protocol for expressing your protein and a method to scale-up production of your protein for large-scale purification. You may already have a method to purify your protein. Since every protein is different, it is difficult to recommend specific techniques for purification. For an overview of methods for purification see (Deutscher, 1990) or (Ausubel *et al.*, 1994).

Some Protein Purification Techniques

Some techniques are listed below and are discussed thoroughly in *Guide to Protein Purification* (Deutscher, 1990). Be sure to perform all steps from cell lysis to purified protein at 4°C.

Ion-Exchange Chromatography Gel Filtration
Affinity Chromatography Chromatofocusing
Isoelectric Focusing Immunoprecipitation
Solubilization (Membrane Proteins) Lectin Affinity Chromatography

Procedure for Cell Lysis

Prepare Breaking Buffer (BB) as described in **Appendix**, page 66.

- 1. Wash cells once in BB by resuspending them and centrifuging 5–10 minutes at $3,000 \times g$ at $4^{\circ}C$.
- 2. Resuspend the cells to an OD_{600} of 50–100 in BB.
- 3. Add an equal volume of acid-washed glass beads (0.5 mm). Estimate volume by displacement.
- 4. Vortex the mixture 30 seconds, then incubate on ice for 30 seconds. Repeat 7 more times. Alternating vortexing with cooling keeps the cell extracts cold and reduces denaturation of your protein.
- 5. Centrifuge the sample at 4° C for 5–10 minutes at $12,000 \times g$.
- 6. Transfer the clear supernatant to a fresh container and analyze for your protein. The total protein concentration should be around 5–10 mg/mL.
- 7. Save the pellet and extract with 6 M urea or 1% Triton X-100 to check for insoluble protein.



Biospec (Bartlesville, OK) makes a Bead Beater $^{\text{\tiny TM}}$ that can handle 5–200 mL volumes of cell suspension.

Analysis of Glycoproteins

When expressing and purifying a glycosylated protein in a heterologous expression system, it is desirable to quickly determine whether the protein is glycosylated properly. For carbohydrate analysis of proteins to characterize glycosylated proteins, see (Ausubel *et al.*, 1994), Unit 17. Further information about glycosylation in eukaryotes is available in a review (Varki & Freeze, 1994).

Protein Purification and Glycosylation, continued

Enzymes for Analyzing Glycoproteins

These are just a few of the enzymes available for carbohydrate analysis. Abbreviations are as follows: Asn - Asparagine, Gal - Galactose, GlcNAc - N-acetylglucosamine, GalNAc - N-acetylgalactosamine, and NeuAc - N-acetylneuraminic acid.

Enzyme	Type of enzyme	Specificity
Endoglycosidase D	Endo	Cleaves various high mannose glycans
Endoglycosidase F	Endo	Cleaves various high mannose glycans
Endoglycosidase H	Endo	Cleaves various high mannose glycans
β-galactosidase	Exo	Removes terminal galactosides from Gal- β 1,3-GlcNAc, Gal- β 1,4-GlcNAc or Gal- β 1,3-GalNAc.
Peptide:N-Glycosidase F	Endo	Glycoproteins between Asn and GlcNAc (removes oligosaccharides)
Sialidases (Neuraminidases) Vibrio cholerae Clostridium perfringens Arthobacter ureafaciens Newcastle disease virus	Exo	NeuAc-α2,6-Gal, NeuAc-α2,6-GlcNAc or NeuAc-α2,3-Gal

Commercial Carbohydrate Analysis

There are a number of commercial vendors who will contract to analyze proteins for glycosylation. A number of companies also supply kits and reagents for researchers to do carbohydrate analysis in their own laboratories. A partial list is provided below:

Company	Type of Service	Contact
Glyko	Kits for Carbohydrate	1-800-334-5956
	Analysis	www.prozyme.com
	Reagents	
	Contract Services	
New England BioLabs	Reagents	1-800-632-5227
		www.neb.com

Appendix

E. coli Media Recipes

Introduction

You select transformants with pAO815, pPIC3K, or pPIC9K on LB agar containing 50 to $100 \mu g/mL$ ampicillin.

LB (Luria-Bertani) Medium

1% Tryptone 0.5% Yeast Extract 1% NaCl pH 7.0

1. For 1 liter, dissolve the following in 950 mL deionized water

10 g tryptone 5 g yeast extract 10 g NaCl

- 2. Adjust the pH of the solution to 7.0 with NaOH and bring the volume up to 1 liter.
- 3. Autoclave the medium for 20 minutes at 15 lbs/sq. in. Let it cool to ~55°C and add desired antibiotics at this point.
- 4. Store the medium at room temperature or at 4°C.

LB agar plates

- 1. Make LB Medium as described above and add to it 15 g/liter agar before autoclaving.
- 2. Autoclave the medium for 20 minutes at 15 lbs/sq. in.
- 3. Let the medium cool to \sim 55°C and add the desired antibiotics. Pour into 10 cm Petri plates. Let the plates harden, invert, and store at 4°C.

Pichia Media Recipes

Introduction

Expressing recombinant proteins in *Pichia pastoris* requires the preparation of several different media. Recipes for these media are included in this section. In addition, Yeast Nitrogen Base is available from Invitrogen (see below for ordering information).

Item	Amount	Cat. no.
Yeast Nitrogen Base	67 g pouch	Q300-07
-with ammonium sulfate -without amino acids	Each pouch contains reagents to prepare 500 mL of a 10X YNB solution	
	500 g	Q300-09

Stock Solutions

10X YNB (13.4% Yeast Nitrogen Base with Ammonium Sulfate without amino acids)

- 1. Dissolve 134 g of yeast nitrogen base (YNB) with ammonium sulfate and without amino acids in 1,000 mL of water. Heat the solution to dissolve YNB completely in water.
- 2. Alternatively, use 34 g of YNB without ammonium sulfate and amino acids and 100 g of ammonium sulfate.
- 3. Filter sterilize and store at 4°C. The shelf life of this solution is approximately one year. If you are using the YNB pouch included in the kit, follow the directions on the pouch.

Note: *Pichia* cells exhibit optimal growth with higher YNB concentrations, therefore, the amount of YNB used in this kit is twice as concentrated as YNB formulations for *Saccharomyces*.

500X B (0.02% Biotin)

- 1. Dissolve 20 mg biotin in 100 mL of water and filter sterilize.
- 2. Store at 4°C. The shelf life of this solution is approximately one year.

100X H (0.4% Histidine)

- 1. Dissolve 400 mg of L-histidine in 100 mL of water. Heat the solution, if necessary, to no greater than 50°C in order to dissolve.
- 2. Filter sterilize and store at 4°C. The shelf life of this solution is approximately one year.

10X D (20% Dextrose)

- 1. Dissolve 200 g of D-glucose in 1,000 mL of water.
- 2. Autoclave for 15 minutes or filter sterilize. The shelf life of this solution is approximately one year.

Stock Solutions, continued

10X M (5% Methanol)

- 1. Mix 5 mL of methanol with 95 mL of water.
- 2. Filter sterilize and store at 4°C. The shelf life of this solution is approximately two months.

10X GY (10% Glycerol)

- 1. Mix 100 mL of glycerol with 900 mL of water.
- 2. Sterilize either by filtering or autoclaving. Store at room temperature. The shelf life of this solution is greater than one year.

100X AA (0.5% of each Amino Acid)

- 1. Dissolve 500 mg each of L-glutamic acid, L-methionine, L-lysine, L-leucine, and L-isoleucine in 100 mL of water.
- 2. Filter sterilize and store at 4°C. The shelf life of this solution is approximately one year.

1 M potassium phosphate buffer, pH 6.0:

- 1. Combine 132 mL of 1 M K_2HPO_4 , 868 mL of 1 M KH_2PO_4 and confirm that the pH = 6.0 ± 0.1 (if the pH needs to be adjusted, use phosphoric acid or KOH).
- 2. Sterilize by autoclaving and store at room temperature. The shelf life of this solution is greater than one year.

YPD or YEPD

Yeast Extract Peptone Dextrose Medium (1 liter)

1% yeast extract

2% peptone

2% dextrose (glucose)

Note: If you are using the YP Base Medium or the YP Base Agar medium pouches included with the Multi-Copy *Pichia* Expression Kit, follow the directions on the pouch.

- 1. Dissolve 10 g yeast extract and 20 g of peptone in 900 mL of water. **Note**: Add 20 g of agar if making YPD slants or plates.
- 2. Autoclave for 20 minutes on liquid cycle.
- 3. Add 100 mL of 10X D.
- 4. Store the liquid medium at room temperature. Store the YPD slants or plates at 4°C. The shelf life is several months.

YPD-Geneticin[®] plates

Yeast Extract Peptone Dextrose Medium

1% yeast extract

2% peptone

2% dextrose (glucose)

2% agar

Variable amounts of Geneticin® (see page 88 for ordering information)

Use 50 mg/mL Geneticin® stock solution to make YPD plates containing Geneticin® at final concentrations of 0.25, 0.5, 0.75, 1.0, 1.5, 1.75, 2.0, 3.0, and 4.0 mg/mL.

For 250 mL (8 to 10 plates of a single Geneticin® concentration):

- 1. Combine 2.5 g yeast extract, 5 g peptone, and 5 g agar in 225 mL deionized water.
- 2. Autoclave for 20 minutes on liquid cycle.
- 3. Add 25 mL of 10X D and mix well.
- 4. Cool YPD to approximately 55–60°C and add appropriate volume of Geneticin® stock (see chart below). Remember to also make several YPD plates **without** Geneticin®.
- 5. Mix well by swirling, but take care to minimize bubble formation.
- 6. Pour agar solution into 10 cm Petri plates. Let plates harden, invert, and store bagged at 4°C. Plates are stable for at least 6 months.

Final [Geneticin®] (mg/mL)	mL Geneticin® stock per 250 mL YPD
0.25	1.25
0.50	2.5
0.75	3.75
1.00	5.0
1.50	7.5
1.75	8.75
2.00	10.0
3.00	15.0
4.00	20.0

MGY and MGYH

<u>Minimal Glycerol Medium + Histidine</u> (1 liter)

1.34% YNB 1% glycerol 4×10^{-5} % biotin ± 0.004 % histidine

- 1. Combine aseptically 800 mL autoclaved water with 100 mL of 10X YNB, 2 mL of 500X B, and 100 mL of 10X GY.
- 2. For growth of his4 strains in this medium, a version can be made that contains histidine (called MGYH) by adding 10 mL of 100X H stock solution.
- 3. Store at 4°C. The shelf life of this solution is approximately two months.

RD and RDH Liquid Media

Regeneration Dextrose Medium + Histidine (1 liter)

1 M sorbitol 2% dextrose 1.34% YNB 4×10^{-5} % biotin 0.005% amino acids ± 0.004 % histidine

- 1. Dissolve 186 g of sorbitol in 700 mL of water and proceed to Step 2.
- 2. Autoclave 20 minutes on liquid cycle.
- 3. Cool and maintain the liquid medium in a 45°C water bath.
- 4. Prepare a prewarmed (45°C) mixture of the following stock solutions:

100 mL of 10X D

100 mL of 10X YNB

2 mL of 500X B

10 mL of 100X AA

88 mL of sterile water

Add to sorbitol solution.

5. For growth of *his4* strains you must add histidine to the media. Add 10 mL of 100X H (histidine) to the prewarmed mixture in Step 4. Store liquid medium at 4°C. Media should last for several months.

RDB and RDHB Agar Plates

- 1. Dissolve 186 g of sorbitol in 700 mL of water and add 20 g of agar.
- 2. Autoclave 20 minutes on liquid cycle.
- 3. Place the autoclaved solution in a 60°C water bath prior to addition of prewarmed mixture of stock solutions. This will keep the medium from becoming too thick to mix reagents.
- 4. Prepare the prewarmed (45°C) mixture from **RD** and **RDH** Liquid Media, Step 4, above. Add to sorbitol/agar solution. If you are selecting for His⁺ transformants, do not add histidine.
- 5. Pour the plates immediately after mixing the solutions in Step 4. Store the plates at 4°C. Shelf life is several months.

RD and RDH Top Agar

- 1. Dissolve 186 g of sorbitol in 700 mL of water and add 10 g of agar or agarose.
- 2. Autoclave 20 minutes on liquid cycle.
- 3. Place the autoclaved solution in a 60°C water bath prior to addition of prewarmed mixture of stock solutions. This will keep the medium from becoming too thick to mix reagents.
- 4. Prepare the prewarmed (45°C) mixture from **RD** and **RDH** Liquid Media, Step 4, previous page. Add to sorbitol/agar solution. If you are selecting for His⁺ transformants, do not add histidine.
- 5. Place the solution to 45°C after adding the solutions in Step 4. During transformation, use as a molten solution at 45°C.
- 6. Store top agar at 4°C. Shelf life is several months.

MD and MDH

Minimal Dextrose Medium + Histidine (1 liter)

1.34% YNB

 4×10^{-5} % biotin

2% dextrose

- 1. For medium, autoclave 800 mL of water for 20 minutes on liquid cycle.
- 2. Cool to about 60°C and then add:

100 mL of 10X YNB

2 mL of 500X B

100 mL of 10X D

- 3. To make MDH, add 10 mL of 100X H stock solution. Mix and store at 4°C.
- 4. For plates, add 15 g agar to the water in Step 1 and proceed.
- 5. If preparing plates, pour the plates immediately. MD stores well for several months at 4°C.

MM and MMH

<u>Minimal Methanol + Histidine</u> (1 liter)

1.34% YNB

 $4 \times 10^{-5}\%$ biotin

0.5% methanol

- 1. For medium, autoclave 800 mL of water for 20 minutes on liquid cycle
- 2. Cool autoclaved water to 60°C and add:

100 mL of 10X YNB

2 mL of 500X B

100 mL of 10X M

- 3. To make MMH, add 10 mL of 100X H stock solution. Mix and store at 4°C.
- 4. For plates, add 15 g agar to the water in Step 1 and proceed.
- 5. After mixing, pour the plates immediately. MM and MMH stores well for several months at 4° C.

BMG and **BMM**

Buffered Minimal Glycerol Buffered Minimal Methanol (1 liter)

100 mM potassium phosphate, pH $6.0 \cdot 1.34\%$ YNB

 $4 \times 10^{-5}\%$ biotin

1% glycerol or 0.5% methanol

- 1. Autoclave 700 mL water for 20 minutes on liquid cycle.
- 2. Cool to room temperature, then add the following and mix well:

100 mL 1 M potassium phosphate buffer, pH 6.0

100 mL 10X YNB

2 mL 500X B

100 mL 10X GY

- 3. For BMM, add 100 mL 10X M instead of glycerol.
- 4. Store media at 4°C. The shelf life of this solution is approximately two months.

BMGY and BMMY

<u>Buffered Glycerol-complex Medium</u> <u>Buffered Methanol-complex Medium</u> (1 liter)

1% yeast extract

2% peptone

100 mM potassium phosphate, pH 6.0

1.34% YNB

 $4 \times 10^{-5}\%$ biotin

1% glycerol or 0.5% methanol

- 1. Dissolve 10 g of yeast extract, 20 g peptone in 700 mL water.
- 2. Autoclave 20 minutes on liquid cycle.
- 3. Cool to room temperature, then add the following and mix well:

100 mL 1 M potassium phosphate buffer, pH 6.0

100 mL 10X YNB

2 mL 500X B

100 mL 10X GY

- 4. For BMMY, add 100 mL 10X M instead of glycerol.
- 5. Store the media at 4°C. The shelf life of this solution is approximately two months.

Breaking Buffer

50 mM sodium phosphate, pH 7.4

1 mM PMSF (phenylmethylsulfonyl fluoride or other protease inhibitors) 1 mM EDTA

5% glycerol

- 1. Prepare a stock solution of your desired protease inhibitors and store it appropriately. Follow manufacturer's recommendations.
- 2. For 1 liter, dissolve in 900 mL deionized water: 6 g sodium phosphate (monobasic), 372 mg EDTA, and 50 mL glycerol
- 3. Use NaOH to adjust the pH of the solution and bring up the volume to 1 liter. Store the breaking buffer at 4°C.
- 4. Right before use, add the protease inhibitors.

Vendors for Baffled Flasks

Bellco (1-800-257-7043) has a wide variety of baffled flasks from 50 to 2,000 mL. Wheaton (1-609-825-1100) only sells side baffle flasks.

Proteins Expressed in *Pichia*

Table

The table below provides a partial list of references documenting successful expression of heterologous proteins in *Pichia pastoris*. Note that both Mut⁺ and Mut^S phenotypes were used successfully as well as secreted and intracellular expression.

Protein	Expression Levels grams/liter	Where Expressed How Expressed	Reference
Enzymes			
Invertase	2.3	Secreted Mut ⁺	(Tschopp <i>et al.,</i> 1987b)
Bovine Lysozyme c2	0.55	Secreted Mut ⁺	(Digan et al., 1989)
Streptokinase (active)	0.08	Intracellular	(Hagenson et al., 1989)
Alpha amylase	2.5	Secreted Mut ^S	(Paifer et al., 1994)
Pectate Lyase	0.004	Secreted Mut ^s	(Guo et al., 1995)
Spinach Phosphoribulokinase	0.1	Intracellular Mut ^s	(Brandes <i>et al.</i> , 1996)
Antigens			
Hepatitis B surface antigen	0.4	Intracellular Mut ^s	(Cregg et al., 1987)
Pertussis Antigen P69	3.0	Intracellular Mut ^s	(Romanos et al., 1991)
Tetanus Toxin Fragment C	12.0	Intracellular Mut ⁺ /Mut ^S	(Clare et al., 1991a)
HIV-1 gp120	1.25	Intracellular Mut ⁺	(Scorer et al., 1993)
Tick Anticoagulant protein	1.7	Secreted Mut ^S	(Laroche et al., 1994)
Bm86 Tick Gut Glycoprotein	1.5	Secreted *	(Rodriguez <i>et al.,</i> 1994)
Regulatory Proteins			
Tumor Necrosis Factor (TNF)	10.0	Intracellular Mut ^S	(Sreekrishna <i>et al.</i> , 1989)
Mouse Epidermal Growth Factor (EGF)	0.45	Secreted Mut ^s	(Clare et al., 1991b)
Human Interferon (IFN) α2b	0.4	Intracellular Mut ^S	(Garcia et al., 1995)

Proteins Expressed in *Pichia*, continued

Table, continued

Protein	Expression Levels grams/liter	Where Expressed How Expressed	Reference
Membrane Proteins			
Human CD38 (soluble portion)	0.05	Secreted Mut ^S	(Fryxell et al., 1995)
Mouse Serotonin Receptor	0.001	Secreted Mut ⁺	(Weiss et al., 1995)
Proteases and Protease	Inhibitors		
Carboxypeptidase B	0.8	Secreted Mut ⁺ /Mut ^s	(Despreaux & Manning, 1993)
Enterokinase	0.021	Secreted Mut ⁺	(Vozza et al., 1996)
Ghilanten	0.01	Secreted Mut ⁺	(Brankamp <i>et al.,</i> 1995)
Kunitz protease inhibitor	1.0	Secreted *	(Wagner et al., 1992)
Human Proteinase Inhibitor 6	0.05	Intracellular Mut ⁺	(Sun et al., 1995)
Antibodies			
Rabbit Single Chain Antibody	>0.1	Secreted Mut ^S	(Ridder et al., 1995)

^{*}Mut phenotype was not described in the paper.

Recombination and Integration in Pichia


Introduction

Like *Saccharomyces cerevisiae*, linear DNA can generate stable transformants of *Pichia pastoris* via homologous recombination between the transforming DNA and regions of homology within the genome (Cregg *et al.*, 1985; Cregg *et al.*, 1989). Such integrants show extreme stability in the absence of selective pressure even when present as multiple copies. The most commonly used expression vectors carry the *HIS4* gene for selection. These vectors are designed to be linearized with a restriction enzyme such that His⁺ recombinants are generated by recombination at the *AOX1* locus (see below) or at the *his4* locus (see next page). Note that single crossover events (insertions) are much more likely to happen than double crossover events (replacements). Multiple insertion events occur spontaneously at about 1–10% of the single insertion events.

Gene Insertion at AOX1 or aox1::ARG4

Gene insertion events at the *AOX1* (GS115) or *aox1::ARG4* (KM71) loci arise from a single crossover event between the loci and any of the three *AOX1* regions on the vector: the *AOX1* promoter, the *AOX1* transcription termination region (TT), or sequences even further downstream of *AOX1* (3´ *AOX1*). This results in the insertion of one or more copies of the vector upstream or downstream of the *AOX1* or the *aox1::ARG4* genes. The phenotype of such a transformant is His⁺ Mut⁺ (GS115) or His⁺ Mut^S (KM71). By linearizing the recombinant vector at a restriction enzyme site located in the 5´ or 3´ *AOX1* regions, Mut⁺ or Mut^S recombinants can be conveniently generated depending on the host strain used.


The figure below shows the result of an insertion of the plasmid 3′ to the intact AOX1 locus (Mut⁺) and the gain of P_{AOXI} , your gene of interest, and HIS4 (expression cassette). This event could also happen at the 5′ AOX1 regions of the plasmid and genome with the resulting insertion positioned 5′ to an intact AOX1 locus. This also occurs with non-linearized plasmid and plasmid that religates, although at a lower frequency.



Recombination and Integration in Pichia, continued

Gene Insertion Events at his4


In GS115 (Mut⁺) and KM71 (Mut⁵) gene insertion events at the *his4* locus arise from a single crossover event between the *his4* locus in the chromosome and the *HIS4* gene on the vector. This event results in the insertion of one or more copies of the vector at the *his4* locus. Because the genomic *AOX1* or *aox1::ARG4* loci are not involved in this recombination event, the phenotype of such a His⁺ transformant has the same Mut phenotype as the parent strain. **By linearizing the recombinant vector at a restriction enzyme site located in** *HIS4* **gene, Mut⁺ or Mut⁵ recombinants can be conveniently generated depending on the host strain used.** The figure below shows the result of an insertion of the plasmid between duplicated copies of the *HIS4/his4* genes, one still mutant, and the other wild-type.



Recombination and Integration in Pichia, continued

Multiple Gene Insertion Events

Multiple gene insertion events at a single locus in a cell do occur spontaneously with a low, but detectable frequency–between 1 and 10% of all selected His⁺ transformants. Multi-copy events can occur as gene insertions either at the *AOX1*, *aox1::ARG4*, or *his4* loci. This results in a Mut⁺ phenotype in GS115 and a Mut^S phenotype in KM71. Quantitative dot blot analysis, Southern blot analysis, and differential hybridization can detect multiple gene insertion events. See page 81 for a protocol to screen for multiple inserts.




Recombination and Integration in Pichia, continued

Gene Replacement at *AOX1* in GS115

In a *his4* strain such as GS115, a gene replacement (omega insertion) event arises from a double crossover event between the AOX1 promoter and $3^{\circ}AOX1$ regions of the vector and genome. This results in the complete removal of the AOX1 coding region (i.e., gene replacement). The resulting phenotype is His⁺ Mut⁵. His⁺ transformants can be readily and easily screened for their Mut phenotype, with Mut⁵ serving as a phenotypic indicator of integration via gene replacement at the AOX1 locus. The net result of this type of gene replacement is a loss of the AOX1 locus (Mut⁵) and the gain of an expression cassette containing P_{AOX1} , your gene of interest, and HIS4. The figure below shows a gene replacement event at the AOX1 locus.

Gene replacement (double-crossover event) is less likely to happen than insertions (single-crossover events). In general, we recommend linearizing your plasmid DNA to create *Pichia* recombinants by single-crossover events. By using GS115 or KM71, the Mut phenotype of the recombinant will be the same as the parent strain.



Electroporation of *Pichia*

Introduction

This method does not require the generation and maintenance of spheroplasts, making it a very convenient method for generating *Pichia* transformants. Efficiencies run about the same as spheroplasting (Scorer *et al.*, 1994).

Preparing Cells

- 1. Grow 5 mL of *Pichia pastoris* in YPD (page 61) in a 50 mL conical at 30°C overnight.
- 2. Inoculate 500 mL of fresh medium in a 2 liter flask with 0.1–0.5 mL of the overnight culture. Grow the culture overnight again to an $OD_{600} = 1.3-1.5$.
- 3. Centrifuge the cells at $1,500 \times g$ for 5 minutes at 4°C. Resuspend the pellet with 500 mL of ice-cold, sterile water.
- 4. Centrifuge the cells as in Step 3, and resuspend the pellet with 250 mL of ice-cold, sterile water.
- 5. Centrifuge the cells as in Step 3 and resuspend the pellet in 20 mL of ice-cold 1 M sorbitol.
- 6. Centrifuge the cells as in Step 3 and resuspend the pellet in 1 mL of ice-cold 1 M sorbitol for a final volume of approximately 1.5 mL.

Note: You may freeze the electrocompetent cells in $80 \, \mu L$ aliquots; however the transformation efficiencies will decrease significantly.

Transformation

- 1. Mix 80 μ L of the cells from Step 6 (above) with 5–20 μ g of linearized DNA (in 5–10 μ L TE Buffer) and transfer them to an ice-cold 0.2 cm electroporation cuvette.
- 2. Incubate the cuvette with the cells on ice for 5 minutes.
- 3. Pulse the cells according to the parameters for yeast (*Saccharomyces cerevisiae*) suggested by the manufacturer of the specific electroporation device being used.
- 4. Immediately add 1 mL of ice-cold 1 M sorbitol to the cuvette. Transfer the cuvette contents to a sterile microcentrifuge tube.
- 5. Spread 200–600 µL of aliquots on MD or RDB plates.
- 6. Incubate the plates at 30°C until colonies appear. Screen for Mut⁺/Mut^S phenotypes as indicated on page 38.

PEG 1000 Transformation Method for Pichia

Introduction

PEG procedure is usually better than LiCl, but not as good as spheroplasting or electroporation for transformation. However, it is convenient for people who do not have an electroporation device. The transformation efficiency is 10^2 to 10^3 transformants per μg of DNA.

Required Solutions

- Buffer A: 1.0 M Sorbitol (Fisher), 10 mM Bicine, pH 8.35 (Sigma), 3% (v/v) ethylene glycol (Merck)
- Buffer B: 40% (w/v) Polyethylene glycol 1000 (Sigma), 0.2 M Bicine, pH 8.35
- Buffer C: 0.15 M NaCl, 10 mM Bicine, pH 8.35

Filter sterilize the above solutions and store them at -20°C.

• Fresh, reagent grade DMSO that is from an unopened bottle or made fresh and stored at –70°C until use.



Cell competence decreases very rapidly after the cells thaw even when held on ice. It is critical to add DNA to frozen cell samples. To perform multiple transformations, it is recommended to process them in groups of six at a time.

Preparing Competent Cells

- 1. Streak *Pichia pastoris* strain for single colonies on a YPD plate and incubate the plate at 30°C for two days.
- 2. Inoculate a 10 mL YPD culture with a single colony from the plate and grow the culture overnight at 30°C with shaking.
- 3. In the morning, use an aliquot of the overnight culture to inoculate a 100 mL YPD culture to a starting OD_{600} of 0.1 and grow at 30°C to an OD_{600} of 0.5 to 0.8.
- 4. Harvest the culture by centrifugation at $3,000 \times g$ at room temperature and wash cells once in 50 mL of Buffer A.
- 5. Resuspend the cells in 4 mL of Buffer A and distribute them in 0.2 mL aliquots to sterile 1.5 mL microcentrifuge tubes. Add 11 μ L of DMSO to each tube, mix, and quickly freeze cells in a bath of liquid nitrogen.
- 6. Store the frozen tubes at -70° C.

PEG 1000 Transformation Method for Pichia, continued

Transformation

- 1. Use up to 50 μ g of each DNA sample in no more than 20 μ L total volume. Add the DNA directly to a still-frozen tube of competent cells. Include carrier DNA (40 μ g of denatured and sonicated salmon sperm DNA) with < 1 μ g DNA samples for maximum transformation frequencies.
- 2. Incubate all sample tubes in a 37°C water bath for five minutes. Mix the samples once or twice during this incubation period.
- 3. Remove the tubes from the bath and add 1.5 mL of Buffer B to each. Mix their contents thoroughly.
- 4. Incubate the tubes in a 30°C water bath for 1 hour.
- 5. Centrifuge the sample tubes at $2,000 \times g$ for 10 minutes at room temperature. Decant the supernatant and resuspend the cells in 1.5 mL of Buffer C.
- 6. Centrifuge the samples a second time and resuspend the cell pellet gently in 0.2 mL of Buffer C.
- 7. Spread the entire contents of each tube on an agar plate containing selective growth medium and incubate plates at 30°C for 3 to 4 days. Screen for Mut phenotype (page 38) or select hyper-resistant Geneticin® colonies (page 41).

Lithium Chloride Transformation Method

Introduction

This is a modified version of the procedure described for *S. cerevisiae* (Gietz & Schiestl, 1996). This protocol is provided as an alternative to transformation by electroporation. Transformation efficiency is between 10^2 to 10^3 cfu/µg linearized DNA.

Preparing Solutions

Lithium acetate does not work with *Pichia pastoris*. Use only lithium chloride.

1 M LiCl in distilled, deionized water. Filter sterilize. Dilute as needed with sterile water.

50% polyethylene glycol (PEG-3350) in distilled, deionized water. Filter sterilize. Store in a tightly capped bottle.

2 mg/mL denatured, fragmented salmon sperm DNA in TE (10 mM Tris-HCl, pH 8.0, 1.0 mM EDTA). Store at -20° C.

Preparing Cells

- 1. Grow a 50 mL culture of *Pichia pastoris* in YPD at 30°C with shaking to an OD_{600} of 0.8 to 1.0 (approximately 10^8 cells/mL).
- 2. Harvest the cells and wash with 25 mL of sterile water and centrifuge at $1,500 \times g$ for 10 minutes at room temperature.
- 3. Decant the water and resuspend the cells in 1 mL of 100 mM LiCl.
- 4. Transfer the cell suspension to a 1.5 mL microcentrifuge tube.
- 5. Pellet the cells at maximum speed for 15 seconds and remove the LiCl with a pipet.
- 6. Resuspend the cells in 400 µL of 100 mM LiCl.
- 7. Dispense $50 \,\mu\text{L}$ of the cell suspension into a $1.5 \,\text{mL}$ microcentrifuge tube for each transformation and use immediately.

Do not store on ice or freeze at -20°C.

Lithium Chloride Transformation Method, continued

Transformation

- 1. Boil a 1 mL sample of single-stranded DNA for five minutes, then quickly chill in ice water. Keep on ice. **Note**: It is neither necessary nor desirable to boil the carrier DNA prior to each use. Store a small aliquot at –20°C and boil every 3–4 times the DNA is thawed.
- 2. Centrifuge the LiCl-cell solution from Step 7, previous page. Remove the LiCl with a pipet.
- 3. For each transformation sample, add the following reagents in the order given to the cells. PEG shields the cells from the detrimental effects of the high concentration of LiCl.

```
240 μL 50% PEG
36 μL 1 M LiCl
25 μLl 2 mg/mL single-stranded DNA
Plasmid DNA (5–10 μg) in 50 μL sterile water
```

- 4. Vortex each tube vigorously until the cell pellet is completely mixed (~1 minute).
- 5. Incubate the tube at 30°C for 30 minutes without shaking.
- 6. Heat shock in a water bath at 42°C for 20–25 minutes.
- 7. Centrifuge the tubes at 6,000 to 8,000 rpm and remove the transformation solution with a pipet.
- 8. Gently resuspend the pellet in 1 mL of sterile water.
- 9. Plate 25 to 100 μ L on RDB or MD plates. Incubate the plates for 2–4 days at 30°C. Screen for Mut phenotype (page 38) or select hyper-resistant Geneticin® colonies (page 41).

Direct PCR Screening of Pichia Clones

Introduction

A simple protocol has been reported in the literature to directly test *Pichia* clones for insertion of your gene by PCR (Linder *et al.*, 1996). Briefly, the cells are lysed by a combined enzyme, freezing, and heating treatment. You may use the genomic DNA directly as a PCR template.

Materials Needed

- A culture or single colony of a *Pichia* transformant
- 1.5 mL microcentrifuge tube
- 5 U/μL solution of Lyticase (Sigma)
- 30°C water bath or heat block
- Liquid nitrogen
- Reagents for PCR

Procedure

- 1. Place 10 μ L of a *Pichia pastoris* culture into a 1.5 mL microcentrifuge tube. For relatively dense cultures, dilute 1 μ L of the culture into 9 μ L water. Alternatively, pick a single colony and resuspend in 10 μ L of water.
- 2. Add 5 μ L of a 5 U/ μ L solution of lyticase and incubate at 30°C for 10 minutes.
- 3. Freeze the sample at –80°C for 10 minutes or immerse in liquid nitrogen for 1 minute.
- 4. Set up a 50 μL PCR for a hot start:

10X Reaction Buffer	5 µL
25 mM MgCl2	5 µL
25 mM dNTPs	1 μL
5´ AOX1 primer (10 pmol/μL)	1 μL
3´ AOX1 primer (10 pmol/μL)	1 μL
Sterile water	27 μL
Cell lysate	5 µL
Total Volume	45 µL

- 5. Place the solution in the thermocycler and incubate at 95°C for 5 minutes.
- 6. Add 5 μ L of a 0.16 U/ μ L solution of *Taq* polymerase (0.8 units).
- 7. Cycle 30 times using the following parameters:

Step	Temperature	Time
Denaturation	95°C	1 minute
Annealing	54°C	1 minute
Extension	72°C	1 minute

Include a final extension of 7 minutes at 72°C.

8. Analyze a 10 µL aliquot by agarose gel electrophoresis.

Isolating Total DNA from Pichia

Introduction

The protocol below allows you to isolate DNA from the desired His⁺ recombinant and the untransformed GS115 or KM71 which is suitable for Southern blot analysis, dot/slot blot analysis or genomic PCR. See *Current Protocols in Molecular Biology*, pages 13.11.1 to 13.11.4 (Ausubel *et al.*, 1994), *Guide to Yeast Genetics and Molecular Biology*, pages 322–323 (Strathern and Higgins, 1991), or (Holm *et al.*, 1986).

Solutions

Prepare the following solutions, because there is not enough of some of these reagents in the kit to perform this experiment.

Minimal Medium (MD, MGYH)

Sterile water

SCED (1 M sorbitol, 10 mM sodium citrate, pH 7.5, 10 mM EDTA, 10 mM DTT) Zymolyase, 3 mg/mL stock solution in water (Seikagaku America, Inc., 1-800-237-4512)

1% SDS in water

5 M potassium acetate, pH 8.9

TE buffer, pH 7.4 (10 mM Tris-HCl, pH 7.4, 1 mM EDTA, pH 8.0)

7.5 M ammonium acetate, pH 7.5

Phenol:chloroform (1:1 v/v)

Preparing Cells

- 1. Grow at 30° C the recombinant strain and the parent strain to an OD_{600} of 5–10 in 10 mL of minimal media such as MD or MGY (recombinant) or MDH or MGYH (GS115 or KM71).
- 2. Collect the cells by centrifugation at $1,500 \times g$ for 5-10 minutes at room temperature.
- 3. Wash the cells with 10 mL sterile water by centrifugation as in Step 2.

Spheroplasting and Lysis

- 1. Resuspend the cells in 2 mL of SCED buffer, pH 7.5. Make this solution fresh.
- 2. Add 0.1–0.3 mg of Zymolyase (mix well before adding to the cells). Incubate at 37°C for 50 minutes to achieve < 80% spheroplasting (monitor the percent spheroplasting using the procedure on pages 33–34).
- 3. Add 2 mL of 1% SDS, mix **gently** and set on ice (0 to 4°C) for 5 minutes.
- 4. Add 1.5 mL of 5 M potassium acetate, pH 8.9, and mix gently.
- 5. Centrifuge at $10,000 \times g$ for 5–10 minutes at 4°C and save the supernatant.

Isolating Total DNA from Pichia, continued

DNA Precipitation

- 1. Transfer the supernatant from Step 5, the previous page, and add 2 volumes of ethanol to the supernatant. Incubate at room temperature for 15 minutes.
- 2. Centrifuge at $10,000 \times g$ for 20 minutes at 4° C.
- 3. Resuspend the pellet **gently** in 0.7 mL of TE buffer, pH 7.4 and transfer to a microcentrifuge tube.
- 4. **Gently** extract with an equal volume of phenol:chloroform (1:1 v/v) followed by an equal volume of chloroform:isoamyl alcohol (24:1). Split the aqueous layer into two microcentrifuge tubes.
- 5. Add 1/2 volume of 7.5 M ammonium acetate, pH 7.5, and 2 volumes of ethanol to each tube. Place on dry ice for 10 minutes or at –20°C for 60 minutes.
- 6. Centrifuge at $10,000 \times g$ for 20 minutes at 4°C and wash the pellets once with 1 mL of 70% ethanol. Briefly air dry the pellets and resuspend each one in 50 μ L of TE buffer, pH 7.5.
- 7. Determine the concentration of the DNA sample. You can store the samples at -20°C separately or combined until ready for use.

Detecting Multiple Integration Events

Introduction

It has been demonstrated in a number of papers (Brierley *et al.*, 1994; Clare *et al.*, 1991a; Romanos *et al.*, 1991; Scorer *et al.*, 1993; Scorer *et al.*, 1994) that multiple integration events may increase the levels of protein expressed. If the expression of your protein is low, you may wish to isolate multicopy integrants. Using the protocol detailed on the previous page, isolate genomic DNA from the His⁺ recombinants that are to be analyzed, as well as from the untransformed GS115 or KM71 to control for any background hybridization. You can use the DNA in either of the following analyses to detect multicopy integration.

Southern Blot Analysis

For a detailed description of this technique as applied to *Pichia pastoris*, see (Clare, *et al.*, 1991a). It is very important to digest your DNA with the right restriction enzyme(s) to generate a blot of digested and gel-separated genomic DNA. It is also important to understand that your strategy will be different if you use pPIC3.5K versus pAO815 to generate your multiple copies. Digesting DNA from *Pichia* recombinants containing multiple copies produces a band that varies in intensity depending on the number of copies of your gene. It is very important to include a control to show the intensity of a single copy gene. You can quantify the band intensities using densitometry to estimate relative gene dosage.

Materials Needed for Quantitative Dot Blot

Prepare 10–15 mL of each following of solutions for each dot blot:

- 50 mM EDTA, 2.5% β-mercaptoethanol pH 9
- 1 mg/mL Zymolyase 100T in water (Seikagaku America, Inc., 1-800-237-4512)
- 0.1 N NaOH, 1.5 M NaCl
- 2X SSC
- 3MM paper.

Detecting Multiple Integration Events, continued

Quantitative Dot Blot Procedure

The following protocol is a summary of a rapid DNA dot blot technique to detect multiple integrants (Romanos, *et al.*, 1991). It is very important to spot equivalent numbers of cells onto filters to quantify the copy number. Alternatively, you may isolate genomic DNA and spot it directly onto nitrocellulose or nylon, fix it, and analyze it.

- 1. Grow Mut⁺ or Mut^S transformants in individual wells of a 96-well microtiter plate in 200 µL of YPD broth at 30°C until all wells have approximately the same density. This may necessitate several passages.
 - Alternatively, you may grow individual transformants in culture tubes and normalize the absorbance at 600 nm by adding medium.
- 2. Filter $50~\mu L$ of each sample onto a nitrocellulose or nylon filter placed into a dot (slot) blot apparatus using multi-channel pipettor. Air dry filters.
- 3. To lyse the cells on the filter, treat the filter with four solutions as follows: place two sheets of 3 MM paper in a tray and soak with 10–15 mL of 50 mM EDTA, 2.5% β -mercaptoethanol pH 9. Make sure that the paper is uniformly soaked and that there are no puddles. Place the nitrocellulose filter face down on the treated 3MM paper. Incubate the filter for 15 minutes at room temperature.
- 4. Remove the nitrocellulose filter from the 3MM paper and replace the 3MM paper with two new sheets. Soak them with 10–15 mL of 1 mg/mL Zymolyase 100T as described in Step 3. Place the nitrocellulose filter face down on the 3MM paper and incubate it for 4 hours at 37°C.
- 5. Remove the nitrocellulose filter from the paper and replace the paper with two new sheets. Soak with 10–15 mL of 0.1 N NaOH, 1.5 M NaCl. Place the nitrocellulose filter face down on the paper and incubate for 5 minutes at room temperature.
- 6. Remove the nitrocellulose filter from the 3MM paper and replace the 3MM paper with two new sheets. Soak them with 10–15 mL of 2X SSC. Place the nitrocellulose filter face down on the 3MM paper and incubate it for 5 minutes at room temperature. Repeat.
- 7. Bake the nitrocellulose filters at 80°C or UV-crosslink the DNA to nylon. The filters may be probed with a nonradioactive-labeled or random-primed, ³²P-labeled probe complementary to your gene.

You can identify multi-copy integrants by a strong hybridization signal relative to the single copy control. You can then quantify dot blots for copy number by densitometry of the film or blot, or by using a β -scanner (if radiolabeled).

Determining the Copy Number of Multiple Integrants, continued

General Guidelines

- Use standard procedures and solutions for Southern blotting as outlined in *Molecular Cloning: A Laboratory Manual* (Sambrook, *et al.*, 1989), pages 9.31–9.58.
- Isolate genomic DNA and quantify it using fluorometry. Be sure to eliminate RNA. It is very important to load the same amount of DNA into each lane to accurately determine the copy number.
- Probe your Southern blot with probes to *HIS4* and your gene. The point mutation in the *his4* gene in the host strain will not interfere with hybridization if you make the probe complementary to the wild-type gene.
- If you used pPIC3.5K to generate multimers, use *Bgl* II to digest your DNA (Clare, *et al.*, 1991a). If you used pPIC3.5K all multimers are NOT necessarily in a head-to-tail configuration. Some multimers may be head-to-head and others tail-to-tail. We recommend that you think about what products may be produced. An expression cassette in the opposite orientation may produce a different band. The number of multiple copies will cause one or two bands (depending on orientation) in the Southern blot to increase in intensity once you have >2 copies.
- If you used pAO815 to generate multimers, use *Bgl* II and *BamH* I to digest the genomic DNA and release the multimer. The molecular weight of the band allows you to determine the number of multimers. If this multimer is too large, you may wish to digest with an enzyme like *Sac* I. This collapses the multimer into single fragments containing your gene. These fragments will produce a band that will be quite intense. The relative intensity of this band versus a band containing a single copy of your gene allows you to determine the copy number.
- Bgl II digested DNA from GS115 and GS115 transformed with pPIC3.5K or pAO815 produces bands of 2.8 kb (the genomic copy of HIS4), and ~6.7 kb (the vector derived copy of HIS4), respectively, when probed with a complementary fragment to HIS4.

Procedure for Total RNA Isolation from Pichia

Introduction

This protocol is designed to isolate 60– $300 \, \mu g$ total RNA (Schmitt $et \, al., 1990$) from Pichia which is suitable for mRNA isolation using Invitrogen's FastTrack® 2.0 or Micro FastTrack™ 2.0 mRNA Isolation Kit. If you wish to use another protocol, scale-up the reaction to yield about 2 mg of total RNA per time point. The mRNA is for Northern blot analysis of Pichia recombinants to determine if the gene of interest is being induced and transcribed. RNA isolation should be done from induced cultures using an uninduced culture as a negative control.

Solutions

Use DEPC-treated water and equipment free of RNase.

MGYH or BMGY medium

3 M sodium acetate, pH 5.3

Buffered phenol

10% SDS in DEPC treated water

Chloroform:isoamyl alcohol (24:1)

DEPC-treated water

Phenol:chloroform (1:1)

65°C water bath

AE buffer (50 mM sodium acetate, pH 5.3, 1 mM EDTA)

Growing Cells

- 1. Grow up two cultures (100–200 mL in MGY or BMGY), but induce only one of them. Use the same protocol for induction that you used in the **Expression** section.
- 2. Take 10 mL time points at 1, 2, 3, 4, and 6 days.
- 3. Harvest the cells from each time point by centrifugation at $1500 \times g$ for 10 minutes at room temperature.
- 4. Resuspend the cell pellet in 400 μL of AE buffer and transfer it to a microcentrifuge tube.

Procedure for Total RNA Isolation from Pichia, continued

Lysing Cells

- 1. Add $40~\mu L$ of 10% SDS to the microcentrifuge tube containing the cell suspension and vortex it for ~20 seconds.
- 2. Add an equal volume (450–500 μ L) of buffer saturated phenol to the tube and vortex it for ~20 seconds.
- 3. Incubate the cell suspension at 65°C for 4 minutes.
- 4. Incubate the cell suspension in a dry ice/ethanol bath until crystals show (~1 minute). Centrifuge the tube at maximum speed for 2 minutes at 4°C.
- 5. Transfer the aqueous phase to a new centrifuge tube, add an equal volume of phenol/chloroform, and vortex it for ~20 seconds. Centrifuge the tube at maximum speed for 2 minutes at 4°C.
- 6. Remove the upper phase to a new tube and add $40~\mu L$ of 3~M sodium acetate, pH 5.3 and 2.5 volumes of 100% ethanol ($-20^{\circ}C$). Centrifuge the tube at maximum speed for 15~minutes at $4^{\circ}C$. Remove the supernatant.
- 7. Wash the pellet with 80% ethanol and air dry it briefly. Resuspend total RNA in 20 μ L of DEPC-treated water and store it at –80°C. The average yield is 60–300 μ g of total RNA.

mRNA Isolation and Northern Analysis

See (Ausubel *et al.*, 1994) for a protocol for mRNA isolation and Northern analysis. The FastTrack® 2.0 mRNA Kit is designed to isolate mRNA from 0.2 to 1 mg total RNA and the Micro-FastTrack™ 2.0 Kit is designed to isolate mRNA from ~100 µg total RNA (see page 88 for ordering information. You will need ~1–5 µg mRNA per time point.

β-Galactosidase Assay

Introduction

The GS115 β -Gal strain is provided as a His⁺ Mut⁺ intracellular expression control. Growth of the strain during Mut⁺ expression provides a positive control for expression conditions. The cell-free β -galactosidase assay described below can also be found in (Miller, 1972), page 403, and can be used to evaluate expression of β -galactosidase.

Required Solutions

- A fresh crude cell lysate of GS115 β-Gal (see page 57)
- Z buffer
- ONPG solution
- 1 M sodium carbonate solution

Recipes for the solutions are below.

Z Buffer

60 mM Na₂HPO₄-7H₂O

40 mM NaH₂PO₄-H₂O

10 mM KCl

1 mM MgSO₄-7H₂O

50 mM β-mercaptoethanol

pH 7.0

1. Dissolve the following in 950 mL deionized water:

Na ₂ HPO ₄ -7H ₂ O	16.1 g
NaH ₂ PO ₄ -H ₂ O	5.5 g
KCl	0.75 g
MgSO ₄ -7H ₂ O	0.246 g
β-mercaptoethanol	2.7 mL

- 2. Adjust pH to 7.0 with either NaOH or HCl and bring the volume up to 1 liter with water.
- 3. **Do not autoclave!** Store at 4°C.

ONPG Solution

4 mg/mL in 100 mM phosphate buffer, pH 7.0

1. Dissolve the following in 90 mL deionized water:

Na ₂ HPO ₄ -7H ₂ O	1.61 g
NaH ₂ PO ₄ -H ₂ O	0.55 g

- 2. Adjust pH to 7.0 with either NaOH or HCl
- 3. Add 400 mg ONPG. Stir to dissolve and bring the volume up to 100 mL with water.
- 4 Store at 4°C away from light.

1 M Sodium Carbonate

Dissolve 12.4 g sodium carbonate in 100 mL of deionized water. Store at room temperature.

β-Galactosidase Assay, continued

Procedure

- 1. Determine protein concentration of your lysate by Lowry, Bradford, or BCA assay.
- 2. Equilibrate Z buffer, ONPG solution, and sodium carbonate solution to 28°C.
- 3. Add 10–50 μ L of your crude assay to 1 mL of Z buffer and equilibrate at 28°C. As a control for spontaneous hydrolysis of ONPG, add an aliquot of your lysis buffer to 1 mL of Z buffer.
- 4. To initiate the reaction, add 0.2~mL~4~mg/mL ONPG to each of the tubes in Step 3.
- 5. Incubate the samples and the control at 28°C until a faint yellow color develops. This should occur at least 10 minutes after the start of the assay to ensure accurate data. Note that the tube with no lysate may not change color.
- 6. Stop the reaction by adding 0.5 mL of 1 M sodium carbonate to each tube. Record the length of incubation for each sample.
- 7. Read the OD_{420} against the control containing buffer alone.
- 8. Determine the protein concentration of your lysate in mg/mL.



If the reaction turns yellow too quickly, you need to dilute your lysate. Try successive 10-fold dilutions of the lysate using your lysis buffer until the reaction starts turning yellow after 10 minutes. This is to ensure that you are measuring a true initial rate.

Determining Specific Activity

Use the following formula to determine the specific activity of the β -galactosidase in units/mg total protein:

β-galactosidase units/mg total protein = $\frac{OD_{420} \times 380}{\text{minutes at } 28^{\circ}\text{C} \times \text{mg protein in reaction}}$

Sample Calculation

Here is a sample calculation: Extract concentration = 10 mg/mL

Assay 10 µL of a 1/100 dilution

Time = 10 minutes

 $OD_{420} = 0.4$

The amount of protein in the reaction = $0.01 \text{ mL} \times 0.01$ (dilution factor) $\times 10 \text{ mg/mL}$

= 0.001 mg protein in the reaction

The specific activity = $\frac{0.400 \times 380}{10 \times 0.001}$ = 15,200 units/mg protein

Pure β -galactosidase has an activity of 300,000 units/mg protein.

Accessory Products

Accessory Products

Many of the reagents supplied in *Pichia* Expression Kit, as well as other reagents and kits used in *Pichia* expression experiments are available separately from Invitrogen. Ordering information is provided below. For more information, refer to www.invitrogen.com or contact Technical Support (see page 89).

Product	Amount	Cat. no.
Pichia Spheroplast Module	10 preparations	K1720-01
(Box 1 and Box 2)	(50 transformations)	
Geneticin [®] , powder	1 g	11811-023
	5 g	11811-031
	25 g	11811-098
Geneticin [®] , liquid (50mg/mL)	20 mL	10131-035
	100 mL	10131-027
Ampicillin Sodium Salt, irradiated	200 mg	11593-027
TA Cloning® Kit (with pCR®2.1 vector)	20 reactions	K2000-01
PureLink™ HiPure Plasmid Miniprep Kit	25 preparations	K2100-02
	100 preparations	K2100-03
Platinum® Taq DNA Polymerase	100 reactions	10966-018
	250 reactions	10966-026
	500 reactions	10966-034
Taq DNA Polymerase, Recombinant	100 units	10342-053
	500 units	10342-020
Easy-DNA™ Kit	1 kit	K1800-01
FastTrack® 2.0 mRNA Isolation Kit	1 kit	K1593-02
Micro-FastTrack [™] 2.0 mRNA Isolation Kit	1 kit	K1520-02
pAO815	20 μg	V180-20
pPIC3.5K	20 μg	V173-20
pPIC9K	20 μg	V175-20

Technical Support

Web Resources



Visit the Invitrogen website at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.
- Complete technical support contact information
- Access to the Invitrogen Online Catalog
- Additional product information and special offers

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our website (www.invitrogen.com).

Corporate Headquarters: 5791 Van Allen Way Carlsbad, CA 92008 USA Tel: 1 760 603 7200

Tel (Toll Free): 1 800 955 6288

Fax: 1760 602 6500

E-mail:

tech_support@invitrogen.com

Japanese Headquarters: LOOP-X Bldg. 6F 3-9-15, Kaigan Minato-ku, Tokyo 108-0022

Tel: 81 3 5730 6509 Fax: 81 3 5730 6519

E-mail: jpinfo@invitrogen.com

European Headquarters: Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF, UK

Tel: +44 (0) 141 814 6100 Tech Fax: +44 (0) 141 814 6117

E-mail:

eurotech@invitrogen.com

SDS

SDSs (Safety Data Sheets) are available at www.invitrogen.com/sds.

Certificate of Analysis

The Certificate of Analysis (CofA) provides detailed quality control information for each product. The CofA for each product is available on our website at www.invitrogen.com/cofa, and is searchable by product lot number, which is printed on each box.

Technical Support, continued

Limited Warranty

Invitrogen (a part of Life Technologies Corporation) is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

All Invitrogen products are warranted to perform according to specifications stated on the certificate of analysis. The Company will replace, free of charge, any product that does not meet those specifications. This warranty limits the Company's liability to only the price of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. The Company reserves the right to select the method(s) used to analyze a product unless the Company agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore the Company makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives. Life Technologies Corporation shall have no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the terms of the most restrictive limited use label license shall control. Life Technologies Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Life Technologies Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Life Technologies is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: outlicensing@invitrogen.com.

Purchaser Notification, continued

Limited Use Label License No. 74: *Pichia Pastoris* Expression System The *Pichia* Expression System is based on the yeast *Pichia pastoris*. *Pichia pastoris* was developed into an expression system by scientists at Salk Institute Biotechnology/ Industry Associates (SIBIA) and Phillips Petroleum for high-level expression of recombinant proteins. All patents for *Pichia pastoris* and licenses for its use as an expression system are owned by Research Corporation Technologies (RCT), Inc., Tucson, Arizona.

Invitrogen has an exclusive license to sell *Pichia* expression kits and vectors to scientists for research purposes only, under the terms described below. Use of *Pichia pastoris* by commercial entities for any commercial purpose requires the user to obtain a commercial license as detailed below. Before using any *Pichia* expression product, please read the following license agreement. If you do not agree to be bound by its terms, contact Invitrogen within 10 days for authorization to return the unused *Pichia* expression products and to receive a full refund. If you do agree to the terms of this license agreement, please complete the User Registration Card and return it to Invitrogen before using the product. *Pichia pastoris* is covered by one or more of the following U.S. patents and corresponding foreign patents owned and licensed by Research Corporation Technologies: 4,683,293; 4,808,537; 4,812,405; 4,818,700; 4,837,148; 4,855,231; 4,857,467; 4,879,231; 4,882,279; 4,885,242; 4,895,800; 4,929,555; 5,002,876; 5,004,688; 5,032,516; 5,122,465; 5,135,868; 5,166,329

Invitrogen Corporation ("Invitrogen") grants you a non-exclusive license to use the enclosed Pichia expression vectors ("Expression Vector") for academic research or for evaluation purposes only. The Expression Vectors are being transferred to you in furtherance of, and reliance on, such license. You may not use the Expression Vectors for any commercial purpose without a license for such purpose from Research Corporation Technologies, Inc., Tucson, Arizona. Commercial purposes include: any use of Expression Products or Expression Vectors in a Commercial Product; any use of Expression Products or Expression Vectors in the manufacture of a Commercial Product; any sale of Expression Products; any use of Expression Products or the Expression Kit to facilitate or advance research or development directed to a Commercial Product; and any use of Expression Products or the Expression Kit to facilitate or advance any research or development program the results of which will be directly applied to the development or manufacture of a Commercial Product. "Expression Products" means products expressed with the Expression Kit, or with the use of any Pichia expression vectors (including the Expression Vector) or host strains. "Commercial Product" means any product intended for sale or commercial use. Commercial entities may conduct their evaluation for one year at which time this license automatically terminates. Commercial entities will be contacted by Research Corporation Technologies during the evaluation period regarding their desire for a commercial license.

Purchaser Notification, continued

Limited Use Label License No. 74: *Pichia Pastoris* Expression System, continued Access to the Expression Kit and Vector must be limited solely to those officers, employees and students of your institution who need access to perform the above-described research or evaluation. You must inform each such officer, employee and student of the provisions of this license agreement and require them to agree, in writing, to be bound by the provisions of this license agreement. You may not distribute any Expression Vector or host strain contained herein or in the Expression Kit to others, even those within your own institution. You may only transfer modified, altered, or original material from the Expression Kit or Vector to a third party following written notification of, and written approval from, Invitrogen so that the recipient can be licensed. You may not assign, sublicense, rent, lease or otherwise transfer this license agreement or any of the rights or obligation there under, except as expressly permitted by Invitrogen and RCT. This license agreement is effective until terminated. You may terminate it at any time by destroying all Pichia Expression products in your control. It will also terminate automatically if you fail to comply with the terms and conditions of the license agreement. You shall, upon termination of the license agreement, destroy all *Pichia* Expression products in your control, and so notify Invitrogen in writing. You may contact Research Corporation Technologies at the following address: David Bramhill, Ph.D., Director, Research Corporation Technologies, 5210 E. Williams Circle, Suite 240, Tucson, Arizona 85711-4410. Tel: 520-748-4443, Fax:520-748-0025, Website: http://rctech.com/licensing/gene-expression.php, Email: dbramhill@rctech.com.

References

- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994) *Current Protocols in Molecular Biology*, Greene Publishing Associates and Wiley-Interscience, New York
- Barr, K. A., Hopkins, S. A., and Sreekrishna, K. (1992) Protocol for Efficient Secretion of HSA Developed from *Pichia pastoris*. Pharm. Eng. 12, 48-51
- Bollag, D. M., M. D., R., and Edelstein, S. J. (1996) Protein Methods, Second Ed., Wiley-Liss, New York
- Brake, A. J., Merryweather, J. P., Coit, D. G., Heberlein, U. A., Masiarz, G. R., Mullenbach, G. T., Urdea, M. S., Valenzuela, P., and Barr, P. J. (1984) a-Factor-Directed Synthesis and Secretion of Mature Foreign Proteins in *Saccharomyces cerevisiae*. Proc. Natl. Acad. Sci. USA *81*, 4642-4646
- Brandes, H. K., Hartman, F. C., Lu, T.-Y. S., and Larimer, F. W. (1996) Efficient Expression of the Gene for Spinach Phosphoribulokinase in *Pichia pastoris* and Utilization of the Recombinant Enzyme to Explore the Role of Regulatory Cysteinyl Residues by Site-directed Mutagenesis. J. Biol. Chem. 271, 6490-6496
- Brankamp, R. G., Sreekrishna, K., Smith, P. L., Blankenship, D. T., and Cardin, A. D. (1995) Expression of a Synthetic Gene Encoding the Anticoagulant-Antimetastatic Protein Ghilanten by the Methylotropic Yeast *Pichia pastoris*. Protein Expression and Purification *6*, 813-820
- Brierley, R. A., Davis, G. R., and Holtz, G. C. (1994) Production of Insulin-Like Growth Factor-1 in Methylotrophic Yeast Cells. United States Patent *5,324,639*
- Buckholz, R. G., and Gleeson, M. A. G. (1991) Yeast Systems for the Commercial Production of Heterologous Protein. Bio/Technology *9*, 1067-1072
- Cavener, D. R., and Stuart, C. R. (1991) Eukaryotic Start and Stop Translation Sites. Nucleic Acids Res. 19, 3185-3192
- Clare, J. J., Rayment, F. B., Ballantine, S. P., Sreekrishna, K., and Romanos, M. A. (1991a) High-level Expression of Tetanus Toxin Fragment c in *Pichia pastoris* Strains Containing Multiple Tandem Integrations of the Gene. Bio/Technology *9*, 455-460
- Clare, J. J., Romanos, M. A., Rayment, F. B., Rowedder, J. E., Smith, M. A., Payne, M. M., Sreekrishna, K., and Henwood, C. A. (1991b) Production of Epidermal Growth Factor in Yeast: High-Level Secretion Using *Pichia pastoris* Strains Containing Multiple Gene Copies. Gene 105, 205-212
- Cregg, J. M., Barringer, K. J., and Hessler, A. Y. (1985) *Pichia pastoris* as a Host System for Transformations. Mol. Cell. Biol. *5*, 3376-3385
- Cregg, J. M., and Higgins, D. R. (1995) Production of Foreign Proteins in the Yeast *Pichia pastoris*. Canadian J. Botany Supp. *73*, 5981-5987
- Cregg, J. M., Madden, K. R., Barringer, K. J., Thill, G., and Stillman, C. A. (1989) Functional Characterization of the Two Alcohol Oxidase Genes from the Yeast, *Pichia pastoris*. Mol. Cell. Biol. *9*, 1316-1323
- Cregg, J. M., Tschopp, J. F., Stillman, C., Siegel, R., Akong, M., Craig, W. S., Buckholz, R. G., Madden, K. R., Kellaris, P. A., Davis, G. R., Smiley, B. L., Cruze, J., Torregrossa, R., Velicelebi, G., and Thill, G. P. (1987) High-Level Expression and Efficient Assembly of Hepatitis B Surface Antigen in the Methylotrophic Yeast *Pichia pastoris*. Bio/Technology *5*, 479-485
- Cregg, J. M., Vedvick, T. S., and Raschke, W. C. (1993) Recent Advances in the Expression of Foreign Genes in *Pichia pastoris*. Bio/Technology *11*, 905-910
- Despreaux, C. W., and Manning, R. F. (1993) The *dacA* Gene of *Bacillus stearothermophilus* Coding for D-Alanine carboxypeptidase: Cloning, Structure, and Expression in *Escherichia coli* and *Pichia pastoris*. Gene 131, 35-41

References, continued

- Deutscher, M. P. (ed) (1990) *Guide to Protein Purification* Vol. 182. Methods in Enzymology. Edited by Abelson, J. N., and Simon, M. I., Academic Press, San Diego, CA.
- Digan, M. E., Lair, S. V., Brierley, R. A., Siegel, R. S., Williams, M. E., Ellis, S. B., Kellaris, P. A., Provow, S. A., Craig, W. S., Velicelebi, G., Harpold, M. M., and Thill, G. P. (1989) Continuous Production of a Novel Lysozyme via Secretion from the Yeast *Pichia pastoris*. Bio/Technology 7, 160-164
- Ellis, S. B., Brust, P. F., Koutz, P. J., Waters, A. F., Harpold, M. M., and Gingeras, T. R. (1985) Isolation of Alcohol Oxidase and Two other Methanol Regulatable Genes from the Yeast, *Pichia pastoris*. Mol. Cell. Biol. *5*, 1111-1121
- Fryxell, K. B., O'Donoghue, K., Graeff, R. M., Lee, H. C., and Branton, W. D. (1995) Functional Expression of Soluble Forms of Human CD38 in *Escherichia coli* and *Pichia pastoris*. Protein Expression and Purification *6*, 329-336
- Garcia, J. N., Aguiar, J. A., Gill, M., Alvarez, A., Morales, J., Ferrero, J., Gonzalez, B., Padron, G., and Menendez, A. (1995) High Level Expression of Human IFN-a2b in *Pichia pastoris*. Biotecnologia Aplicada 12, 152-155
- Gietz, R. D., and Schiestl, R. H. (1996) in *Methods in Molecular Biology* (Evans, I. H., ed), Humana Press, Totowa, NJ
- Guo, W., Gonzalez-Candelas, L., and Kolattukudy, P. E. (1995) Cloning of a New Pectate Lyase Gene *pelC* from *Fusarium solani* f. sp. *pisi* (*Nectria haematococca*, Mating Type VI) and Characterization of the Gene Product Expressed in *Pichia pastoris*. Arch. Biochem. Biophys *323*, 352-360
- Guthrie, C., and Fink, G. R. (eds) (1991) *Guide to Yeast Genetics and Molecular Biology* Vol. 194. Methods in Enzymology. Edited by Abelson, J. N., and Simon, M. I., Academic Press, San Diego, CA
- Hagenson, M. J., Holden, K. A., Parker, K. A., Wood, P. J., Cruze, J. A., Fuke, M., Hopkins, T. R., and Stroman, D. W. (1989) Expression of Streptokinase in *Pichia pastoris* Yeast. Enzyme Microbiol. Technol. *11*, 650-656
- Henikoff, S., and Cohen, E. H. (1984) Sequences Responsible for Transcription Termination on a Gene Segment in *Saccharomyces cerevisiae*. Mol. Cell. Biol. 4, 1515-1520
- Hinnen, A., Hicks, J. B., and Fink, G. R. (1978) Transformation of Yeast Cells. Proc. Natl. Acad. Sci. USA 75, 1292-1293
- Irniger, S., Egli, C. M., and Braus, G. H. (1991) Different Classes of Polyadenylation Sites in the Yeast *Saccharomyces cerevisiae*. Mol. Cell. Bio. *11*, 3060-3069
- Johnston, M. (1987) A Model Fungal Gene Regulatory Mechanism: the *GAL* Genes of *Saccharomyces Cerevisiae*. Microbiol. Rev 51, 458-476
- Koutz, P. J., Davis, G. R., Stillman, C., Barringer, K., Cregg, J. M., and Thill, G. (1989) Structural Comparison of the *Pichia pastoris* Alcohol Oxidase Genes. Yeast 5, 167-177
- Kozak, M. (1987) An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. Nucleic Acids Res. 15, 8125-8148
- Kozak, M. (1990) Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. Proc. Natl. Acad. Sci. USA 87, 8301-8305
- Laroche, Y., Storme, V., Meutter, J. D., Messens, J., and Lauwereys, M. (1994) High-Level Secretion and Very Efficient Isotopic Labeling of Tick Anticoagulant Peptide (TAP) Expressed in the Methylotrophic Yeast, *Pichia pastoris*. Bio/Technology 12, 1119-1124
- Linder, S., Schliwa, M., and Kube-Granderath, E. (1996) Direct PCR Screening of *Pichia pastoris* Clones. BioTechniques 20, 980-982
- Nico-Farber, K., Harder, W., Ab, G., and Veenhuis, M. (1995) Review: Methylotrophic Yeasts as Factories for the Production of Foreign Proteins. Yeast *11*, 1331-1344

References, continued

- Paifer, E., Margolles, E., Cremata, J., Montesino, R., Herrera, L., and Delgado, J. M. (1994) Efficient Expression and Secretion of Recombinant Alpha Amylase in *Pichia pastoris* Using Two Different Signal Sequences. Yeast *10*, 1415-1419
- Ridder, R., Schmitz, R., Legay, F., and Gram, H. (1995) Generation of Rabbit Monoclonal Antibody Fragments from a Combinatorial Phage Display Library and Their Production in the Yeast *Pichia pastoris*. Bio/Technology *13*, 255-260
- Rodriguez, M., Rubiera, R., Penichet, M., Montesinos, R., Cremata, J., Falcon, V., Sanchez, G., Bringas, R., Cordoves, C., Valdes, M., Lleonart, R., Herrera, L., and delaFuente, J. (1994) High Level Expression of the *B. microplus* Bm86 Antigen in the Yeast *Pichia pastoris* Forming Highly Immunogenic Particles for Cattle. J. Biotechnol. *33*, 135-146
- Romanos, M. (1995) Advances in the Use of Pichia pastoris for High-Level Expression. Curr. Opin. Biotech. *6*, 527-533
- Romanos, M. A., Clare, J. J., Beesley, K. M., Rayment, F. B., Ballantine, S. P., Makoff, A. J., Dougan, G., Fairweather, N. F., and Charles, I. G. (1991) Recombinant *Bordetella pertussis* Pertactin p69 from the Yeast *Pichia pastoris* High Level Production and Immunological Properties. Vaccine *9*, 901-906
- Romanos, M. A., Scorer, C. A., and Clare, J. J. (1992) Foreign Gene Expression in Yeast: A Review. Yeast 8, 423-488
- Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Second Ed., Cold Spring Harbor Laboratory Press, Plainview, New York
- Scorer, C. A., Buckholz, R. G., Clare, J. J., and Romanos, M. A. (1993) The Intracellular Production and Secretion of HIV-1 Envelope Protein in the Methylotrophic Yeast *Pichia pastoris*. Gene *136*, 111-119
- Scorer, C. A., Clare, J. J., McCombie, W. R., Romanos, M. A., and Sreekrishna, K. (1994) Rapid Selection Using G418 of High Copy Number Transformants of *Pichia pastoris* for High-level Foreign Gene Expression. Bio/Technology 12, 181-184
- Sreekrishna, K., Nelles, L., Potenz, R., Cruse, J., Mazzaferro, P., Fish, W., Fuke, M., Holden, K., Phelps, D., Wood, P., and Parker, K. (1989) High Level Expression, Purification, and Characterization of Recombinant Human Tumor Necrosis Factor Synthesized in the Methylotrophic Yeast *Pichia pastoris*. Biochemistry *28*, 4117-4125
- Sreekrishna, K., Potenz, R. H. B., Cruze, J. A., McCombie, W. R., Parker, K. A., Nelles, L., Mazzaferro, P. K., Holden, K. A., Harrison, R. G., Wood, P. J., Phelps, D. A., Hubbard, C. E., and Fuke, M. (1988) High Level Expression of Heterologous Proteins in Methylotrophic Yeast *Pichia pastoris*. J. Basic Microbiol. 28, 265-278
- Sun, J., Coughlin, P., Salem, H. H., and Bird, P. (1995) Production and Characterization of Recombinant Human Proteinase Inhibitor 6 Expressed in *Pichia pastoris*. Biochim. Biophys. ACTA 1252, 28-34
- Thill, G. P., Davis, G. R., Stillman, C., Holtz, G., Brierley, R., Engel, M., Buckholz, R., Kinney, J., Provow, S., Vedvick, T., and Siegel, R. S. (1990) in *International Symposium on the Genetics of Microorganisms* Vol. 2, pp. 477-490
- Tschopp, J. F., Brust, P. F., Cregg, J. M., Stillman, C., and Gingeras, T. R. (1987a) Expression of the *lacZ* Gene from Two Methanol Regulated Promoters in *Pichia pastoris*. Nucleic Acids Res. *15*, 3859-3876
- Tschopp, J. F., Sverlow, G., Kosson, R., Craig, W., and Grinna, L. (1987b) High Level Secretion of Glycosylated Invertase in the Methylotrophic Yeast *Pichia pastoris*. Bio/Technology *5*, 1305-1308

References, continued

- Varki, A., and Freeze, H. H. (1994) The Major Glycosylation Pathways of Mammalian Membranes. A Summary. Subcell. Biochem. 22, 71-100
- Vedvick, T., Buckholz, R. G., Engel, M., Urcan, M., Kinney, J., Provow, S., Siegel, R. S., and Thill, G. P. (1991) High-level Secretion of Biologically Active Aprotonin from the Yeast *Pichia pastoris*. J. Ind. Microbiol. 7, 197-201
- Vozza, L. A., Wittwer, L., Higgins, D. R., Purcell, T. J., Bergseid, M., Collins-Racie, L. A., LaVallie, E. R., and Hoeffler, J. P. (1996) Production of a Recombinant Bovine Enterokinase Catalytic Subunit in the Methylotrophic Yeast *Pichia pastoris*. Bio/Technology *14*, 77-81
- Wagner, S. L., Siegel, R. S., Vedvick, T. S., Raschke, W. C., and VanNostrand, W. E. (1992) High-level Expression, Purification, and Characterization of the Kunitz-type Protease Inhibitor Domain of Protease Nixin-2/amyloid b-Protein Precursor. Biochem. Biophys. Res. Commun. 186, 1138-1145
- Wegner, G. H. (1990) Emerging Applications of the Methylotrophic Yeasts. FEMS Microbiology Reviews 87, 279-284
- Weiss, H. M., Haase, W., Michel, H., and Reilander, H. (1995) Expression of Functional Mouse 5-HT_{5A} Serotonin Receptor in the Methylotrophic Yeast *Pichia pastoris*: Pharmacological Characterization and Localization. FEBS 377, 451-456
- Zaret, K. S., and Sherman, F. (1984) Mutationally Altered 3´ Ends of Yeast *CYC1* mRNA Affect Transcript Stability and Translational Efficiency. J. Mol. Biol. 177, 107-136

©2009 Life Technologies Corporation. All rights reserved.

For research use only. Not intended for any animal or human therapeutic or diagnostic use.



Corporate Headquarters

5791 Van Allen Way Carlsbad, CA 92008 T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information, visit our web site at **www.invitrogen.com**