

EtOH Precipitation of DNA

- Primarily used to clean-up DNA that has been stored in phenol
- 1. In the 2ml epi tube containing the DNA sample to clean, add 1ml of COLD 100% ethanol and 100ul sodium acetate (NaOAc) pH 5.2
 - a. Mix by inverting
- 2. Place in -20C freezer for at least 3 hours, overnight is fine
- 3. Spin at 14,000rpm in a chilled centrifuge at 4C for 10-15min
 - a. Room temp centrifugation may not lower yield (I suspect)
- 4. Discard supernatant
 - a. Careful not to disturb pellet!
 - b. Pellet may be invisible
- 5. Re-suspend in 1ml fresh 70% ethanol, centrifuge for 5min at 14,000rpm
- 6. Discard supernatant
- 7. Dry pellet in SpeedVac for up to 10min at low or medium heat
- 8. Re-suspend pellet in 1xTE or Qiagen AE buffer
 - a. Generally use 50-200ul depending on amount of DNA present
 - b. Thermomixer is probably ideal at RT