Small-scale protein prep in S. pombe

- 1. Grow yeast in liquid culture to desired density (we grow to 0.5 at OD_{595})
- 2. Pellet 20 ml by centrifuging at 3000 rpm for 10 minutes at 4°C (increasing the time up to 20 minutes or centrifuging at room temperature won't do any harm).

All subsequent steps should be carried out with cold solutions on ice.

- 3. Wash once in water and resuspend pellet in 200 ul protein lysis buffer (or other appropriate buffer) in a 1.5 ml screwcap tube.
- 4. Add 0.5 mm glass beads (we use acid washed glass beads from Sigma, 425-600 um) to lysis buffer meniscus.
- 5. Vortex aggressively at 4°C for 15 min using a multi-tube holder. The buffer should be yellow and most yeasts appear broken by microscopy (If you vortex for 15 minutes or more you may need not to confirm by microscopy).
- 6. Puncture screwcap tube bottom with a needle. The needle can by heated with a flame to facilitate puncturing.
- 7. Place the punctured tube in a fresh 1.5 ml tube.
- 8. Spin for 10 seconds at maximum speed in a microcentrofuge to drain lysate into the fresh tube. Then spin the collected lysate for 5 min at maximum speed in a microcentrofuge to pellet insoluble material.
- 9. Collect the supernatant in a new 1.5 ml tube. The Supernatant can be used as a clear cell lysate for subsequent biochemistry.
- 10. Quantify the protein concentration using your method of choice. We use the Bio- Rad DC protein assay

If you started with 20 ml of 0.5 OD_{595} cultured cells you may get 150 - 300 ug of protein or more.

Protein Lysis Buffer

50 mM Tris (pH 7.5)

150 mM NaCl

5 mM EDTA

10% Glycerol

1 mM phenylmethylsulphonylfluoride (PMSF) – added fresh

References

Forsburg and Rhind (2006). Yeast 23: 173-183