




## VWR Electrophoresis Systems

Horizontal or vertical

Molded for durability

Safe and reliable



## VWR Horizontal Electrophoresis Systems

- Molded, durable gel boxes and gel trays
- Rubber gates for tapeless gel casting
- EasyLift safety lid with buffer recapture system

Designed for ease of use and safe operation, the VWR Electrophoresis Systems and accessories are built to withstand the rigors of every day use.

All of the gel boxes and gel trays are molded of thick acrylic for leakproof performance. The gel trays are UV transparent for easy visualization of samples. The electrodes are protected by “cassettes” and are easily removed for cleaning and replacement.


Casting is simple and tape-free with the durable rubber casting gates. Slots in the sides of the trays allow for easy comb placement. The two piece design of the combs provides for vertical adjustment, giving the user control over the depth of the well. The combs supplied with the systems are 1.5 mm thick. Multiple combs, including those compatible with multichannel pipettors and different thicknesses, are available.

The “EasyLift” gel box lid is easily removed using the side tabs and pressure pads. The lid is domed, to prevent condensation from dripping onto the gel, and is surrounded by a drip ring, to help recover the condensate and maintain buffer concentration.

A variety of sizes are available, from the compact mini system to the large maxi system which will accommodate as many as 450 samples. A broad range of accessories are also available including various sizes of gel trays and combs.

	<b>Mini7</b>	<b>Mini10</b>	<b>Midi10</b>	<b>MidiPlus10</b>	<b>MidiPlus15</b>	<b>Maxi20</b>
Gel dimensions (W x L)	7 x 7 cm	7 x 10 cm	10 x 10 cm	15 x 10 cm	15 x 15 cm	20 x 20 cm
Unit dimensions (W x L x H)	9 x 21 x 9 cm	9 x 21 x 9 cm	12.5 x 22 x 9 cm	17.5 x 26.5 x 9 cm	17.5 x 26.5 x 9 cm	23 x 39.5 x 9 cm
Buffer volume	225 mL	225 mL	300 mL	500 mL	500 mL	1,200 mL
Maximum sample capacity	32	64	100	140	210	450
Supplied Combs	8 tooth (2 ea)	8 tooth (3 ea)	10* tooth (2 ea)	16* tooth (1 ea) 20 tooth (1 ea)	16* tooth (1 ea) 28* tooth (1 ea) 20 tooth (1 ea)	20 tooth* (2 ea) 40 tooth* (2 ea)
VWR Catalog Number	89032-288	89032-290	89032-292	89032-294	89032-296	89032-298
Price	\$250.00	\$260.00	\$275.00	\$320.00	\$420.00	\$440.00

\*Denotes multichannel compatible combs.



## VWR Mini 7

This true "mini-gel" system has a small footprint and is designed for quick runs of low to medium numbers of samples. Will also accept a 10 cm long tray, allowing for larger numbers of samples or longer runs. Supplied with one 7 x 7 cm gel tray, well visualization strips and two eight tooth combs



89032-288

## VWR Mini 10

Supplied with a 10 cm tray, this minigel systems maximizes the number of samples that can be run in a unit with a small footprint. Ideal for labs with medium throughput or for longer runs of smaller samples. Supplied with one 7 x 10 cm gel tray, well visualization strips and three 8 tooth combs.



89032-290

## VWR Midi 10

Extending only in width from the gel box above, this system will hold more samples with a relatively low increase in buffer volume. The Midi 10 is perfect for processing large numbers of samples in short to medium runs. Supplied with one 10 x 10 cm tray, well visualization strips and two 10 tooth multichannel compatible combs.



89032-292

## VWR Midi Plus 10 and 15

These two versatile systems are ideal for resolving close molecular weight samples and restriction fragments. Their multichannel combs make loading large numbers of samples quick and easy. Supplied with either a 15 x 10 or 15 x 15 cm gel tray and well visualization strips. The 10 cm unit is supplied with one 16 tooth and one 20 tooth comb. The 15 cm unit is supplied with one each of the 16, 20 and 28 tooth combs.



With 15 x 10 cm tray  
89032-294

With 15 x 15 cm tray  
89032-296

## VWR Maxi 20

For laboratories processing large numbers of samples, such as PCR products from 96 well plates, RFLP analysis or plasmid screenings, the Maxi 20 is an ideal choice. The 20 cm width of the gel provides for high throughput screening of samples, while the 20 cm length provides for good resolution. The unit is supplied with one 20 x 20 cm gel tray, well visualization strips and two 20 and 40 tooth combs, both of which are multichannel compatible.



89032-298

## VWR Power Source Power Supplies

- Perfect partner for VWR Gel Boxes
- Designed for safety and performance
- Constant current, voltage or power modes

VWR Power Source Power Supplies meet a wide range of applications including DNA, RNA and protein electrophoresis, blotting and many others. Designed with safety in mind, they are loaded with the features that today's molecular biology applications demand.

Two models are available with voltage ranges to 300V and current ranges to 3000mA. Operating modes include constant voltage or constant current with automatic crossover. In addition, the Model 250V offers a constant power operating mode. Operation may be timed or continuous.


All parameters, including voltage, current, power (where applicable) and time are easily set using the membrane keypad and displayed digitally on the back lit LCD. Programming is straightforward and simple, and allows setting up multiple voltage/current changes within a single program. Each unit includes four sets of color coded output terminals allowing multiple gels to be run from the same power supply.

VWR Power Source Power Supplies offer short circuit and over current/over voltage protection. No load and load change detection systems are also incorporated into the software for safety. Programming also includes an option for automatic recovery after a power failure.

Compact in design with a small footprint, the VWR Power Source Power Supplies are lightweight and stackable to save space. For easy viewing and setting of parameters on the bench or on a shelf, the front of the power supply can be raised with the integral folding feet.

### Specifications

Model	300V	250V
Output Voltage Range/Increments	2-300V/1V	5-250V/1V
Output Current Range/Increments	4-500mA/1mA	10-3000mA/10mA
Output Power Range/Increments	90 watts max/na	1-300W/1W
Operating constant modes	current or voltage	current or voltage or power
Timer	1-999 minutes or continuous	1 minute to 99 h 59 min or continuous
Programmable	No	Yes
Dimensions (W x D x H)	19 x 25 x 8 cm	19 x 25 x 8 cm
Weight	2.2 kg/4.9 lbs	2.5 kg/5.6 lbs
Input Voltage	120V	120V
VWR Catalog Number	93000-744	93000-746
Price	\$378.75	\$756.49





## VWR Modular Mini Vertical PAGE Systems

- PAGE, Electroblotting and 2D Electrophoresis Modules
- Exceptionally easy gel casting
- Customizable accessories

Why purchase three different systems for protein electrophoresis, IEF and blotting when the VWR Modular Mini Vertical PAGE System will do them all. The system features inserts for each application and a common buffer tank in which to run them.

Each insert has its own electrode assembly, which connects to the lid of the buffer tank. For safety, the power leads connect directly to the lid - when the lid is removed, power is disconnected from the system. The inserts also have a small pad that fits into position in the lid. This ensures that the lid is always placed on properly and aids in its removal. The buffer tank is molded to prevent separation and leaking.

Inserts, accessories and running tanks are manufactured from quality materials for a long life. A variety of spacers, combs and other accessories are available to customize the system. All spacers and combs are color coded or clearly labeled to indicate thickness. Each system is supplied with a cooling pack for superior results during longer runs or quick runs at higher voltages. Three packages are available to meet your needs. Components may also be purchased separately.

	Mini PAGE System	Mini PAGE Blotting System	Mini PAGE Tube Gel System
Gel Dimensions	10 x 10 cm	NA	8 cm long
Capacity	2 gels/run*	3 gels/transfer	10 tubes/run
Maximum samples	20 per gel	NA	10 (1/gel)
Buffer tank dim (W x D x H)	19 x 13 x 15 cm	19 x 13 x 15 cm	19 x 13 x 15 cm
Buffer volume	250 mL	1200 mL	500 mL
Included items**	2 ea 10 x 10 cm notched plates 2 ea 10 x 10 cm plates w/spacers 2 ea 1 mm comb, 12 wells PAGE Module and casting base	3 ea blotting cassettes 6 ea fiber pads Blotting module	20 ea capillary tubes 20 ea tube supports 10 ea blocking plugs 2D module
VWR Catalog Number	89032-300	89032-302	89032-304
Price	\$495.00	\$695.00	\$725.00

\*It is possible to run 4 gels by stacking the plates.

\*\*Each system includes buffer tank, lid with leads and cooling pack. The Blotting and 2D Gel Systems also include the PAGE components.

Molded gel box for leakproof performance

Pins to orient lid and aid in removal

Easy to use casting stand

10 gel capacity

Tube Gel Module

PAGE Module

Capacity for 3 gels

Gaskets for easy insertion/removal

Wide variety of combs and spacers

Easy to set pressure pads


Run up to four 10 x 10 gels

High intensity, coiled electrodes

Blotting Module

## VWR Mini Page System

Set up of 10 x 10 cm gels is easy with the VWR Mini PAGE System. Spacers are bonded to the glass plates. Before pouring the gel, the plates are placed in the gel insert and held in place by pressure strips. The insert is then placed in the casting stand. After the gel has polymerized, the insert is simply removed from the casting stand and placed in the buffer chamber. Each side of the insert will hold one set of plates. A third plate may be added to each set, allowing for a total capacity of four gels. The system is supplied with buffer tank, gel insert, casting stand, two notched plates, two plates with bonded spacers, two 12 tooth combs and a cooling pack.



89032-300



89032-302

## VWR Mini PAGE Electroblotting System

Three gels can be blotted quickly and easily with the VWR Mini PAGE Electroblotting System. Gels are sandwiched snugly in the blotting cassettes which are placed into the blotting insert. Coiled electrodes within the insert generate a high intensity current, allowing a standard transfer to take place in about an hour. The system includes the blotting insert, three blotting cassettes, six fiber pads and all of the components of the VWR Mini PAGE system listed above.

## VWR Mini PAGE Tube Gel System

The insert included in the VWR Mini PAGE Tube Gel System has a capacity for 10 first dimension gels. Gel ports are gasketed for easy insertion and removal of tubes. Included in the system are the tube gel insert, 10 blanking caps, 20 disposable capillary tubes, 20 tube gel supports and all of the components of the VWR Mini PAGE system above.



89032-304

**VWR** 

**vwr.com 1.800.932.5000**

**Order from Over 1,000,000 Products**

Prices, product appearance and specifications are current at the time of printing, subject to change without notice. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. All prices are in U.S. dollars unless otherwise noted. Offers valid in USA, void where prohibited by law or company policy, while supplies last. Visit [vwr.com](http://vwr.com) to view our privacy policy and additional disclaimers.

©2006 VWR International, Inc. All rights reserved. Printed in U.S.A.