# A Novel Approach to Optimize Sensitivity in Ultrasound Molecular Imaging

Mehmet Kaya

Esra Talu, Shukui Zhao, Kanaka Hettiarachchi, Abraham. P. Lee, Marjorie L. Longo, and Paul A. Dayton

UNC/NCSU Joint Department of Biomedical Engineering




Department of Chemical Engineering, University of California, Davis Department of Biomedical Engineering, University of California, Irvine

# Background

- Ultrasonic Molecular Imaging
  - Uses targeted contrast agents
  - Attachment of specific ligands
    - Antibodies
    - Peptides
    - Polysaccharides
  - Bind to cellular receptors
  - Recognize disease antigens
  - Targeted Pathologies
    - Tumors
    - Inflammation
    - Thrombus


## Targeted Ultrasound Contrast Agents

- Blood doesn't scatter US well
- Microbubble: A contrast agent
- Filled with gas: High echogenicity compared to liquid based particles or blood cells
- Diameters of the order of 1 to 5 μm
- Smaller than Red Blood Cells (6-8 µm)
- Improve the sensitivity of 2-D and 3-D ultrasound imaging




#### Structure of Microbubbles

- A monolayer amphiphilic phospholipid shell
- The water insoluble hydrocarbons
- Charged phosphate head groups
- Use of PEG:
  - Coalescence
  - Phagocytosis
  - Physical barrier to enzymatic agents, adsorption of blood proteins




## **Current Contrast Agent Production**

- produced by mechanical agitation of lipid solution
- polydisperse size distribution



# Problems with Current Contrast Agents & Ultrasound Systems for Molecular Imaging


- Size distribution
- The resonant frequency




- Limited frequency bandwidth
- Only a small percentage of the contrast agents have diameters which result in resonant frequencies optimized for the limited bandwidth
- Fewer contrast agents retained
- Filtering does not work on lipid-shelled bubbles

The Problem

#### A Solution





bubble size=1.8±1.5 microns
F receiver bandwidth=5.0 +/- 1.5 MHz

detectable bubbles: 18%


bubble size=1.8±0.2 microns


F receiver = 5.0 +/- 1.5 MHz

detectable bubbles: 92%

## Monodisperse Microbubble Production


• Use of Microfluidic Flow Focusing Device







#### Why use a Microfluidic System?

• Relatively cheap, simple, and easy to multiplex







(video playback at 1/100 speed)




**50** μm


#### Tailoring the Size Distribution Microbubbles






Lipid Encapsulation for Stability







Cell Targeting




**Incorporate Targeting Ligand** 

#### Ultrasound Imaging System








Microtube




US Transducers


- Ultrasound Transducers
- Transmitter 2.5 MHZ/Receiver 5 MHz
- Liquid flow rate: 35 ul/min
- PFC flow rate: ~2-4 psi




# Acoustic Studies of Monodisperse Contrast Agents


- Standard deviation of echo amplitude is smaller
- Correlation coefficient 0.95 vs. 0.70
- Improved Signal Response of Monodisperse Contrast Agents







# Detection of Monodisperse Agents



# **Summary & Conclusions**

- Successful production of monodisperse contrast agents
- Monodisperse microbubble contrast agents have the potential to increase sensitivity in ultrasound molecular imaging
- More sensitive detection of monodisperse contrast agents compared to polydisperse agents
- Production rate of contrast agents is one of the challenges (in comparison to ~10<sup>10</sup> µl with polydisperse)

# Funding Acknowledgement: NIH Roadmap for Medical Research, Grant EB005325