


MIT BioMicro Center: Sample Quantification



Katherine Sullivan, Justin Elliott and Stuart Levine


MICROFLUIDICS








Agilent 2100 Bioanalyzer and Caliper LabChip GX



RT-PCR





providing Thing 3.

General Information

- •The MIT Biomicro Center now offers three Roche Lightcycler 480s, a fully-integrated, high throughput real-time PCR platform for highly accurate qualitative and quantitative detection of nucleic acids. •96 well and 384 well plates
- Rapid cycle times

m/m

•Can detect >4 fluors

Melting Curve Genotyping

Temperature

Illumina QC

Cycles


Standard Curve

Amplification Curve

SYBR green.


Applications:






Absolute Quantification Relative Quantification


 Compare expression levels of 2 target sequences in a single




FLUIDIGM BIOMARK

Massively parallel real-time PCR machine





Integrated Fluid Circuit (IFC)





- Dense network of channels able to regulate flow of solution down to picoliter scale
- Nanoflex valves open and close by application of air pressure
- Reaction volumes fixed by chamber volume, consistent.

Digital PCR

- 9,216 qPCR data points
- Requires pre-amplification
- 1/200th of typical reagents:

qPCR	384 well	96.96 Dynamic Array
Master Mix	46,000μL	240µl
Primer-probe	9,200μL	480µl
Plates	24	1
Pipette steps	18,432	192
	Master Mix Primer-probe Plates	Master Mix46,000μLPrimer-probe9,200μLPlates24







Other Applications







Expected Ratios



Allele X Frequency 9,216 data points from a single 96.96 dynamic array

CALIPER v. BIOANALYZER: DNA

The Caliper Labchip GX and Agilent Bioanalyzer both provide a platform for analyzing DNA samples at high sensitivity. We preformed several validation assays to compare the performance of the two platforms on samples prepared for Illumina sequencing.


Manufacturers' Specifications

	BIOANALYZER	LABCHIP
Size range	50-7000bp	50-5000bp
Linear Range	5-500 pg/μl	10-500 pg/μl
Accuracy	+/- 20%	+/- 30%
Batch Size	11 samples	96 samples

ILLUMINA LIBRARY QUALITY CONTROL


Evaluation of libraries prior to clustering to ensure consistent sequencing results.

Bioanalyzer




The Bioanalyzer or Labchip GX confirms there is sample, and provides the average size, approximate concentration and molarity. It also can identify samples with problems such as hyperamplification, "primer-dimer" or low complexity


Cross Comparison





Reproducibility

10 DNA samples run six times on the Lab Chip and twice on the BioAnalyzer. Error bars indicate the standard deviation, a box plot of the errors is shown below.










Lab Chip GX

NEB 50bp ladder and 21 samples were compared on both the Bioanalyzer and the LabChip. Results indicate the two instruments have very similar behaviors as expected.

RT PCR




Provides accurate concentration measurement for clustering. Absolute quantification based on protocols developed at the Broad Institute.



The Ct of library samples is compared to that of a standard curve

Multiplexing

Data obtained from RT PCR is used to pool samples of varying concentrations for equal or biased representation in multiplexed lanes.



Low complexity