




| Laboratory of Integrated Bio Medical | Micro/Nanotechnology & Applications





Interfacing Biology and Engineering at the Micro and Nano Scale

March 8th, 2008 Institute of Biological Engineering 2008 Annual Conference

Rashid Bashir

Micro and Nanotechnology Laboratory
Electrical and Computer Engineering & BioEngineering
University of Illinois, Urbana-Champaign
http://libna.micro.uiuc.edu/



Acknowledgements



Researchers:

- Prof. Demir Akin
- Dr. JeongMi Moon
- Piyush Bajaj
- Vincent Chan
- Brian Dorvell
- Oguz Elibol
- Yi-Shao Liu
- Kidong Park
- Bobby Reddy
- Shuaib Salamat
- Murali Venkatesan
- Nick Watkins











Recent Alums

Prof. J. Jang, Prof. S. Iqbal, Prof. S.
 W. Lee, Prof. L. Yang,

BioVitesse, Inc. Co-Founder




Funding Agencies

- National Institute of Health, NIBIB
- NASA Institute on Nano-electronics and Computing (INAC)
- US Department of Agriculture (Center for Food Safety Engineering at Purdue)
- NSF NSEC at OSU
- NIH Nanomedicine NDC


Faculty Collaborators

- Prof. A. Alam (ECE, Purdue)
- Prof. D. Bergstrom (Med Chem, Purdue)
- Prof. A. Bhunia (Food Science, Purdue)
- Prof. S. Claire (IU-SOM)
- Prof. P. Guo (Vet Medicine, BME, Univ. of Cincinnati)
- Prof. M. Ladisch (Ag& Bio Engr, Purdue)
- Prof. M. Toner (Harvard Med School)
- Prof. J. P. Robinson (BMS, BME, Purdue)
- Prof. W. Rodriguez (Harvard Med School)





Evolution of Technologies





Research Projects

- Petri dish on a chip with electrical detection
- Cantilevers for virus and bacterial detection
- Detection of CD4+ cells from blood
- Nanopores for single DNA molecule characterization
- Label free field effect silicon sensors for DNA and proteins

Micro/ Nano-biosensors and Lab-on-Chip

Nanoscale
Therapeutics and
Nanomedicine

Cell Biology and Tissue Engineering


- Bacterial Mediated Delivery of Nanoparticles
- Macrophage Stealth Delivery of Nanovectors
- Phi-29 Nanomotor for Nanomedicine

- 3-D Tissue Engineering with Stereo-Lithography
- Nanomechanics for Cell Biology



Microfluidic-based technology platforms for bacterial detection


L. Monocytogenes
On a chip




- R. Bashir, ECE, BME
 A. Bhunia, Food Science
 M. Ladisch, BME, ABE
 J. P. Robinson, BMS, BME
- Rapid detection of live bacteria has wide applications
 - Industrial Microbiology (Pharmaceutical facility monitoring, Food Safety)
 - Blood transfusion
 - Human Diagnostics
 - Global Health
- Rapid time to result limited by:
 - Amplification, enrichment, culture
 - 48-72 hours for growth and identification
- Center for Food Safety Engineering at Purdue (<u>www.cfse.purdue.edu</u>) through cooperative agreement with USDA-ARS















Integrated Chips for Study of Microorganisms and Cells




















Silicon Nanowires and






Biotinylated

Antibody (C11E9)

Streptavidin

Biotinylated BSA

Overall Detection Scheme



* Size information are obtained from Biochemistry 2nd 9 edition, R. H. Garrett and C. M. Grisham, 1995.

SiO₂ surface


Gomez, Morisette, Bashir, IEEE/ASME JMEMS, 2005 Yang, Li, Akin, Huang, Bhunia, Ladisch, Bashir, Lab Chip, 2006





Bacterial Growth & Impedance Microbiology!



Invented in late 1800s (Petri and Koch)
Still the most widespread means to
grow and detect the presence of








Yang, L, Banada, P., Bhunia, A. Bashir, R., Biotech Bioengr, 2005




Petri Dish-on-a-Chip






Bacterial Cell Concentration on-Chip


 DEP-based concentration system collects particles from a large flow stream and diverts them to a smaller stream







On-Chip Incubation of *L. Monocytogenes*











BioVitesse, Inc.




- 'Petri dish on a chip' to miniaturize impedance microbiology
- To quickly and reliably detect and identify <u>live bacteria</u> in 2 to 4 hours, instead of 2 to 10 days
- Provides in-process quality control monitoring systems
- To the industrial microbiological market (Bio/Pharma and Food Safety)

Automated System





Chip Cartridge









Sample, Media, Sanitizing Fluid


BioVitesse, Inc. Company Confidential





DNA - the code of life in all Cells!












Molecular Biology of the Cell, Alberts, Bray, et al.





Nanopores for single DNA molecule characterization

$2c/bp \times 3e9bp = $60M!$



- α-hemolysin channel, a biological protein based-pore
- Diameter of 2.6 nm.
- Translocation of ssDNA using electrophoresis
- Nanoscale coulter counter






The model of DNA passing through an α -hemolysin channel.




Nanopore Fabrication

Key steps

- EBL to write dots
- TMAH etch
- Pore in SOI
- Oxidation
- Shrinking in TEM




Pore shrinking temporal profile of a nanopore channel All Images at 1,000,000X. (NPC-2)






Adding DNA Selectivity to Solid State Nanopores

Protein Ion Channels




Nanopores with selectivity towards specific targets




[Kasianowicz et al., 2006]



Selective Nanopores Measurements Overview









Single-Base Mismatch Selectivity

- Channel-Molecular Interaction
 - MM-DNA vs. PC-DNA
- PC-DNA
 - Interactions with Binding Sites
 - Faster and More than MM-DNA
- MM-DNA
 - Electrostatic Friction
 - Mechanical Resistance
 - Inability to open HPL



$$J = \frac{n}{\tau}(c_1 - c_2)$$






Channel Interaction Attractive vs. Repulsive Potential

Key Assumptions

- PC/HPL: Attractive Potential
- MM/HPL: Repulsive Potential
- Magnitudes of Potentials
- ϕ span part of the channel

$$rac{ au_{MM}}{ au_{PC}} = e^{|\phi_{MM}| - |\phi_{PC}|}$$




(a) Attractive potential (b) Repulsive potential, spanning part of the channel



Nanopore Channel Sensor Array for Characterization of Single Molecules

- Miniature DNA analyzer
 - Transduction events at single molecule, nanoscale
 - Interface to the macroscale
 - Array of nanopores serial or parallel ?
 - Site selective functionalization strategy
- Sequencing by digestion
 - Smaller molecules bases!
 - Digest and detect





Microcantilever Mass Sensors

Unloaded Resonant Frequency:

$$f_0 = \frac{1}{2\pi} \sqrt{\frac{k}{m *}}$$

Spring constant for a rectangular shaped cantilever beam: $k = \frac{Et^3w}{4l^3}$


Loaded Resonant frequency : $f_1 = \frac{1}{2\pi} \sqrt{\frac{k}{m*+\delta m}}$ δm is the added mass

$$\Delta m = \frac{k}{4\pi^2} \left(\frac{1}{f_1^2} - \frac{1}{f_o^2} \right)$$

- k = spring constant
- m = mass of cantilever
- f_0 = unloaded resonant frequency
- $f_1 = loaded$ resonant frequency

Mass Change Detection








Nanomechanical Sensors for Viral Detection


Objectives:


To develop technology for the rapid detection of virus particles in fluid and air using Nanomechanical Cantilever Sensors


Frequency Shift, $\Delta f = 60 \text{ kHz}$


- \Rightarrow Mass change, $\Delta m = 9 \text{ fg}$
 - \Rightarrow This corresponds 1 vaccinia virus.















25


A. Gupta, et al. Proc. Nat. Acad. Sci. 2006



Living Cantilever Arrays












Cell Attachment on Cantilevers

DEP Assisted Capture










Hela cells, captured with DEP. Conductive silicon electrode generate DEP forces to capture Hela cells.





Living Cantilever Arrays


















Living Cantilever Arrays














Calculated Mass: 2.06 ng

Volume: 2349 um³ Density³: 1.055g/cc Estimated mass : 2.48 ng





Point of Test Biosensors




Silicon Based BioChips!

- Disposable, one-time-use devices
- Rapid, sensitive, integrated
- Detection and monitoring of disease and state of health










- Focus on translation: Top-down integrated nano-structures and nanosensors platforms – POC sensors for biomolecules
- Chip makers and diagnostic industry
- System level fabrication and modeling/simulations (from 'molecules to systems') with considerations of integrated systems, fluid and biomolecule transport, packaging and interfaces, etc.




BioMEMS and **Bionanotechnology**

Apply micro-systems and nanotechnology to develop novel devices and systems that have a biomedical impact or are bioinspired





Minimum Detectable Mass



The frequency measurement is limited by thermo-mechanical noise on the cantilever beam.

Minimum Detectable Frequency,

$$\Delta f_{,\min} = \frac{1}{A} \sqrt{\frac{f_0 k_B TB}{2 \pi kQ}}$$

Minimum Detectable Mass,

$$\Delta m_{,min} = \frac{1}{A} \sqrt{\frac{4 k_B TB}{Q}} \frac{m_{eff}^{5/4}}{k^{3/4}}$$

k_B = Boltzmann constant

T = Temperature in Kelvin

B = Bandwidth measurement, (= $2\pi f_0/2Q$)


100kDa Protein ~ 10⁻¹⁹ g DNA bp ~ 10⁻²¹ g

Q can increase by 100X by driving the cantilevers

Roukes, et. al.




40nm NPs in KB Cancer Cells






200nm Beads in KB Cells



- Bead incubation for 1 hour
- Confocal images NPs are NOT internalized
- Cross-sectional scans







Cell Concentration and Recovery

- Proprietary technology for cell concentration and recovery in small volumes – CCR™ kit
- One-time use filter cartridge
- Recovery rates ~ 50%











Recovery Results

E. coli

100 ml sample in PBST

Run#	Initial cell conc. (cells/ml)	Final pressure drop (psi)	Recovery collection 1 (viable cells)	Total cells in collection 1	Recovery collection 2 (viable cells)	Recovery collection 3 (viable cells)	Leakage through permeate side	Total recovery (3 collections)
1	52±12	46	72.4±11.6%	3767±603	4.9±0.9%	0.4±0.1%	0.0±0.0%	77.7%
2	52±12	40	62.2±15.5%	3233±808	2.9±1.6%	0.5±0.3%	0.0±0.0%	65.5%
3	52±12	38	62.8±14.7%	3267±1531	4.0±1.1%	0.2±0.1%	0.0±0.0%	67.0%

Listeria innocua

Run#	Initial cell conc. (cells/ml)	Final pressure drop (psi)	Recovery collection 1 (viable cells)	Total cells in collection 1	Recovery collection 2 (viable cells)	Recovery collection 3 (viable cells)	Leakage through permeate side	Total recovery (3 collections)
1	66±22	31	56.6±3.8%	3510±236	12.7±2.4%	0.9±0.4%	0.0±0.0%	70.2%
2	66±22	30	60.0±3.9%	3717±244	13.7±2.4%	2.3±0.7%	0.0±0.0%	75.9%
3	66±22	28	37.2±3.9%	2303±243	25.1±2.6%	6.0±0.6%	0.0±0.0%	68.2%

Streptococcus faecalis

Run#	Initial cell conc. (cells/ml)	Final pressure drop (psi)	Recovery collection 1 (viable cells)	Total cells in collection 1	Recovery collection 2 (viable cells)	Recovery collection 3 (viable cells)	Leakage through permeate side	Total recovery (3 collections)
1	109±24	40	53.7±1.9%	5867±208	5.6±2.4%	1.7±0.3%	0.0±0.0%	61.0%
2	109±24	40	37.3±8.3%	4067±907	17.3±0.9%	2.4±0.3%	0.0±0.0%	56.9%
3	109±24	40	56.2±16.6%	6133±1815	6.4±0.3%	2.9±0.5%	0.0±0.0%	65.5%





Targeted Nano-therapeutics


Objective:

To develop novel means to deliver nanotherapeutics to cancer cells and tumors

- Current drugs are ineffective at penetrating to the tumors due to inadequate vascularization at these tumor sites.
- Nanoparticles have the similar challenge of reaching tumor sites if adequate vascularization does not exist.
- Penetrate multiple barriers
- Need active transport to the cells and tumors







- Current drugs are ineffective at penetrating through the solid tumors due to inadequate vascularization at these tumor sites.
- Nanoparticles have the similar challenge of reaching tumor sites if adequate vascularization does not exist since they are carried by blood flow.

Bacterial Mediated Delivery of Nanoparticles in Cells

D. Akin, J.P. Robinson, S. Muhammad, A. Bhunia, R. Bashir

Objective:


To develop biologically-inspired intelligent devices for intracellular monitoring and therapeutic intervention of signal transduction networks


Bacteria as a carrier!

- Intracellular bacteria may provide means to penetrate and deliver therapeutic compounds to solid tumors that are otherwise inaccessible to drugs delivered by oral and parenteral routes.
- Bacteria can be engineered to have affinity for certain tumors.
- Bacteria can be engineered to express proteins on their surface for therapeutic or diagnostic intervention




Functionalization Strategy








Nanoparticles on Bacteria






Bacterial Mediated Delivery of NP







Expression of GFP in KB Cells

Expression in-vitro in cells







Expression in-vitro in mice










Work featured in Popular Science, EE Times, Nature News, and CNBC TV Show



A. Gupta, D. Akin, R. Bashir, Applied Physics Letters, 2004. A. Gupta, et al. Proc. Nat. Acad. Sci. 2006