

Raw data

Sample Name: 78%BSA by mass	X-ray	: Cu / 30 kV / 15 mA	Counter	: MiniFlex2 counter
File : 78%BSA by mass trial 2 (wh	Goniometer	: MiniFlex2 goniometer		
Comment :	Attachment	: Miniflex 6 sample changer		
Date : April-09-14 10:37:02	Filter	: Unused	Scan mode	: Continuous
Operator : Administrator	I.Monochro	: Unused	Scan speed	: 1.000 deg./min.
	C.Monochro	: Fixed Monochromator	Sampling width	: 0.030 deg.
	DivSlit	: 1.25 deg.	Scan axis	: 2theta/theta
	SctSlit	: 1.25 deg.	Scan range	: 10.000 -> 70.000 deg.
	RecSlit	: 0.3mm	Theta offset	: 0.000 deg.
	Monochro RS	: 0.8mm		

Memo :

